

Integrating the Healthcare Enterprise



5

IHE Patient Care Coordination (PCC) Technical Framework Supplement

CDA Content Modules

10

15

Trial Implementation

20

Date: August 31, 2012
Author: IHE PCC Technical Committee
25 Email: pcc@ihe.net

Foreword

- This is a supplement to the IHE Patient Care Coordination Technical Framework V8.0. Each supplement undergoes a process of public comment and trial implementation before being incorporated into the volumes of the Technical Frameworks.
- 30 This supplement is different than traditional IHE supplements. It serves as the Trial Implementation staging area for content modules. The content modules (section level, entry level) defined during Trial Implementation are gathered in this document to provide a central location for readers of supplements from PCC, QRPH and/or other domains that use the dictionary of content modules first defined by PCC. After individual modules are successfully tested and reviewed, they will be moved to Final Text. At that time, they are removed from this document. Thus, this supplement will continue to exist for some time as new content modules are defined and documented here. Likewise, content modules will be removed as they go to final text.
- 35 40 This supplement describes changes to the existing technical framework documents and where indicated amends text by addition (**bold underline**) or removal (**~~bold strikethrough~~**), as well as addition of large new sections introduced by editor's instructions to "add new text" or similar, which for readability are not bolded or underlined.
- 45 "Boxed" instructions like the sample below indicate to the Volume Editor how to integrate the relevant section(s) into the relevant Technical Framework volume:

Replace Section X.X by the following:

- General information about IHE can be found at: www.ihe.net
- 50 Information about the IHE PCC domain can be found at: <http://www.ihe.net/Domains/index.cfm>
Information about the structure of IHE Technical Frameworks and Supplements can be found at: <http://www.ihe.net/About/process.cfm> and <http://www.ihe.net/profiles/index.cfm>
- The current version of the IHE Technical Framework can be found at:
http://www.ihe.net/Technical_Framework/index.cfm
- 55

CONTENTS

	INTRODUCTION	15
60	PROFILE ABSTRACT	15
	OPEN ISSUES AND QUESTIONS	15
	CLOSED ISSUES.....	15
	VOLUME 1 – INTEGRATION PROFILES	16
	GLOSSARY	17
	2.5 HISTORY OF ANNUAL CHANGES	17
65	VOLUME 2 – TRANSACTIONS AND CONTENT MODULES	18
	6.1 CONVENTIONS	18
	6.2 FOLDER CONTENT MODULES.....	19
	<i>6.2.1 EDES Folder Specification</i>	19
	<i>6.2.2 APR Folder Specification.....</i>	19
70	<i>6.2.3 LDR Folder Specification</i>	19
	6.3 HL7 VERSION 3.0 CONTENT MODULES	19
	<i>6.3.1 CDA Document Content Modules</i>	19
	6.3.1.X History and Physical Specification 1.3.6.1.4.1.19376.1.5.3.1.1.16.1.4.....	19
	6.3.1.x.1 Format Code	19
75	6.3.1.x.2 LOINC Code.....	19
	6.3.1.x.3 Standards.....	19
	6.3.1.x.4 Specification	20
	6.3.1.x.5 Conformance.....	20
	<i>6.3.2 CDA Header Content Modules</i>	22
80	6.3.2.1 Language Communication 1.3.6.1.4.1.19376.1.5.3.1.2.1	22
	6.3.2.2 Employer and School Contacts 1.3.6.1.4.1.19376.1.5.3.1.2.2	22
	6.3.2.3 Healthcare Providers and Pharmacies 1.3.6.1.4.1.19376.1.5.3.1.2.3	22
	6.3.2.4 Patient Contacts 1.3.6.1.4.1.19376.1.5.3.1.2.4	22
	6.3.2.5 Spouse 1.3.6.1.4.1.19376.1.5.3.1.2.4.1.....	23
85	6.3.2.5.1 Parent Template	23
	6.3.2.5.2 Specification	23
	6.3.2.5.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4/'><templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4.1'>.....	23
	6.3.2.5.4 <associatedEntity classCode='PRS'>	23
90	6.3.2.5.5 <code code='127848009 184142008' displayName=' ' codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT'>.....	23
	6.3.2.5.6 Completed Example.....	24
	6.3.2.6 Natural Father of Fetus 1.3.6.1.4.1.19376.1.5.3.1.2.4.2	24
	6.3.2.6.1 Parent Template	24
95	6.3.2.6.2 Specification	24
	6.3.2.6.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4/'><templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4.2'>....	25
	6.3.2.6.4 <associatedEntity classCode='PRS'>	25
	6.3.2.6.5 <code code='xx-fatherofbaby' displayName=' ' codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT'>	25
	6.3.2.6.6 Completed Example.....	25
100	6.3.2.7 Authorization 1.3.6.1.4.1.19376.1.5.3.1.2.5	26
	6.3.2.7.1 Parent Template	26
	6.3.2.7.2 Specification	26
	6.3.2.7.3 <authorization typeCode='AUTH'>	26
	6.3.2.7.4 <consent classCode='CONS' moodCode='EVN'>	26
105	6.3.2.7.5 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.5'>	26
	6.3.2.7.6 <id root=' '/>	26

	6.3.2.7.7 <code code=' ' codeSystem=' ' codeSystemName=' ' displayName=' '/>	27
	6.3.3 CDA Section Content Modules.....	27
110	6.3.3.1 Reasons for Care	27
	6.3.3.1.1 Reason for Referral	27
	6.3.3.1.2 Coded Reason for Referral	27
	6.3.3.1.3 Chief Complaint.....	27
	6.3.3.1.4 Hospital Admission Diagnosis	27
	6.3.3.1.5 Proposed Procedure Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.1	27
115	6.3.3.1.6 EBS Estimated Blood Loss Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.2	28
	6.3.3.1.7 Proposed Anesthesia Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.3	29
	6.3.3.1.8 Reason for Procedure Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.4	30
	6.3.3.1.9 Reason for Visit Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.1.1	31
	6.3.3.1.10 Injury Incident Description Section 1.3.6.1.4.1.19376.1.5.3.1.1.19.2.1	31
120	6.3.3.2 Other Condition Histories	32
	6.3.3.2.1 History of Present Illness	32
	6.3.3.2.2 Hospital Course.....	32
	6.3.3.2.3 Active Problems.....	32
125	6.3.3.2.4 Discharge Diagnosis	32
	6.3.3.2.5 History of Past Illness	32
	6.3.3.2.6 Encounter Histories.....	32
	6.3.3.2.7 History of Outpatient Visits	32
	6.3.3.2.8 History of Inpatient Visits	32
	6.3.3.2.9 List of Surgeries	32
130	6.3.3.2.10 Coded List of Surgeries.....	32
	6.3.3.2.11 Allergies and Other Adverse Reactions	32
	6.3.3.2.12 Family medical History	33
	6.3.3.2.13 Coded Family Medical History	33
	6.3.3.2.14 Social History Section	33
135	6.3.3.2.15 Functional Status.....	33
	6.3.3.2.16 Review of Systems.....	33
	6.3.3.2.17 Hazardous Working Conditions	33
	6.3.3.2.18 Pregnancy History.....	33
	6.3.3.2.19 Medical Devices.....	33
140	6.3.3.2.20 Foreign Travel.....	33
	6.3.3.2.21 Pre-procedure Family Medical History Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.5 (Deprecated)	33
	6.3.3.2.22 Coded Functional Status Assessment Section 1.3.6.1.4.1.19376.1.5.3.1.1.12.2.1	33
	6.3.3.2.22.1 Standards	34
	6.3.3.2.22.2 Parent Template	34
145	6.3.3.2.23 Pain Scale Assessment Section 1.3.6.1.4.1.19376.1.5.3.1.1.12.2.2	35
	6.3.3.2.24 Braden Score Section 1.3.6.1.4.1.19376.1.5.3.1.1.12.2.3	36
	6.3.3.2.25 Geriatric Depression Scale Section 1.3.6.1.4.1.19376.1.5.3.1.1.12.2.4	37
	6.3.3.2.26 Physical Function Section 1.3.6.1.4.1.19376.1.5.3.1.1.12.2.5	38
	6.3.3.2.26.1 Constraints	38
150	6.3.3.2.27 Preprocedure Review of Systems Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.13	42
	6.3.3.2.28 Estimated Delivery Date Section 1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.1	43
	6.3.3.2.29 History of Tobacco Use Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.8	44
	6.3.3.2.30 Current Alcohol/Substance Abuse Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.10	45
	6.3.3.2.31 History of Blood Transfusion Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.12	45
155	6.3.3.2.32 Anesthesia Risk Review of Systems Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.14	46
	6.3.3.2.33 Implanted Medical Device Review Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.46	47
	6.3.3.2.34 Pregnancy Status Review Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.47	47
	6.3.3.2.35 History of Infection Section 1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1	48
	6.3.3.2.36 Coded Social History Section 1.3.6.1.4.1.19376.1.5.3.1.3.16.1	48
160	6.3.3.2.37 Coded History of Infection Section 1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1	49
	6.3.3.2.38 Prenatal Events Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.2	49
	6.3.3.2.39 Labor and Delivery Events Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.3	51
	6.3.3.2.40 Newborn Delivery Information Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.4	52
	6.3.3.2.41 Postpartum Hospitalization Treatment Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.7	55

165	6.3.3.2.42 Event Outcomes Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.9	57
	6.3.3.2.43 Newborn Status at Maternal Discharge 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.8	57
	6.3.3.2.44 History of Surgical Procedures Section 1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2	58
	6.3.3.2.45 Operative Note Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.6	58
170	6.3.3.2.46 Child Functional Status Assessment 1.3.6.1.4.1.19376.1.7.3.1.1.13.3	59
	6.3.3.2.47 Psychomotor Development Section 1.3.6.1.4.1.19376.1.7.3.1.1.13.4	59
	6.3.3.2.48 Eating and Sleeping Assessment Section 1.3.6.1.4.1.19376.1.7.3.1.1.13.5	60
	6.3.3.2.49 Coded Event Outcomes 1.3.6.1.4.1.19376.1.7.3.1.1.13.7	61
	6.3.3.2.50 Intentionally blank	62
175	6.3.3.2.51 Intentionally blank	62
	6.3.3.2.52 Intentionally blank	62
	6.3.3.2.53 Notifications, Alerts, and Reminders Section 1.3.6.1.4.1.19376.1.5.3.1.1.20.2.1.x	62
	6.3.3.2.54 Pain Assessment Panel Section 1.3.6.1.4.1.19376.1.5.3.1.1.20.2.4	63
	6.3.3.2.55 History of Cognitive Function Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.11	63
180	6.3.3.2.56 Isolation Status Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.8	64
	6.3.3.2.57 Restraints Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.10	64
	6.3.3.2.58 Risk Indicators for Hearing Loss	65
	6.3.3.2.59 Cancer Diagnosis Section 1.3.6.1.4.1.19376.1.7.3.1.3.14.1	65
	6.3.3.3 Medications	66
185	6.3.3.3.1 Medications Section	66
	6.3.3.3.2 Admission Medication History Section	66
	6.3.3.3.3 Medications Administered Section	66
	6.3.3.3.4 Hospital Discharge Medications Section	66
	6.3.3.3.5 Immunizations Section	67
190	6.3.3.4 Physical Exams	67
	6.3.3.4.1 Physical Examination Section	67
	6.3.3.4.2 Detailed Physical Examination Section	67
	6.3.3.4.3 Hospital Discharge Physical Exam Section	67
	6.3.3.4.4 Vital Signs Section	67
195	6.3.3.4.5 Coded Vital Signs Section	67
	6.3.3.4.29 Extremities	67
	6.3.3.4.30 Coded Detailed Physical Examination Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.15.1	67
	6.3.3.4.31 Pelvis Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.10	69
	6.3.3.4.32 Admission Physical Exam Section 1.3.6.1.4.1.19376.1.5.3.1.1.22.1.1.2.1	70
200	6.3.3.4.33 Discharge Status 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.12	70
	6.3.3.5 Relevant Studies	71
	6.3.3.5.1 Results	71
	6.3.3.5.2 Coded Results	71
	6.3.3.5.3 Hospital Studies Summary	71
	6.3.3.5.4 Coded Hospital Studies Summary	71
205	6.3.3.5.5 Consultations 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.8	71
	6.3.3.5.6 Antenatal Testing and Surveillance Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5	71
	6.3.3.5.7 Coded Antenatal Testing and Surveillance Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5.1	72
	6.3.3.5.9 Intentionally blank	73
210	6.3.3.5.10 Intentionally blank	73
	6.3.3.5.11 Hearing Screening Coded Results	73
	6.3.3.5.11.1 Parent Template	74
	6.3.3.6 Plans of Care	74
	6.3.3.6.1 Care Plan	75
	6.3.3.6.2 Assessment and Plan	75
215	6.3.3.6.3 Discharge Disposition	75
	6.3.3.6.4 Discharge Diet	75
	6.3.3.6.5 Advance Directives	75
	6.3.3.6.6 Coded Advance Directives	75
	6.3.3.6.7 Transport Mode	75
220	6.3.3.6.8 Procedure Care Plan Status Report Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.45	75
	6.3.3.6.9 Health Maintenance Care Plan Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.50	76
	6.3.3.6.10 Health Maintenance Care Plan Status Report Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.41	77

	6.3.3.6.11 Provider Orders Section 1.3.6.1.4.1.19376.1.5.3.1.1.20.2.1	77
225	6.3.3.6.12 Birth Plan Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.1	78
	6.3.3.6.13 Immunization Recommendations 1.3.6.1.4.1.19376.1.5.3.1.1.18.3.1	79
	6.3.3.6.14 Patient Education Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.38	79
	6.3.3.6.15 Care Plan Section 1.3.6.1.4.1.19376.1.5.3.1.3.31	80
	6.3.3.6.16 Diet and Nutrition Section 1.3.6.1.4.1.19376.1.5.3.1.1.20.2.2	81
230	6.3.3.6.17 Intake and Output Section 1.3.6.1.4.1.19376.1.5.3.1.1.20.2.3	82
	6.3.3.6.18 Intentionally blank	82
	6.3.3.6.19 Intentionally blank	82
	6.3.3.6.20 Procedure Care Plan Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.40	82
235	6.3.3.7 Administrative and Other Information	85
	6.3.3.7.1 Payers	85
	6.3.3.7.2 Referral Source	85
	6.3.3.7.3 Transport Mode	85
	6.3.3.7.4 ED Disposition	85
	6.3.3.7.5 Intentionally blank	85
240	6.3.3.7.6 Sending Facility Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.1	85
	6.3.3.7.7 Receiving Facility Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.2	86
	6.3.3.7.8 Mass Casualty Incident Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.3	86
	6.3.3.7.9 Unit Response Level Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.4	87
	6.3.3.7.10 Extra Attendants Information Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.6	87
	6.3.3.7.11 Provider Level Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.9	88
245	6.3.3.8 Interventions	88
	6.3.3.8.3 Procedures and Interventions Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11	89
	6.3.3.8.4 Intravenous Fluids Administered Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.6	89
	6.3.3.9 Impressions	90
	6.3.3.9.1 Pre-procedure Impressions Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.42 (Deprecated)	90
250	6.3.3.9.2 Pre-procedure Risk Assessment Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.44	90
	6.3.3.9.3 Antepartum Visit Summary Flowsheet Section 1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.2	91
	6.3.3.9.4 Progress Note Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.7	92
	6.3.3.9.5 ED Diagnosis Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.9	93
	6.3.3.9.6 Acuity Assessment Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.2	93
	6.3.3.9.7 Assessments Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.4	94
255	6.3.4 CDA Entry Content Modules	95
	6.3.4.25 Family History Observation 1.3.6.1.4.1.19376.1.5.3.1.4.13.3	95
	6.3.4.25.1 Standards	95
260	6.3.4.25.2 Parent Template	95
	6.3.4.25.3 Specification	95
	6.3.4.25.4 <templateId root='2.16.840.1.113883.10.20.1.22'> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13.3'>	95
	6.3.4.25.5 <code code=' ' displayName=' ' codeSystem=' ' codeSystemName=' '/>	95
	6.3.4.25.6 <value xsi:type='CD' code=' ' displayName=' ' codeSystem=' ' codeSystemName=' '/>	96
265	6.3.4.27 EDD Observation 1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.1	96
	6.3.4.27.1 Specification	97
	6.3.4.27.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.1'>	97
	6.3.4.27.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13'>	98
	6.3.4.27.4 <code code='11778-8' codeSystem='2.16.840.1.113883.6.1'>	98
270	6.3.4.27.5 <value xsi:type='TS' value=' '>	98
	6.3.4.27.6 <author typeCode='AUT'><assignedAuthor><id root=' ' extension=' '/></assignedAuthor></author>	98
	6.3.4.27.7 <author typeCode='AUT'><time value=' '/></author>	98
	6.3.4.27.8 <entryRelationship typeCode='SPRT'>	98
275	6.3.4.27.9 <observation> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13'> : </observation> [1st nesting] 99	99
	6.3.4.27.10 <code code=' ' codeSystem='2.16.840.1.113883.6.1'> [1st nesting]	99
	6.3.4.27.11 <entryRelationship typeCode='DRIV'>	99
	6.3.4.27.12 <observation> <templateId root=' ' /> : </observation> [2st nesting]	99
	6.3.4.27.13 <code code=' ' codeSystem='2.16.840.1.113883.6.1'> [2nd nesting]	100

280	6.3.4.27.14 <repeatNumber value=''/> <interpretationCode code=' ' codeSystem=' '/> <targetSiteCode code=' ' codeSystem=' '/>	100
	6.3.4.28 Antepartum Visit Summary Battery 1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.2.....	100
	6.3.4.28.1 Specification	101
285	6.3.4.28.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.2'>.....	101
	6.3.4.28.3 <organizer classCode='BATTERY' moodCode='EVN'>.....	101
	6.3.4.28.4 <id root=' ' extension=' '/>.....	101
	6.3.4.28.6 <author><time><assignedAuthor><id>/<assignedAuthor></author>.....	102
	6.3.4.28.7 <statusCode code='completed'>.....	102
	6.3.4.28.8 <component>	102
290	6.3.4.29 Advance Directive Observation 1.3.6.1.4.1.19376.1.5.3.1.4.13.7.....	104
	6.3.4.29.1 Standards.....	104
	6.3.4.29.2 Specification	104
	6.3.4.29.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13'> <templateId root='2.16.840.1.113883.10.20.1.17'> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13.7'>.....	104
295	6.3.4.29.4 <code code=' ' codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT'>	105
	6.3.4.29.5 <value xsi:type='BL' value='true false'>	105
	6.3.4.29.6 <reference typeCode='REFR'> <templateId root='2.16.840.1.113883.10.20.1.36'> <externalDocument classCode='DOC' moodCode='EVN'> <id root=' ' extension=' '/> <text><reference value=' '/></text> 105	
300	6.3.4.30 Blood Type Observation 1.3.6.1.4.1.19376.1.5.3.1.4.13.6.....	106
	6.3.4.30.1 Standards.....	106
	6.3.4.30.2 Specification	106
	6.3.4.30.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13'> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13.6'> <templateId root='2.16.840.1.113883.10.20.1.31'>.....	106
305	6.3.4.30.4 <code code='882-1' displayName='ABO+RH GROUP' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'> 106	
	6.3.4.31 Encounters 1.3.6.1.4.1.19376.1.5.3.1.4.14	107
	6.3.4.31.1 Standards.....	107
	6.3.4.31.2 Specification	107
310	6.3.4.31.2.1 <encounter classCode='ENC' moodCode='APT ARQ EVN'>	108
	6.3.4.31.2.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.14'>	108
	6.3.4.31.2.3 <id root=' ' extension=' '/>	108
	6.3.4.31.2.4 <code code=' ' codeSystem='2.16.840.1.113883.5.4' codeSystemName='ActEncounterCode' /> ..	108
	6.3.4.31.2.5 <text><reference value="#xxx"/></text>	109
315	6.3.4.31.2.6 <effectiveTime><low value='/'><high value='/'></effectiveTime>	109
	6.3.4.31.2.7 <priorityCode code='CS'>	109
	6.3.4.31.2.8 <performer>	109
	6.3.4.31.2.9 <participant typeCode='LOC'> <participantRole classCode='SDLOC'>	109
	6.3.4.31.2.10 <id/>	109
320	6.3.4.31.2.11 <code/>	109
	6.3.4.31.2.12 <addr>...</addr>	109
	6.3.4.31.2.13 <telecom value=' ' use=' '/>	109
	6.3.4.31.2.14 <playingEntity classCode='PLC'> <name>...</name> </playingEntity>	110
325	6.3.4.32 Update Entry 1.3.6.1.4.1.19376.1.5.3.1.4.16	110
	6.3.4.32.1 Specification	110
	6.3.4.32.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.16'>	110
	6.3.4.32.3 <reference typeCode='RPLC'>	110
	6.3.4.32.4 <externalAct classCode='ACT' moodCode='EVN'>	111
	6.3.4.32.5 <id root=' ' extension=' '/>	111
330	6.3.4.33 Procedure Entry 1.3.6.1.4.1.19376.1.5.3.1.4.19	111
	6.3.4.33.1 Standards.....	111
	6.3.4.33.2 Specification	112
	6.3.4.33.2.1 <procedure classCode='PROC' moodCode='EVN INT'>	112
	6.3.4.33.2.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.19'>	112
335	6.3.4.33.2.3 <id root=' ' extension=' '/>	112
	6.3.4.33.2.4 <code code=' ' displayName=' ' codeSystem=' ' codeSystemName=' '/>	112
	6.3.4.33.2.5 <text><reference value="#xxx"/></text>	112

	6.3.4.33.2.6 <statusCode code='completed active aborted cancelled' />	113
340	6.3.4.33.2.7 <effectiveTime><low value="/" /><high value="/" /></effectiveTime>	113
	6.3.4.33.2.8 <priorityCode code="/" />	113
	6.3.4.33.2.9 <approachSiteCode code=" displayName=" codeSystem=" codeSystemName="/" />	113
	6.3.4.33.2.10 <targetSiteCode code=" displayName=" codeSystem=" codeSystemName="/" />	113
	6.3.4.33.2.11 <entryRelationship typeCode='COMP' inversionInd='true' />	113
	6.3.4.33.2.12 <entryRelationship typeCode='RSON' />	113
345	6.3.4.34 Transport 1.3.6.1.4.1.19376.1.5.3.1.1.10.4.1	113
	6.3.4.34.1 Specification	114
	6.3.4.34.1.1 <act classCode='ACT' moodCode='INT EVN' />	114
	6.3.4.34.1.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.10.4.1' />	114
350	6.3.4.34.1.3 <id root=" extension="/" />	114
	6.3.4.34.1.4 <code code=" displayName=" codeSystem='2.16.840.1.113883.6.102.4.2' codeSystemName='DEEDS4.02' />	114
	6.3.4.34.1.5 <originalText><reference value="#xxx" /><originalText>	114
	6.3.4.34.1.6 <text><reference value="#text" /></text>	114
355	6.3.4.34.1.7 <effectiveTime>	115
	6.3.4.34.1.8 <low value="/" />	115
	6.3.4.34.1.9 <high value="/" />	115
	6.3.4.35 Encounter Disposition 1.3.6.1.4.1.19376.1.5.3.1.1.10.4.2	115
	6.3.4.35.1 Specification	116
360	6.3.4.35.1.1 <act classCode='ACT' moodCode='INT EVN' />	116
	6.3.4.35.1.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.10.4.2' />	116
	6.3.4.35.1.3 <id root=" extension="/" />	116
	6.3.4.35.1.4 <code code=" displayName=" codeSystem=" codeSystemName=" />	117
	6.3.4.35.1.5 <text><reference value="#xxx" /></text>	117
	6.3.4.35.1.6 <effectiveTime><low value="/" /><high value="/" /><effectiveTime>	117
365	6.3.4.35.1.7 <performer typeCode='PRF' />	117
	6.3.4.35.1.8 <assignedEntity />	117
	6.3.4.35.1.9 <id root=" extension="/" />	117
	6.3.4.35.1.10 <addr></addr>	117
370	6.3.4.35.1.11 <telecom value=" use="/" />	117
	6.3.4.35.1.12 <assignedPerson><name/" /></assignedPerson>	118
	6.3.4.35.1.13 <participant typeCode='RCV' ><time value="/" /> <participantRole classCode='ROL' ><id root=" extension="/" /> <addr></addr> <telecom value=" use="/" /> <playingEntity><name/" /></playingEntity>	118
	6.3.4.36 Coverage Entry 1.3.6.1.4.1.19376.1.5.3.1.4.17	118
375	6.3.4.36.1 Standards	118
	6.3.4.36.2 Specification	118
	6.3.4.36.3 <act classCode='ACT' moodCode='DEF' />	119
	6.3.4.36.4 <templateId root='2.16.840.1.113883.10.20.1.20' /> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.17' /> 119	119
380	6.3.4.36.5 <id root=" extension="/" />	119
	6.3.4.36.6 <code code='48768-6' displayName='PAYMENT SOURCES' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />	119
	6.3.4.36.7 <statusCode code='completed' />	119
	6.3.4.36.8 <entryRelationship typeCode='COMP' />	119
385	6.3.4.36.9 <sequenceNumber value=' />	119
	6.3.4.37 Payer Entry 1.3.6.1.4.1.19376.1.5.3.1.4.18	119
	6.3.4.37.1 Standards	120
	6.3.4.37.2 Specification	120
	6.3.4.37.3 <act classCode='ACT' moodCode='EVN' />	120
390	6.3.4.37.4 <templateId root='2.16.840.1.113883.10.20.1.26' /> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.18' /> 121	121
	6.3.4.37.5 <id root=" extension="/" />	121
	6.3.4.37.6 <code code=' ' displayName=' ' codeSystem=' ' codeSystemName=' ' />	121
	6.3.4.37.7 <statusCode code='completed' />	121
395	6.3.4.37.8 <performer typeCode='PRF' > <assignedEntity classCode='ASSIGNED' />	121
	6.3.4.37.9 <id root=" extension="/" />	121

	6.3.4.37.10 <code code='PAYOR GUAR PAT' displayName=' ' codeSystem='2.16.840.1.113883.5.110' codeSystemName='RoleClass'/>	122
400	6.3.4.37.11 <addr></addr>	122
	6.3.4.37.12 <telecom value=' ' use=' '/>	122
	6.3.4.37.13 <representedOrganization typeCode='ORG'> <name></name>	122
	6.3.4.37.14 <participant typeCode='COV'> <participantRole classCode='PAT'>	122
	6.3.4.37.15 <id root=' ' extension=' '/>	122
	6.3.4.37.16 <code code= displayName= codeSystem=2.16.840.1.113883.5.111' codeSystemName='RoleCode' />	122
405	6.3.4.37.17 <addr></addr>	123
	6.3.4.37.18 <telecom value=' ' use=' '/>	123
	6.3.4.37.19 <playingEntity><name></name></playingEntity>	123
	6.3.4.37.20 <participant typeCode='HLD'> <participantRole classCode='IND'>	123
410	6.3.4.37.21 <id root=' ' extension=' '/>	123
	6.3.4.37.22 <addr></addr>	123
	6.3.4.37.23 <telecom value=' ' use=' '/>	124
	6.3.4.37.24 <playingEntity><name></name></playingEntity>	124
	6.3.4.37.25 <entryRelationship typeCode='REFR'> <act classCode='ACT' moodCode='DEF'>	124
415	6.3.4.37.26 <id root=' ' extension=' '/>	124
	6.3.4.37.27 <text><reference value=' '/></text>	124
	6.3.4.38 Pain Score Observation 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.1	124
	6.3.4.38.1 Parent Template	124
	6.3.4.38.2 Specification	125
420	6.3.4.38.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />	125
	6.3.4.38.4 <code code='38208 5' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'> <translation code='406127006' displayName='Pain intensity' codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT' />	125
425	6.3.4.38.5 <value xsi:type='CO' value=' ' />	125
	6.3.4.38.6 <interpretationCode code='301379001 40196000 76948002 67849003' codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT' />	126
	6.3.4.38.7 <methodCode code=' ' codeSystem=' ' codeSystemName=' '/>	126
	6.3.4.38.8 <targetSiteCode code=' ' codeSystem=' ' codeSystemName=' '/>	126
	6.3.4.39 Braden Score Observation 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.2	126
430	6.3.4.40 Braden Score Component 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.3	126
	6.3.4.41 Geriatric Depression Score Observation 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.4	126
	6.3.4.42 Geriatric Depression Score Component 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.5	127
	6.3.4.43 Survey Panel 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.7	127
	6.3.4.43.1 Parent Template	127
435	6.3.4.43.2 Uses	127
	6.3.4.43.3 Specification	127
	6.3.4.43.3.1 <organizer classCode='CLUSTER' moodCode='EVN'>	127
	6.3.4.43.3.2 <templateId root='2.16.840.1.113883.10.20.1.32' /> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.7' />	127
440	6.3.4.43.3.3 <id root=' ' extension=' '/>	128
	6.3.4.43.3.4 <code code=' ' displayName=' ' codeSystem=' ' codeSystemName=' '/>	128
	6.3.4.43.3.5 <statusCode code='completed' />	128
	6.3.4.43.3.6 <effectiveTime value=' '/>	128
	6.3.4.43.3.7 <!-- one or more survey observations --> <component typeCode='COMP'>	128
445	6.3.4.44 Survey Observation 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.6	128
	6.3.4.44.1 Parent Template	128
	6.3.4.44.2 Uses	128
450	6.3.4.44.3 Specification	128
	6.3.4.44.3.1 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' /> <templateId root='2.16.840.1.113883.10.20.1.31' /> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.6' />	129
	6.3.4.44.3.2 <code code=' ' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />	129
	6.3.4.44.3.3 <value xsi:type='CO CD INT PQ' ... />	129
	6.3.4.44.3.4 <interpretationCode code=' ' codeSystem=' ' codeSystemName=' '/>	129

	6.3.4.44.3.5 <methodCode code=' ' codeSystem=' ' codeSystemName=' '><targetSiteCode code=' ' codeSystem=' ' codeSystemName=' '>.....	129
455	6.3.4.45 Acuity 1.3.6.1.4.1.19376.1.5.3.1.1.13.3.1	129
	6.3.4.45.1 Specification	130
	6.3.4.45.1.1 <observation classCode='OBS' moodCode='EVN'>.....	130
	6.3.4.45.1.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.1'>	130
	6.3.4.45.1.3 <id root=" extension="/">	130
460	6.3.4.45.1.4 <code code=" displayName=" codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT'>.....	130
	6.3.4.45.1.5 <originalText><reference value="#xxx"/><orginalText>.....	130
	6.3.4.45.1.6 <text><reference value="#text"/></text>.....	130
	6.3.4.45.1.7 <effectiveTime>	130
465	6.3.4.45.1.8 <high value="/".....	131
	6.3.4.46 Intravenous Fluids 1.3.6.1.4.1.19376.1.5.3.1.1.13.3.2.....	131
	6.3.4.46.1 Specification	132
	6.3.4.46.1.1 Medication Fields	132
470	6.3.4.46.1.2 <substanceAdministration classCode='SBADM' moodCode='INT EVN'>	133
	6.3.4.46.1.3 <templateId root='2.16.840.1.113883.10.20.1.24'><templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.7.1' />	133
	6.3.4.46.1.4 <id root=" extension="/">	133
	6.3.4.46.1.5 <code code=" displayName=" codeSystem=" codeSystemName=">	133
	6.3.4.46.1.6 <text><reference value="/"></text>	134
475	6.3.4.46.1.7 <statusCode code='completed active'/>	134
	6.3.4.46.1.8 <effectiveTime xsi:type='IVL_TS'>	134
	6.3.4.46.1.9 <low value="/"><high value="/".....	134
	6.3.4.46.1.10 <approachSiteCode code=" codeSystem="> originalText><reference value="/"></originalText></approachSiteCode>	134
480	6.3.4.46.1.11 <doseQuantity><low value=" unit="/"><high value=" unit="/" ></doseQuantity>	134
	6.3.4.46.1.12 <low high value="/" > <translation> <originalText><reference value="/"></originalText></translation></low high >	135
	6.3.4.46.1.13 <rateQuantity><low value=" unit="/"><high value=" unit="/"></rateQuantity>	135
	6.3.4.46.1.14 <consumable>	135
485	6.3.4.47 Nursing Assessments Battery 1.3.6.1.4.1.19376.1.5.3.1.1.13.3.4.....	135
	6.3.4.47.1 Specification	136
	6.3.4.47.1.1 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.4'>	136
	6.3.4.47.1.2 <organizer classCode='BATTERY' moodCode='EVN'>	136
490	6.3.4.47.1.3 <id root=" extension="/">	136
	6.3.4.47.1.4 <code code='X-ASSESS' codeSystem='2.16.840.1.113883.6.1'>	136
	6.3.4.47.1.5 <author/><time/><assignedAuthor><id/></assignedAuthor></author>	137
	6.3.4.47.1.6 <statusCode code='completed' />	137
	6.3.4.47.1.7 <component>	137
	6.3.4.48 Antenatal Testing and Surveillance Battery 1.3.6.1.4.1.19376.1.5.3.1.1.21.3.10.....	137
495	6.3.4.48.1 Specification	138
	6.3.4.48.1.1 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.3.10'>	138
	6.3.4.48.1.2 <organizer classCode='BATTERY' moodCode='EVN'>	138
	6.3.4.48.1.3 <id root=" extension="/">	138
500	6.3.4.48.1.4 <code code='XX- XX-ANTENATALTESTINGBATTERY' codeSystem='2.16.840.1.113883.6.1'>	138
	6.3.4.48.1.5 <author/><time/><assignedAuthor><id/></assignedAuthor></author>	139
	6.3.4.48.1.6 <statusCode code='completed' />	139
	6.3.4.48.1.7 <component>	139
505	6.3.4.49 Immunization Recommendation 1.3.6.1.4.1.19376.1.5.3.1.4.12.2	139
	6.3.4.49.1 <substanceAdministrationIntent typeCode='SBADM' moodCode='INT PRP' negationInd='true false'>140	140
	6.3.4.49.2 <templateId root='2.16.840.1.113883.10.20.1.25'>	140
	6.3.4.49.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.12.2'>	140
	6.3.4.49.4 <id root=" extension="/">	140
510	6.3.4.49.5 <code code='IMMUNIZ' codeSystem='2.16.840.1.113883.5.4' codeSystemName='ActCode'>	141
	6.3.4.49.6 <text><reference value='#xxx'></text>	141

	6.3.4.49.7 <statusCode code='active'>.....	141
	6.3.4.49.8 <effectiveTime><low value=''/><high value=''/></effectiveTime>.....	141
	6.3.4.49.9 <routeCode code="" codeSystem="" codeSystemName='RouteOfAdministration'> <approachSiteCode code="" codeSystem="" codeSystemName='HumanSubstanceAdministrationSite'> <doseQuantity value="" units="">.....	141
515	141	
	6.3.4.49.10 <consumable typeCode='CSM'>.....	141
	6.3.4.49.11 <entryRelationship typeCode='SUBJ'> <observation classCode='OBS' moodCode='EVN'> <templateId root='2.16.840.1.113883.10.20.1.46'> <code code='30973-2' displayName='Dose Number' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'> <statusCode code='completed'> <value xsi:type='INT' value=''/> </observation> </entryRelationship>	141
520	141	
	6.3.4.49.12 <!-- Optional <entryRelationship> element referencing guidelines -->	142
	6.3.4.50 Alert Entry 1.3.6.1.4.1.19376.1.5.3.1.4.12.3	142
	6.3.4.50.1 <entryRelationship typeCode='SUBJ' inversionInd='true'>.....	142
	6.3.4.50.2 <templateId root='TBD'>.....	142
525	142	
	6.3.4.50.3 <observation classCode='OBS ALRT' moodCode='EVN'>.....	142
	6.3.4.50.4 <id root="" extension="">.....	143
	6.3.4.50.5 <code code="" displayName="" codeSystem='2.16.840.1.113883.5.4' codeSystemName='ActDetectedIssueCode'>	143
	6.3.4.50.6 <text><reference value="#ref-1"/></text>	143
530	143	
	6.3.4.51 Antigen Dose 1.3.6.1.4.1.19376.1.5.3.1.4.12.1	143
	6.3.4.51.1 Specification	143
	6.3.4.51.2 <substanceAdministration typeCode='SBADM' moodCode='EVN' negationInd='false'>.....	144
	6.3.4.51.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.12.1'>	144
	6.3.4.51.4 <id root="" extension="">	144
535	144	
	6.3.4.51.5 <code/><text/><statusCode/><effectiveTime value=''/><routeCode code=''' codeSystem=''' codeSystemName='RouteOfAdministration'> <approachSiteCode code=''' codeSystem=''' codeSystemName='HumanSubstanceAdministrationSite'>	144
	6.3.4.51.6 <doseQuantity value=' ' units=' '/>	144
540	144	
	6.3.4.51.7 <consumable typeCode='CSM'>	144
	6.3.4.52 Intentionally blank.....	144
	6.3.4.53 Intentionally blank.....	144
	6.3.4.54 Observation Request 1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1	145
	6.3.4.54.1 Uses	145
545	145	
	6.3.4.54.2 Specification	145
	6.3.4.54.2.1 <observation classCode='OBS' moodCode='INT PRP GOL'>	145
	6.3.4.54.2.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1'>	146
	6.3.4.54.2.3 <templateId root='2.16.840.1.113883.10.20.1.25'>	146
	6.3.4.54.2.4 <id root="" extension="">	146
550	146	
	6.3.4.54.2.5 <code code=''' displayName=''' codeSystem=''' codeSystemName='''>	146
	6.3.4.54.2.6 <text><reference value="#xxx"/></text> -OR- <text>text</text>	146
	6.3.4.54.2.7 <statusCode code='active'>	146
	6.3.4.54.2.8 <effectiveTime value=''/>	146
	6.3.4.54.2.9 <value xsi:type=' ' .../>	146
555	146	
	6.3.4.54.2.10 <methodCode code=''' codeSystem=''' codeSystemName='''>	146
	6.3.4.54.2.11 <targetSiteCode code=''' codeSystem=''' codeSystemName='''>	147
	6.3.4.54.2.12 <author><assignedAuthor classCode='ASSIGNED'>...<assignedAuthor></author>	147
	6.3.4.55 Risk Indicators for Hearing Loss Entry 1.3.6.1.4.1.19376.1.7.3.1.1.15.5.1	147
	6.3.4.55.1 Specification	147
560	147	
	6.3.4.55.2 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.15.5.1'>	148
	6.3.4.55.3 <organizer classCode='BATTERY' moodCode='EVN'>	148
	6.3.4.55.4 <id root="" extension="">	148
	6.3.4.55.5 <code code='58232-0' codeSystem='2.16.840.1.113883.6.1'>	149
	6.3.4.55.6 <author><time/><assignedAuthor><id/></assignedAuthor></author>	149
565	149	
	6.3.4.55.7 <statusCode code='completed'>	149
	6.4.4.55.8 <component>	149
	6.3.4.56 Cancer Diagnosis Entry 1.3.6.1.4.1.19376.1.7.3.1.4.14.1	149
	6.3.4.56.1 Parent Template	149
	6.3.4.56.2 Specification	150

570	6.3.4.56.3 <act classCode='ACT' moodCode='EVN'>.....	152
	6.3.4.56.4 <templateId root='2.16.840.1.113883.10.20.1.27'> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5.2'>.....	152
	6.3.4.56.5 <!-- 1..* entry relationships identifying problems of concern --><entryRelationship type='SUBJ'><observation classCode='OBS' moodCode='EVN'><templateIDroot='1.3.6.1.4.1.19376.1.5.3.1.4.5.'>...</observation>	152
575	6.3.4.56.6 <observation classCode="OBS" moodCode="EVN">	152
	6.3.4.56.7 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5.'> <templateId root='1.3.6.1.4.1.19376.1.7.3.1.4.14.1'>.....	152
	6.3.4.56.8 <code code="282291009" codeSystem=" 2.16.840.1.113883.6.96" codeSystemName="SNOMED CT" displayName="Diagnosis"/>.....	153
580	6.3.4.56.9 <statusCode code="completed"/>.....	153
	6.3.4.56.10 <effectiveTime value="xxx"/>.....	153
	6.3.4.56.11 <value xsi:type='CD' code=' ' codeSystem=' ' codeSystemName=' ' displayName=' '>	153
	6.3.4.56.12 <qualifier><name code="31206-6" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" displayName=" Behavior ICD-O-3 Cancer"/><value code="" codeSystem="" codeSystemName=" " displayName=" "/></qualifier>	153
585	6.3.4.56.13 <qualifier><name code="21861-0" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" displayName="Dx confirmed by Cancer"/><value xsi:type="CD" code="" codeSystem="" codeSystemName=" " displayName=" "/></qualifier>	153
	6.3.4.56.14 <targetSiteCode code=" " codeSystem=" " codeSystemName=" " displayName=" ">	154
590	6.3.4.56.15 <qualifier><name code="20228-3" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" displayName="Anatomic part Laterality"/> <value code="" codeSystem="" codeSystemName=" " displayName=" "/></qualifier>	154
	6.3.4.56.16 <entryRelationship typeCode="SUBJ" inversionInd="false">	154
595	6.3.4.56.17 <observation classCode="OBS" moodCode="EVN"> <templateId root="1.3.6.1.4.1.19376.1.7.3.1.4.14.2"/> [1 st nesting]	154
	6.3.4.56.18 <code code="xxxxx-x" displayName="TNM Clinical Stage Information" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/> [1 st nesting]	154
	6.3.4.56.19 <statusCode code="completed"/> [1 st nesting]	154
600	6.3.4.56.20 <value xsi:type="CD" code="" codeSystem="" codeSystemName="" displayName=" "> [1 st nesting]	154
	6.3.4.56.21 <qualifier><name code="21909-7" displayName=" Descriptor.clinical Cancer" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/> <value xsi:type="CD" code="" codeSystem="" codeSystemName=" " displayName=" "/></qualifier> [1 st nesting]	155
605	6.3.4.56.22 <qualifier><name code="21917-0" displayName="Version TNM Classification" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/><value xsi:type="CD" code="" codeSystem="" codeSystemName=" " displayName=" "/></qualifier> [1 st nesting]	155
	6.3.4.56.23 <participant typeCode="PPRF"> <participantRole> <code code="21910-5" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" displayName="Stager.clinical Cancer"/><playingEntity nullFlavor="NA"> <code xsi:type="CE" code="" codeSystem="" codeSystemName=" " displayName=" "/> [1 st nesting]	155
610	6.3.4.56.24 <!-- 0..3 entryRelationships identifying simple observations for TNM Clinic Tumor, TNM Clinical Nodes, and TNM Clinical Metastases--><entryRelationship typeCode="COMP" inversionInd="false"><observation classCode='OBS'moodCode='EVN'><templateIDroot='1.3.6.1.4.1.19376.1.5.3.1.4.13'>...</observation>[2nd nesting]	155
	6.3.4.56.25 <code code="" displayName=" " codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/> [2 nd nesting]	156
	6.3.4.56.26 <value xsi:type="CD" code="" codeSystem="" codeSystemName=" " displayName=" "/>	156
	6.3.4.56.1 Parent Template	156
	6.3.4.56.2 Specification	157
620	6.4 HL7 VERSION 2.0 CONTENT MODULES	159
	6.5 PCC VALUE SETS	159
	<i>6.5.A Antepartum History of Past Illness Value Set 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.1</i>	159
	<i>6.5.C Antepartum Family History and Genetic Screening Value Set 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.4</i>	161
	<i>6.5.D Antepartum Review of Systems Menstrual History Value Set 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.5</i>	162
625	<i>6.5.E Antepartum History of Infection Value Set 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.6</i>	162
	<i>6.5.F Antepartum Laboratory Value Set 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.7</i>	163

	<i>6.5.G Antepartum Education Value Set 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.8.....</i>	166
	<i>6.5.H JCIH-EHDI Risk Indicators for Hearing Loss (LOINC®) Value Set</i>	168
	6.5.H.1 Metadata.....	168
630	6.5.H.2 JCIH-EHDI Risk Indicators for Hearing Loss (LOINC®) Value Set Value Set	169
	<i>6.5.I JCIH-EHDI Risk Indicators for Hearing Loss Codes</i>	170
	6.5.I.1 Metadata	170
	6.5.I.2 JCIH-EHDI Risk Indicators for Hearing Loss Value Set.....	171
	6.5.I.3 Pending Codes for SNOMED-CT Findings/Situation to support Risk Indicators for Hearing Loss.....	171
635	<i>6.5.J JCIH-EHDI Risk Indicators for Hearing Loss - Procedures Codes</i>	171
	6.5.J.1 Metadata.....	171
	6.5.J.2 JCIH-EHDI Risk Indicators for Hearing Loss - Procedures Value.....	172
	<i>6.5.K Newborn Hearing Procedure Codes</i>	173
640	6.5.K.1 Metadata.....	173
	6.5.K.2 JCIH-EHDI Newborn Hearing Procedure Value Set.....	174
	<i>6.5.L JCIH-EHDI Newborn Hearing Screening Method Codes.....</i>	174
	6.5.L.1 Metadata	174
	6.5.L.2 JCIH-EHDI Newborn Hearing Screening Method Value Set.....	175
645	<i>6.5.M JCIH-EHDI Hearing Screen Right Codes– Right.....</i>	175
	6.5.M.1 Metadata	175
	6.5.M.2 JCIH-EHDI Hearing Screen Right Value Set.....	176
	<i>6.5.N JCIH-EHDI Hearing Screen Left Codes</i>	176
	6.5.N.1 Metadata	176
	6.5.N.2 JCIH-EHDI Hearing Screen Left Value Set.....	177
650	<i>6.5.O JCIH-EHDI Reason for no Hearing Loss Diagnosis or Screening Codes(SNOMED)</i>	177
	6.5.O.1 Metadata	177
	6.5.O.2 JCIH-EHDI Reason for no Hearing Loss Diagnosis or Screening Value Set.....	178
	<i>6.5.P JCIH-EHDI Newborn Hearing Loss Referrals Codes.....</i>	179
	6.5.P.1 Metadata.....	179
655	6.5.P.2 JCIH-EHDI Newborn Hearing Loss Referrals Value Set	180
	<i>6.5.Q JCIH-EHDI Newborn Hearing Loss Reason for no Follow-up – Patient Reason Codes</i>	180
	6.5.Q.1 Metadata	180
	6.5.Q.2 JCIH-EHDI Newborn Hearing Loss Reason for no Follow-up – Patient Reason Value Set	181
	<i>6.5.R Joint Commission Medical Reason Codes.....</i>	183
660	6.5.R.1 Metadata	183
	6.5.R.2 Joint Commission Medical Reason Value Set	183
	<i>6.5.S JCIH-EHDI Inpatient Screening Results not Performed Codes</i>	184
	6.5.S.1 Metadata	184
	6.5.S.2 JCIH-EHDI Inpatient Screening Results not Performed Value Set	184
665	<i>6.5.T JCIH-EHDI Evidence of Hearing Screening Performed Codes</i>	185
	6.5.T.1 Metadata	185
	6.5.T.2 JCIH-EHDI Evidence of Hearing Screening Performed Value Set	186
	<i>6.5.U JCIH-EHDI Procedure Declined Value Set Codes</i>	186
	6.5.U.1 Metadata	186
670	6.5.U.2 JCIH-EHDI JCIH-EHDI Procedure Declined Value Set Value Set	187
	<i>6.5.V JCIH-EHDI Newborn Hearing Screening Abnormal Results Value Set Codes.....</i>	187
	6.5.V.1 Metadata	187
	6.5.V.2 JCIH-EHDI Newborn Hearing Screening Abnormal Results Value Set	188
	<i>6.5.W Primary Site Value Set.....</i>	188
675	<i>6.5.X Histologic Type Value Set.....</i>	189
	<i>6.5.Y Derived AJCC Descriptor (T,N,M) Value Set.....</i>	189
	<i>6.5.Z TNM Edition Value Set</i>	189
	<i>6.5.AA TNM Stage Group Value Set</i>	190
	<i>6.5.BB TNM Stage Descriptor Value Set</i>	191
680	<i>6.5.CC TNM Tumor Value Set</i>	191
	<i>6.5.DD TNM Node Value Set.....</i>	192
	<i>6.5.DD TNM Metastasis Value Set</i>	193

APPENDIX Q: DOCUMENT CONSTRUCTION	194
--	------------

Introduction

690 This supplement is written for Trial Implementation. It is written as changes to the latest revision of the documents listed below. The reader should have already read and understood these documents:

1. [PCC Technical Framework Volume 1](#)
2. [PCC Technical Framework Volume 2](#)

This supplement also references other documents¹. The reader should have already read and understood these documents:

- 695
1. [IT Infrastructure Technical Framework Volume 1](#)
 2. [IT Infrastructure Technical Framework Volume 2](#)
 3. [IT Infrastructure Technical Framework Volume 3](#)
 4. HL7 and other standards documents referenced in Volume 1 and Volume 2

700 This supplement defines a number of PCC content modules that are shared between various content documents. These are provided for trial implementation and will be published in the same format for Trial Implementation. Upon completion, some content modules will be moved to Final Text; others may remain in Trial Implementation.

Profile Abstract

This supplement does not describe a profile

705

Open Issues and Questions

Closed Issues

¹ The first three documents can be located on the IHE Website at http://www.ihe.net/Technical_Framework/index.cfm#IT. The remaining document can be obtained from its respective publisher.

Volume 1 – Integration Profiles

710 **Glossary**

<i>Add the following terms to the Glossary:</i>

2.5 History of Annual Changes

715

<i>Add the following bullet to the end of the bullet list in Section 2.5</i>
--

- Added a set of CDA Content Modules shared across several Integration Profiles for the 2010-2011 documentation cycle.
- In the 2011-2012 documentation cycle, the following CDA Section Content Modules were added as well as various Entry Content Modules and Value Sets:

720

- PCC Transport Summary Profiles supplement

- Sending Facility
- Receiving Facility
- Mass Causality Incident
- Unit Response Level

725

- Protocols Used
- Extra Attendants Information
- Invasive Airway

- Isolation Status
- Restraints

730

- Ventilator Usage
- Provider Level

- QRPH EHCP Profile
 - Risk Indicators for Hearing Loss
 - Hearing Screening Coded Results

735

- QRPH PRPH-Ca Profile
 - Cancer Diagnosis

- In the 2012-2013 documentation cycle, edits were made based on CPs. No new content modules were added and none were removed.

740

Volume 2 – Transactions and Content Modules

Add Section 6.1

745 6.1 Conventions

Various tables used in this section will further constrain the content. Within this volume, the follow conventions are used.

R

750 A "Required" data element is one that shall always be provided. If there is information available, the data element must be present. If there is no information available, or it cannot be transmitted, the data element must contain a value indicating the reason for omission of the data. (See PCC TF-2: 5.3.4.2 for a list of appropriate statements).

R2

755 A "Required if data present" data element is one that shall be provided when a value exists. If the information cannot be transmitted, the data element shall contain a value indicating the reason for omission of the data. If no such information is available to the creator or if such information is not available in a well identified manner (e.g., buried in a free form narrative that contains additional information relevant to other sections) or if the creator requires that information be absent, the R2 section shall be entirely absent. (See Section PCC TF-2: 5.3.4.2 for a list of appropriate statements).

O

An optional data element is one that may be provided, irrespective of whether the information is available or not. If the implementation elects to support this optional section, then its support shall meet the requirement set forth for the "Required if data present" or R2.

760 C

A conditional data element is one that is required, required if known, or optional depending upon other conditions. These will have further notes explaining when the data element is required, et cetera.

770 Note: The definitions of R, R2, and O differ slightly from other IHE profiles. This is due in part to the fact that local regulations and policies may in fact prohibit the transmission of certain information, and that a human decision to transmit the information may be required in many cases.

Add Section 6.2

6.2 Folder Content Modules

775 This section contains modules that describe the content requirements of Folders used with XDS, XDM or XDR. When workflows are completed normally, the folders will contain documents with the optionality specified in the tables shown below. Under certain circumstances, the folders will not meet the optionality requirements described below, for example, when the patient leaves before treatment is completed.

780 **6.2.1 EDES Folder Specification**

This section intentionally left blank.

6.2.2 APR Folder Specification

This section intentionally left blank.

6.2.3 LDR Folder Specification

785 This section intentionally left blank.

6.3 HL7 Version 3.0 Content Modules

This section contains content modules based upon the HL7 CDA Release 3.0 Standard, and related standards and/or implementation guides.

6.3.1 CDA Document Content Modules

790

Add Section 6.3.1.X

6.3.1.X History and Physical Specification 1.3.6.1.4.1.19376.1.5.3.1.1.16.1.4

795 The History and Physical document content module is a Medical Summary and inherits all header constraints from Medical Summary (1.3.6.1.4.1.19376.1.5.3.1.1.2). The intention of this document content module is to provide a base from which other document content modules may be derived. Future work may also result in a content profile for History and Physical.

6.3.1.x.1 Format Code

The XDSDocumentEntry format code for this content is **urn:ihe:pcc:hp:2008**

6.3.1.x.2 LOINC Code

800 The LOINC code for this document is **34117-2 HISTORY AND PHYSICAL**

6.3.1.x.3 Standards

CDAR2 [HL7 CDA Release 2.0](#)

CDTHP [CDA for Common Document Types History and Physical Notes \(DSTU\)](#)

6.3.1.x.4 Specification

This section references content modules using Template ID as the key identifier. Definitions of the modules are found in either:

- 805 • IHE Patient Care Coordination Volume 2: Final Text
 • IHE PCC CDA Content Modules Supplement (this document, for Trial Implementation)

Table 6.3.1.x.4-1: History and Physical Data Elements

Data Element Name	Opt	Template ID
Chief Complaint	R	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.1
History of Present Illness	R	1.3.6.1.4.1.19376.1.5.3.1.3.4
History of Past Illness	R	1.3.6.1.4.1.19376.1.5.3.1.3.8
Medications	R	1.3.6.1.4.1.19376.1.5.3.1.3.19
Allergies and Other Adverse Reactions Section	R	1.3.6.1.4.1.19376.1.5.3.1.3.13
Social History	R	1.3.6.1.4.1.19376.1.5.3.1.3.16
Family History	R	1.3.6.1.4.1.19376.1.5.3.1.3.14
Review of Systems	R	1.3.6.1.4.1.19376.1.5.3.1.3.18
Detailed Physical Examination This section SHALL include Vital Signs (1.3.6.1.4.1.19376.1.5.3.1.3.25) as a subsection.	R	1.3.6.1.4.1.19376.1.5.3.1.1.9.15
Results Diagnostic Findings; use this OR Coded Results	R	1.3.6.1.4.1.19376.1.5.3.1.3.27
Coded Results Diagnostic Findings; use this OR Results	R	1.3.6.1.4.1.19376.1.5.3.1.3.28
Assessment and Plan	R	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.5

810 **6.3.1.x.5 Conformance**

CDA Release 2.0 documents that conform to the requirements of this content module shall indicate their conformance by the inclusion of the appropriate <templateId> elements in the header of the document. This is shown in the sample document below. A CDA Document may conform to more than one template. This content module inherits from the [Medical Summaries](#) content module, and so must conform to the requirements of that template as well, thus all <templateId> elements shown in the example below shall be included.

```

<ClinicalDocument xmlns='urn:hl7-org:v3'>
  <typeId extension="POCD_HD000040" root="2.16.840.1.113883.1.3"/>
  <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.2'/>
  <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.1.4'/>
  <id root=' ' extension=' ' />
  <code code='34117-2' displayName='HISTORY AND PHYSICAL'
    codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
  <title>History and Physical</title>
  <effectiveTime value='20080601012005' />
  <confidentialityCode code='N' displayName='Normal'
    codeSystem='2.16.840.1.113883.5.25' codeSystemName='Confidentiality' />
  <languageCode code='en-US' />
  :
  <component><structuredBody>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.1' />
        <!-- Required Chief Complaint Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.4' />
        <!-- Required History of Present Illness Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.8' />
        <!-- Required History of Past Illness Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.19' />
        <!-- Required Medications Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.13' />
        <!-- Required Allergies and Other Adverse Reactions Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.16' />
        <!-- Required Social History Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.14' />
        <!-- Required Family History Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.18' />
        <!-- Required Review of Systems Section content -->
      </section>
    </component>
    <component>
      <section>
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.15' />
        <!-- Required Detailed Physical Examination Section content -->
      </section>
    </component>
  </structuredBody>
</component>

```

```
885    </component>
886    <component>
887        <section>
888            <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.27' />
889            <!-- Required Results Section content -->
890        </section>
891    </component>
892    <component>
893        <section>
894            <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.28' />
895            <!-- Required Coded Results Section content -->
896        </section>
897    </component>
898    <component>
899        <section>
900            <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.5' />
901            <!-- Required Assessment and Plan Section content -->
902        </section>
903    </component>
904    </structuredBody></component>
905 </ClinicalDocument>
```

Figure 6.3.1.x.5-1: Sample History and Physical Document

Add Section 6.3.2

6.3.2 CDA Header Content Modules

Add Section 6.3.2.1

910 6.3.2.1 Language Communication 1.3.6.1.4.1.19376.1.5.3.1.2.1

Add Section 6.3.2.2

6.3.2.2 Employer and School Contacts 1.3.6.1.4.1.19376.1.5.3.1.2.2

Add Section 6.3.2.3

6.3.2.3 Healthcare Providers and Pharmacies 1.3.6.1.4.1.19376.1.5.3.1.2.3

915 Add Section 6.3.2.4

6.3.2.4 Patient Contacts 1.3.6.1.4.1.19376.1.5.3.1.2.4

920 Add Section 6.3.2.5

6.3.2.5 Spouse 1.3.6.1.4.1.19376.1.5.3.1.2.4.1

The spouse header element records the spouse of a patient, and inherits other constraints from the [Patient Contacts](#) entry. Items in bold in the example below show the additional constraints on this element.

- 925 This element SHALL be included as a participant in the header of the CDA document in the event of the pregnancy. If this does not apply to the patient this element SHALL use a null flavor.

6.3.2.5.1 Parent Template

The parent of this template is [Patient Contacts](#).

930 **6.3.2.5.2 Specification**

```
<participant typeCode='IND'>
  <templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4' />
  <templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4.1' />
  <time value='20070213' />
  <associatedEntity classCode='PRS'>
    <code code='xx-spouse|184142008' displayName=' ' codeSystem='2.16.840.1.113883.6.96'
codeSystemName='SNOMED CT' />
    <addr></addr>
    <telecom value=' ' use=' ' />
    <assignedPerson><name></name></assignedPerson>
  </associatedEntity>
</participant>
```

6.3.2.5.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4.1'><templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4.1'>

- 945 The <templateId> element identifies this person as a spouse and must be recorded exactly as shown above.

```
<rule context='h17:templateId/@root="1.3.6.1.4.1.19376.1.5.3.1.2.4.1"]'>
  <assert test='h17:templateId/@root="1.3.6.1.4.1.19376.1.5.3.1.2.4"'>
    A participant using template 1.3.6.1.4.1.19376.1.5.3.1.2.4.1 must also use template
    1.3.6.1.4.1.19376.1.5.3.1.2.4.
  </assert>
</rule>
```

6.3.2.5.4 <associatedEntity classCode='PRS'>

The classCode attribute of the <associatedEntity> element shall be PRS.

955

```
<rule context='h17:templateId/@root="1.3.6.1.4.1.19376.1.5.3.1.2.4.1"]'>
  <assert test='..h17:associatedEntity/@classCode = "PRS"'>
    The classCode attribute of the associated entity shall be PRS.
  </assert>
</rule>
```

960 **6.3.2.5.5 <code code='127848009|184142008' displayName=' '
codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT'>**

This element SHALL use 127848009 to represent the patient's spouse or 184142008 to represent the patient's next of kin. The code system name is SNOMED CT.

6.3.2.5.6 Completed Example

```

965 <!-- Husband/Domestic Partner -->
<participant typeCode="IND">
  <associatedEntity classCode="NOK">
    <code code="184142008" displayName="patient's next of kin"
      codeSystem="2.16.840.1.113883.6.96" codeSystemName="SNOMED CT"/>
970   <addr>
    <streetAddressLine>45 Chunn Dr.</streetAddressLine>
    <city>Spring Hill</city>
    <state>TN</state>
    <postalCode>37174</postalCode>
    <country>USA</country>
975   </addr>
   <telecom value="tel: (999) 555-1212" use="WP"/>
   <associatedPerson>
     <name>
       <prefix>Mr.</prefix>
       <given>John</given>
       <family>Youngston</family>
     </name>
980   </associatedPerson>
   </associatedEntity>
985 </participant>

```

Add Section 6.3.2.6

6.3.2.6 Natural Father of Fetus 1.3.6.1.4.1.19376.1.5.3.1.2.4.2

990 This header element records the natural father of the fetus, and inherits other constraints from the Patient Contacts (1.3.6.1.4.1.19376.1.5.3.1.2.4) entry. Items in bold in the example below show the additional constraints on this element.

995 This element SHALL be included as a participant in the header of the CDA document in the event of the pregnancy. If the father of the baby is unknown this element SHALL use a null flavor.

6.3.2.6.1 Parent Template

The parent of this template is Patient Contacts (1.3.6.1.4.1.19376.1.5.3.1.2.4).

6.3.2.6.2 Specification

```

1000 <participant typeCode='IND'>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4' />
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4.1' />
<time value='20070213' />
1005 <associatedEntity classCode='PRS'>
  <code code='xx-fatherofbaby' displayName=' ' codeSystem='2.16.840.1.113883.6.96'
    codeSystemName='SNOMED CT' />
    <addr></addr>
    <telecom value=' ' use=' ' />
    <assignedPerson><name></name></assignedPerson>
  </associatedEntity>
1010 </participant>

```

6.3.2.6.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4'><templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.4.2'>

The <templateId> element identifies this person as the natural father and must be recorded exactly as shown above.

```

<rule context='h17:templateId/@root="1.3.6.1.4.1.19376.1.5.3.1.2.4.2"]'>
  <assert test='h17:templateId/@root="1.3.6.1.4.1.19376.1.5.3.1.2.4"'>
    A participant using template 1.3.6.1.4.1.19376.1.5.3.1.2.4.2 must also use template
    1.3.6.1.4.1.19376.1.5.3.1.2.4.
  </assert>
</rule>

```

6.3.2.6.4 <associatedEntity classCode='PRS'>

The classCode attribute of the <associatedEntity> element SHALL be PRS.

```

<rule context='h17:templateId/@root="1.3.6.1.4.1.19376.1.5.3.1.2.4.2"]'>
  <assert test='..../h17:associatedEntity/@classCode = "PRS"'>
    The classCode attribute of the associated entity shall be PRS.
  </assert>
</rule>

```

6.3.2.6.5 <code code='xx-fatherofbaby' displayName=' '

codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT'>

For father of baby the code SHALL be xx-fatherofbaby (requested). The code system name is SNOMED CT.

6.3.2.6.6 Completed Example

```

<!-- Father of baby -->
<participant typeCode="IND">
  <associatedEntity classCode="NOK">
    <code code="xx-fatherofbaby" displayName="Father of Baby"
      codeSystem="2.16.840.1.113883.6.96" codeSystemName="SNOMED CT"/>
    <addr>
      <streetAddressLine>18 Oak Valley Dr.</streetAddressLine>
      <city>Monteagle</city>
      <state>TN</state>
      <postalCode>37205</postalCode>
      <country>USA</country>
    </addr>
    <telecom value="tel:(999) 555-1212" use="WP"/>
    <associatedPerson>
      <name>
        <prefix>Mr.</prefix>
        <given>Thomas</given>
        <family>Caster</family>
      </name>
    </associatedPerson>
  </associatedEntity>
</participant>

```

Add Section 6.3.2.7

6.3.2.7 Authorization 1.3.6.1.4.1.19376.1.5.3.1.2.5

- 1060 Each <authorization> element in the CDA Header represents an informed consent. When the document being shared represents the informed consent to a policy expressed by the XDS Affinity Domain within the document, it shall do so in an <authorization> element. More than one <authorization> element may be present. The consent to share information shall have a unique identifier contained in the <id> element, representing the patient consent to that policy.
- 1065 The policy being consented to shall be represented in the <code> element. Note that other <authorization> elements may be present representing other sorts of consents associated with the document.

6.3.2.7.1 Parent Template

None.

6.3.2.7.2 Specification

- 1070
- ```
<authorization typeCode='AUTH'>
 <consent classCode='CONS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.5' />
 <id root='/'>
 <code code='' codeSystem='' codeSystemName='' displayName='' />
 <statusCode code='completed' />
 </consent>
</authorization>
```
- 1075

1080 **6.3.2.7.3 <authorization typeCode='AUTH'>**

At least one <authorization> element must be present in a consent medical document in documents shared by Document Source actors that implement the privacy option. The typeCode attribute shall be present and be valued with AUTH, indicating that this is an authorization act related to the document.

1085 **6.3.2.7.4 <consent classCode='CONS' moodCode='EVN'>**

Each authorization element shall have one <consent> element. The classCode shall be present and be valued with CONS, indicating that the related act is an informed consent. The moodCode shall be EVN, indicating that this element represents an act that has occurred.

#### **6.3.2.7.5 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.2.5' />**

- 1090 The <templateId> element shall be recorded as shown above and identifies this consent as an authorization entry.

#### **6.3.2.7.6 <id root=' />**

The <consent> element shall have one identifier that is used to uniquely identify the consent act. This identifier shall contain a root attribute, and shall not contain an extension attribute.

**1095 6.3.2.7.7 <code code=' ' codeSystem=' ' codeSystemName=' ' displayName=' '/>**

The <consent> element shall have one <code> element that is used to identify the consent policy that was agreed to by the patient.

*Add Section 6.3.3*

**6.3.3 CDA Section Content Modules****1100 Add Section 6.3.3.1****6.3.3.1 Reasons for Care**

*Add Section 6.3.3.1.1*

**6.3.3.1.1 Reason for Referral**

*Add Section 6.3.3.1.2*

**1105 6.3.3.1.2 Coded Reason for Referral**

*Add Section 6.3.3.1.3*

**6.3.3.1.3 Chief Complaint**

*Add Section 6.3.3.1.4*

**6.3.3.1.4 Hospital Admission Diagnosis****1110 Add Section 6.3.3.1.5****6.3.3.1.5 Proposed Procedure Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.1	
<b>General Description</b>	The proposed procedure section shall contain a description of the procedures for which a risk assessment is required including procedure names and codes, patient position, dates, and names of surgeons. It shall include entries for procedures as described in the Entry Content Modules and the required and optional subsections.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
29554-3	R	PROCEDURE
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.19	R	Procedure Entry
<b>Subsections</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.9.4	R	Reason for Procedure
1.3.6.1.4.1.19376.1.5.3.1.1.9.3	R	Proposed Anesthesia
1.3.6.1.4.1.19376.1.5.3.1.1.9.2	R	Estimated Blood Loss

1.3.6.1.4.1.19376.1.5.3.1.1.9.40	R	Procedure Care Plan
----------------------------------	---	---------------------

```

1115 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.1' />
 <id root=' ' extension=' ' />
 <code code='29554-3' displayName='PROCEDURE'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Required Procedure Entry element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.19' />
 :
 </entry>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.4' />
 <!-- Required Reason for Procedure Section content -->
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.3' />
 <!-- Required Proposed Anesthesia Section content -->
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.2' />
 <!-- Required if known Estimated Blood Loss Section content -->
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.40' />
 <!-- Required if known Procedure Care Plan Section content -->
 </section>
</component>
</section>
</component>

```

**Figure 6.3.3.1.5-1: Specification for Proposed Procedure Section**

1155

Add Section 6.3.3.1.6
-----------------------

**6.3.3.1.6 EBS Estimated Blood Loss Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.2**

Template ID	1.3.6.1.4.1.19376.1.5.3.1.1.9.2		
General Description	The estimated blood loss section shall contain a description of the blood loss for the procedure.		
LOINC Code	Opt	Description	
8717-1	R	OPERATIVE NOTE ESTIMATED BLOOD LOSS	
Entries	Opt	Description	

1.3.6.1.4.1.19376.1.5.3.1.4.13	R	Simple Observation
--------------------------------	---	--------------------

```

1160 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.2' />
<id root=' ' extension=' '/>
<code code='8717-1' displayName='OPERATIVE NOTE ESTIMATED BLOOD LOSS'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
<!-- Required Simple Observation element -->
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
</entry>
</section>
</component>

```

**Figure 6.3.3.1.6-1: EBS Specification for Estimated Blood Loss Section**

Add Section 6.3.3.1.7
-----------------------

**6.3.3.1.7 Proposed Anesthesia Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.3**

Template ID	1.3.6.1.4.1.19376.1.5.3.1.1.9.3	
General Description	The proposed anesthesia section shall contain a description of the anesthetic techniques for which a risk assessment is required. It shall include entries for anesthetic procedures as described in the Entry Content Modules.	
LOINC Code	Opt	Description
10213-7	R	Surgical operation note anesthesia
Entries	Opt	Description
1.3.6.1.4.1.19376.1.5.3.1.4.19	R	Procedure Entry The procedure entries shall be in INT mood.

1180

```

1185 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.3' />
<id root=' ' extension=' ' />
<code code='10213-7' displayName='Surgical operation note anesthesia'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
 <!-- Required Procedure Entry element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.19' />
 :
</entry>
</section>
</component>

```

**Figure 6.3.3.1.7-1: Specification for Anesthesia Administered Section**

1200     *Add Section 6.3.3.1.8*

**6.3.3.1.8 Reason for Procedure Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.4**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.4	
<b>General Description</b>	The reason for procedure section shall contain a description of the reason that the patient is receiving the procedure. It shall include entries for conditions as described in the Entry Content Module.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
10217-8	R	OPERATIVE NOTE INDICATIONS
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.5	R2	Problem Entry

```

1205 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.4' />
<id root=' ' extension=' ' />
<code code='10217-8' displayName='OPERATIVE NOTE INDICATIONS'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
 <!-- Required if known Problem Entry element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5' />
 :
</entry>
</section>
</component>

```

**Figure 6.3.3.1.8-1: Specification for Reason for Procedure Section**

1220     *Add Section 6.3.3.1.9*

**6.3.3.1.9 Reason for Visit Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.1.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.1.1	
<b>General Description</b>	This contains a narrative description of the patient's reason for visit.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
29299-5	R	REASON FOR VISIT

1225

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.1.1' />
 <id root=' ' extension=' ' />
 <code code='29299-5' displayName='REASON FOR VISIT'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

1230

1235

**Figure 6.3.3.1.9-1: Specification for Reason for Visit Section***Add Section 6.3.3.1.10*

1240

**6.3.3.1.10 Injury Incident Description Section 1.3.6.1.4.1.19376.1.5.3.1.1.19.2.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.19.2.1	
<b>General Description</b>	This section shall include a description of the incident leading to the injury, including status of relevant safety equipment in use (e.g., safety belts, air bag, helmet).	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
11374-6	R	Injury incident description

1245

1250

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.19.2.1' />
 <id root=' ' extension=' ' />
 <code code='11374-6' displayName='Injury incident description'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

**Figure 6.3.3.1.10-1: Sample Injury Incident Description Section**

1255

*Add Section 6.3.3.2*

**6.3.3.2 Other Condition Histories**

*Add Section 6.3.3.2.1*

**6.3.3.2.1 History of Present Illness**

1260 *Add Section 6.3.3.2.2*

**6.3.3.2.2 Hospital Course**

*Add Section 6.3.3.2.3*

**6.3.3.2.3 Active Problems**

*Add Section 6.3.3.2.4*

1265 **6.3.3.2.4 Discharge Diagnosis**

*Add Section 6.3.3.2.5*

**6.3.3.2.5 History of Past Illness**

*Add Section 6.3.3.2.6*

**6.3.3.2.6 Encounter Histories**

1270 *Add Section 6.3.3.2.7*

**6.3.3.2.7 History of Outpatient Visits**

*Add Section 6.3.3.2.8*

**6.3.3.2.8 History of Inpatient Visits**

*Add Section 6.3.3.2.9*

1275 **6.3.3.2.9 List of Surgeries**

*Add Section 6.3.3.2.10*

**6.3.3.2.10 Coded List of Surgeries**

*Add Section 6.3.3.2.11*

**6.3.3.2.11 Allergies and Other Adverse Reactions**

1280 *Add Section 6.3.3.2.12*

**6.3.3.2.12 Family medical History**

*Add Section 6.3.3.2.13*

**6.3.3.2.13 Coded Family Medical History**

*Add Section 6.3.3.2.14*

1285 **6.3.3.2.14 Social History Section**

*Add Section 6.3.3.2.15*

**6.3.3.2.15 Functional Status**

*Add Section 6.3.3.2.16*

**6.3.3.2.16 Review of Systems**

1290 *Add Section 6.3.3.2.17*

**6.3.3.2.17 Hazardous Working Conditions**

*Add Section 6.3.3.2.18*

**6.3.3.2.18 Pregnancy History**

*Add Section 6.3.3.2.19*

1295 **6.3.3.2.19 Medical Devices**

*Add Section 6.3.3.2.20*

**6.3.3.2.20 Foreign Travel**

*Add Section 6.3.3.2.21*

**6.3.3.2.21 Pre-procedure Family Medical History Section**

1300 **1.3.6.1.4.1.19376.1.5.3.1.1.9.5 (Deprecated)**

*Add Section 6.3.3.2.22*

**6.3.3.2.22 Coded Functional Status Assessment Section**

**1.3.6.1.4.1.19376.1.5.3.1.1.12.2.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.12.2.1
<b>Parent Template</b>	Functional Status (1.3.6.1.4.1.19376.1.5.3.1.3.17, see PCC TF-2: 6.3.3.2.15)
<b>General Description</b>	The coded functional status assessment section provides a machine readable and narrative description of the patient's status of normal functioning at the time the document was created.

	Functional status includes information concerning: Ambulatory ability Mental status or competency Activities of Daily Living (ADL's) including bathing, dressing, feeding, grooming Home/living situation having an effect on the health status of the patient Ability to care for self Social activity, including issues with social cognition, participation with friends and acquaintances other than family members Occupation activity, including activities partly or directly related to working, housework or volunteering, family and home responsibilities or activities related to home and family Communication ability, including issues with speech, writing or cognition required for communication Perception, including sight, hearing, taste, skin sensation, kinesthetic sense, proprioception, or balance	
LOINC Code	Opt	Description
47420-5	R	Functional Status Assessment
Subsections	Opt	Description
1.3.6.1.4.1.19376.1.5.3.1.1.12.2.2	R	<a href="#">Pain Scale Assessment</a>
1.3.6.1.4.1.19376.1.5.3.1.1.12.2.3	O <sup>Note 1</sup>	<a href="#">Braden Score Assessment</a>
1.3.6.1.4.1.19376.1.5.3.1.1.12.2.4	O <sup>Note 1</sup>	<a href="#">Geriatric Depression Scale</a>
1.3.6.1.4.1.19376.1.5.3.1.1.12.2.5	O <sup>Note 1</sup>	<a href="#">Minimum Data Set</a>

Note 1: At least one of the above optional subsections shall be present

1305

### 6.3.3.2.22.1 Standards

**CDAR2** [HL7 CDA Release 2.0](#)

**CRS** [HL7 Care Record Summary](#)

**CCD** [ASTM/HL7 Continuity of Care Document](#)

**LOINC** [Logical Observation Identifier Names and Codes](#)

**SNOMED** [Systemitized Nomenclature of Medicine Clinical Terminology](#)

### 6.3.3.2.22.2 Parent Template

The parent of this template is Functional Status (see PCC TF-2: 6.3.3.2.15).

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.17' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.2.1' />
 <id root=' ' extension=' ' />
 <code code='47420-5' displayName='Functional Status Assessment'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.2.2' />
 <!-- Required Pain Scale Assessment Section content -->
 </section>
 </component>
 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.2.3' />
 <!-- Optional Braden Score Assessment Section content -->
 </section>
 </component>
 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.2.4' />
 <!-- Optional Geriatric Depression Scale Section content -->
 </section>
 </component>
 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.2.5' />
 <!-- Optional Minimum Data Set Section content -->
 </section>
 </component>
 </section>
</component>

```

**Figure 6.3.3.2.22.2-1: Specification for Coded Functional Status Assessment Section**

1350

Add Section 6.3.3.2.23

**6.3.3.2.23 Pain Scale Assessment Section 1.3.6.1.4.1.19376.1.5.3.1.1.12.2.2**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.12.2.2	
<b>General Description</b>	The Pain Scale Assessment contains a coded observation reflecting the patient's reported intensity of pain on a scale from 0 to 10.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
38208-5	R	Pain severity
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.12.3.1	R	<a href="#">Pain Score Observation</a>

```

1355 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.2.2' />
<id root=' ' extension=' ' />
<code code='38208-5' displayName='Pain severity'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
 <!-- Required Pain Score Observation element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.1' />
 :
</entry>
</section>
</component>

```

**Figure 6.3.3.2.23-1: Specification for Pain Scale Assessment Section****Add Section 6.3.3.2.24****1375 6.3.3.2.24 Braden Score Section 1.3.6.1.4.1.19376.1.5.3.1.1.12.2.3**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.12.2.3	
<b>General Description</b>	This section reports the Braden score and its related assessments in machine and human readable form.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
38228-3	R	BRADEN SCALE SKIN ASSESSMENT PANEL
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.12.3.2	R	<a href="#">Braden Score Observation</a>

1380

```

1385 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.2.3' />
<id root=' ' extension=' ' />
<code code='38228-3' displayName='BRADEN SCALE SKIN ASSESSMENT PANEL'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
 <!-- Required Braden Score Observation element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.2' />
 :
</entry>
</section>
</component>

```

**Figure 6.3.3.2.24-1: Specification for Braden Score Section**

1405 *Add Section 6.3.3.2.25*

**6.3.3.2.25 Geriatric Depression Scale Section 1.3.6.1.4.1.19376.1.5.3.1.1.12.2.4**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.12.2.4	
<b>General Description</b>	This section reports the Geriatric Depression Scale score and its related assessments in machine and human readable form.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
48542-5	R	Geriatric Depression Scale (GDS) Panel
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.12.3.4	R	<a href="#">Geriatric Depression Score Observation</a>

```

1410 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.2.4' />
<id root=' ' extension=' ' />
<code code='48542-5' displayName='Geriatric Depression Scale (GDS) Panel'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
 <!-- Required Geriatric Depression Score Observation element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.4' />
 :
</entry>
</section>
</component>

```

**Figure 6.3.3.2.25-1: Specification for Geriatric Depression Scale Section**

Add Section 6.3.3.2.26
------------------------

### 6.3.3.2.26 Physical Function Section 1.3.6.1.4.1.19376.1.5.3.1.1.12.2.5

1430

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.12.2.5	
<b>General Description</b>	This section reports scores from section G of the Minimum Data Set.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
46006-3	R	Physical functioning and structural problems
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.12.3.7	O	<u><a href="#">Survey Panel</a></u> At least one Survey Panel or Survey Observation shall be present.
1.3.6.1.4.1.19376.1.5.3.1.1.12.3.6	O	<u><a href="#">Survey Observation</a></u> At least one Survey Panel or Survey Observation shall be present.

1435

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.2.5' />
 <id root=' ' extension=' ' />
 <code code='46006-3' displayName='Physical functioning and structural problems'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Optional Survey Panel element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.7' />
 :
 </entry>
 <entry>
 :
 <!-- Optional Survey Observation element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.6' />
 :
 </entry>
 </section>
</component>

```

1440

1445

1450

1455

**Figure 6.3.3.2.26-1: Specification for Physical Function Section**

#### 6.3.3.2.26.1 Constraints

1460

[Survey Panel](#) found in this section shall be identified using the panel codes found in the table below, and shall contain one or more survey observations from that panel.

[Survey Observation](#) found in this section shall use the LOINC codes from table 6.3.3.2.26.1 to express the answer to one or more questions from the Minimum Data Set Section G. The Survey Observations shall not contain a <methodCode> or <targetSiteCode> element, as these are not appropriate to the MDS Survey instrument.

1465

**Table 6.3.3.2.26.1-1: Panel Codes**

<b>Panel Code</b>	<b>Observation Code</b>	<b>Description</b>	<b>Data Type</b>	<b>Value Set</b>
46007-1	Panel	ADL self performance or support		
	45588-1	Bed mobility - self-performance	CO	2.16.840.1.113883.6.257.755
	45589-9	Bed mobility - support provided	CO	2.16.840.1.113883.6.257.768
	45590-7	Transfer - self-performance	CO	2.16.840.1.113883.6.257.755
	45591-5	Transfer - support provided	CO	2.16.840.1.113883.6.257.768
	45592-3	Walk in room - self-performance	CO	2.16.840.1.113883.6.257.755
	45593-1	Walk in room - support provided	CO	2.16.840.1.113883.6.257.768
	45594-9	Walk in corridor - self-performance	CO	2.16.840.1.113883.6.257.755
	45595-6	Walk in corridor - support provided	CO	2.16.840.1.113883.6.257.768
	45596-4	Locomotion on unit - self-performance	CO	2.16.840.1.113883.6.257.755
	45597-2	Locomotion on unit - support provided	CO	2.16.840.1.113883.6.257.768
	45598-0	Locomotion off unit - self-performance	CO	2.16.840.1.113883.6.257.755
	45599-8	Locomotion off unit - support provided	CO	2.16.840.1.113883.6.257.768
	45600-4	Dressing - self-performance	CO	2.16.840.1.113883.6.257.755
	45601-2	Dressing - support provided	CO	2.16.840.1.113883.6.257.768
	45602-0	Eating - self-performance	CO	2.16.840.1.113883.6.257.755
	45603-8	Eating - support provided	CO	2.16.840.1.113883.6.257.768
	45604-6	Toilet use - self-performance	CO	2.16.840.1.113883.6.257.755
	45605-3	Toilet use - support provided	CO	2.16.840.1.113883.6.257.768
	45606-1	Personal hygiene - self-performance	CO	2.16.840.1.113883.6.257.755
	45607-9	Personal hygiene - support provided	CO	2.16.840.1.113883.6.257.768
46008-9	Panel	Bathing		
	45608-7	Bathing - self-performance	CO	2.16.840.1.113883.6.257.860
	45609-5	Bathing - support provided	CO	2.16.840.1.113883.6.257.768
46009-7	Panel	Test for balance		
	45610-3	Balance while standing	CO	2.16.840.1.113883.6.257.876
	45523-8	Balance while sitting	CO	2.16.840.1.113883.6.257.876
46010-5	Panel	Functional limitation in range of motion		
	45524-6	Range of motion^Neck	CO	2.16.840.1.113883.6.257.889
	45525-3	Voluntary movement^Neck	CO	2.16.840.1.113883.6.257.898
	45526-1	Range of motion^Upper Extremity	CO	2.16.840.1.113883.6.257.889
	45527-9	Voluntary movement^Upper Extremity	CO	2.16.840.1.113883.6.257.898
	45528-7	Range of motion^Hand	CO	2.16.840.1.113883.6.257.889

<b>Panel Code</b>	<b>Observation Code</b>	<b>Description</b>	<b>Data Type</b>	<b>Value Set</b>
	45529-5	Voluntary movement^Hand	CO	2.16.840.1.113883.6.257.898
	45530-3	Range of motion^Lower Extremity	CO	2.16.840.1.113883.6.257.889
	45531-1	Voluntary movement^Lower Extremity	CO	2.16.840.1.113883.6.257.898
	45532-9	Range of motion^Foot	CO	2.16.840.1.113883.6.257.889
	45533-7	Voluntary movement^Foot	CO	2.16.840.1.113883.6.257.898
	45534-5	Other - range of motion	CO	2.16.840.1.113883.6.257.889
	45535-2	Other - voluntary movement	CO	2.16.840.1.113883.6.257.898
46011-3	Panel	Modes of locomotion		
	45536-0	Uses cane, walker or crutch	CO	2.16.840.1.113883.6.257.117
	45537-8	Wheeled self	CO	2.16.840.1.113883.6.257.117
	45538-6	Other person wheeled	CO	2.16.840.1.113883.6.257.117
	45539-4	Uses wheelchair for primary locomotion	CO	2.16.840.1.113883.6.257.117
	45540-2	No modes of locomotion	CO	2.16.840.1.113883.6.257.117
46012-1	Panel	Modes of transfer		
	45541-0	Bedfast all or most of the time	CO	2.16.840.1.113883.6.257.117
	45542-8	Bed rails for bed mobility or transfer	CO	2.16.840.1.113883.6.257.117
	45543-6	Lifted manually	CO	2.16.840.1.113883.6.257.117
	45544-4	Lifted mechanically	CO	2.16.840.1.113883.6.257.117
	45545-1	Transfer aid	CO	2.16.840.1.113883.6.257.117
	45546-9	No mode of transfer	CO	2.16.840.1.113883.6.257.117
No Panel	45611-1	Task segmentation	CO	2.16.840.1.113883.6.257.117
46013-9	Panel	ADL functional rehabilitation potential		
	45612-9	Resident sees increased independence capability	CO	2.16.840.1.113883.6.257.117
	45613-7	Staff sees increased independence capability	CO	2.16.840.1.113883.6.257.117
	45614-5	Resident slow performing tasks or activity	CO	2.16.840.1.113883.6.257.117
	45615-2	Difference in morning to evening activities of daily living	CO	2.16.840.1.113883.6.257.117
	45616-0	Activities of daily living rehabilitation potential - none of above	CO	2.16.840.1.113883.6.257.117
	45617-8	Change in activities of daily living function	CO	2.16.840.1.113883.6.257.464

The coded original values used in the observations above are described in more detail in the table below.

1470

<b>Explanation</b>	<b>Coded Value</b>
2.16.840.1.113883.6.257.755	
INDEPENDENT-No help or oversight -OR- Help/oversight provided only 1 or 2 times during last 7 days	0
SUPERVISION-Oversight, encouragement or cueing provided 3 or more times during last 7 days -OR- Supervision (3 or more times) plus physical assistance provided only 1 or 2 times during last 7 days	1
LIMITED ASSISTANCE-Resident highly involved in activity; received physical help in guided maneuvering of limbs or other nonweight bearing assistance 3 or more times - OR-More help provided only 1 or 2 times during last 7 days	2
EXTENSIVE ASSISTANCE-While resident performed part of activity, over last 7-day period, help of following type(s) provided 3 or more times: - Weight-bearing support - Full staff performance during part (but not all) of last 7 days	3
TOTAL DEPENDENCE-Full staff performance of activity during entire 7 days	4
ACTIVITY DID NOT OCCUR during entire 7 days	8
2.16.840.1.113883.6.257.768	
No setup or physical help from staff	0
Setup help only	1
One person physical assist	2
ADL activity itself did not occur during entire 7 days	8
2.16.840.1.113883.6.257.860	
Independent-No help provided	0
Supervision-Oversight help only	1
Physical help limited to transfer only	2
Physical help in part of bathing activity	3
Total dependence	4
Activity itself did not occur during entire 7 days	8
2.16.840.1.113883.6.257.876	
Maintained position as required in test	0
Unsteady, but able to rebalance self without physical support	1
Partial physical support during test; or stands (sits) but does not follow directions for test	2
Not able to attempt test without physical help	3
2.16.840.1.113883.6.257.889	
No limitation	0
Limitation on one side	1
Limitation on both sides	2
2.16.840.1.113883.6.257.898	
No loss	0
Partial loss	1
Full loss	2

Explanation	Coded Value
2.16.840.1.113883.6.257.117	
No	0
Yes	1
UTD	-
2.16.840.1.113883.6.257.464	
No change	0
Improved	1
Deteriorated	2

*Add Section 6.3.3.2.27***6.3.3.2.27 Preprocedure Review of Systems Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.13**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.13	
<b>Parent Template</b>	Review of Systems (1.3.6.1.4.1.19376.1.5.3.1.3.18)	
<b>General Description</b>	The pre-procedure review of systems section shall contain only required and optional subsections dealing with the responses the patient gave to a set of routine questions on body systems in general and specific risks of anesthesia not covered in general review of systems.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
10187-3	R	REVIEW OF SYSTEMS
<b>Subsections</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.9.46	R	Implanted Medical Device Review
1.3.6.1.4.1.19376.1.5.3.1.1.9.47	R2	Pregnancy Status Review
1.3.6.1.4.1.19376.1.5.3.1.1.9.14	R	Anesthesia Risk Review of Systems

```

1475 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.18' />
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.13' />
<id root=' ' extension=' ' />
<code code='10187-3' displayName='REVIEW OF SYSTEMS'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.46' />
 <!-- Required Implanted Medical Device Review Section content -->
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.47' />
 <!-- Required if known Pregnancy Status Review Section content -->
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.14' />
 <!-- Required Anesthesia Risk Review of Systems Section content -->
 </section>
</component>

1500 </section>
</component>

1505 </section>
</component>

```

**Figure 6.3.3.2.27-1: Specification for Preprocedure Review of Systems Section****Add Section 6.3.3.2.28****6.3.3.2.28 Estimated Delivery Date Section 1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.1	
<b>General Description</b>	This section contains the physician's best estimate of the patient's due date. This is generally done both on an initial evaluation, and later confirmed at 18-20 weeks. The date is supported by evidence such as the patient's history of last menstrual period, a physical examination, or ultrasound measurements.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57060-6	R	Estimated date of delivery
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.1	R	<u><a href="#">Estimated Delivery Date Observation</a></u> This is a simple observation to represent the estimated due date with a supporting observation or observations that state the method used and date implied by that method. If one observation is present, then it is to be interpreted as the initial EDD. If the initial observation dates indicate the EDD is within the 18 to 20 weeks completed gestation, that observation will also populate the 18-20 week update. If the initial observation indicates an EDD of more than 20 weeks EGA, then no value will be placed in the 18-20 week update field.

```

1515 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.1'>
 <id root=' ' extension=' ' />
 <code code='57060-6' displayName='Estimated date of delivery'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Required Estimated Due Date Observation element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.1'>
 :
 </entry>
 </section>
</component>

```

**Figure 6.3.3.2.28-1: Specification for Estimated Delivery Dates Section***Add Section 6.3.3.2.29***6.3.3.2.29 History of Tobacco Use Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.8**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.8	
<b>General Description</b>	The history of tobacco use section shall contain a description of the responses the patient gave to a set of routine questions on the history of tobacco use.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
11366-2	R	HISTORY OF TOBACCO USE

```

1535 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.8'>
 <id root=' ' extension=' ' />
 <code code='11366-2' displayName='HISTORY OF TOBACCO USE'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>

```

**Figure 6.3.3.2.29-1: Specification for History of Tobacco Use Section***Add Section 6.3.3.2.30*

**6.3.3.2.30 Current Alcohol/Substance Abuse Section****1.3.6.1.4.1.19376.1.5.3.1.1.9.10**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.10	
<b>General Description</b>	The history of alcohol/substance abuse section shall contain a description of the responses the patient gave to a set of routine questions on the current abuse of alcohol or other substances.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
18663-5	R	HISTORY OF PRESENT ALCOHOL AND/OR SUBSTANCE ABUSE

1555

1560

1565

1570

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.10'>
 <id root=' ' extension=' '/>
 <code code='18663-5' displayName='HISTORY OF PRESENT ALCOHOL AND/OR SUBSTANCE ABUSE'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'>
 <text>
 Text as described above
 </text>
 </code>
 </section>
 </component>

```

**Figure 6.3.3.2.30-1: Specification for Current Alcohol/Substance Abuse Section***Add Section 6.3.3.2.31***6.3.3.2.31 History of Blood Transfusion Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.12**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.12	
<b>General Description</b>	The History of Blood Transfusion section shall contain a narrative description of the blood products the patient has received in the past, including any reactions to blood products.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
56836-0	R	History of blood transfusion

1575

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.12' />
 <id root=' ' extension=' ' />
 <code code='56836-0' displayName='History of blood transfusion'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>

```

**Figure 6.3.3.2.31-1: Specification for History of Blood Transfusion Section***Add Section 6.3.3.2.32***6.3.3.2.32 Anesthesia Risk Review of Systems Section****1.3.6.1.4.1.19376.1.5.3.1.1.9.14**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.14	
<b>Parent Template</b>	1.3.6.1.4.1.19376.1.5.3.1.3.18	
<b>General Description</b>	The anesthesia review of systems section shall contain a description of the responses the patient gave to a set of routine questions on specific risks of anesthesia not covered in general review of systems such as broken teeth, airway limitations, positioning limitations, recent infections, and history of personal anesthesia problems..	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57081-2	R	Anesthesia Risk Review of Systems

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.18' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.14' />
 <id root=' ' extension=' ' />
 <code code='57081-2' displayName='Anesthesia Risk Review of Systems'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>

```

**Figure 6.3.3.2.32-1: Specification for Anesthesia Risk Review of Systems Section***Add Section 6.3.3.2.33*

1615 **6.3.3.2.33 Implanted Medical Device Review Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.46**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.46	
<b>General Description</b>	The implanted medical device review section shall contain a description of the medical devices that are inserted into the patient, whether internal or partially external.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57080-4	R	Implanted medical device

1620 <component>

```

<section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.46' />
 <id root=' ' extension=' ' />
 <code code='57080-4' displayName='Implanted medical device'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
</section>
</component>
```

1625

1630

**Figure 6.3.3.2.33-1: Specification for Implanted Medical Device Review Section**

Add Section 6.3.3.2.34
------------------------

1635 **6.3.3.2.34 Pregnancy Status Review Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.47**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.47	
<b>General Description</b>	The pregnancy status review section shall contain a description of the responses the patient gave to a set of routine questions regarding potential pregnancy in females of child-bearing-age.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
11449-6	R	Pregnancy Status-Reported

1640 <component>

```

<section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.47' />
 <id root=' ' extension=' ' />
 <code code='11449-6' displayName='Pregnancy Status-Reported'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
</section>
</component>
```

1645

1650

**Figure 6.3.3.2.34-1: Specification for Pregnancy Status Review Section**

Add Section 6.3.3.2.35
------------------------

### 6.3.3.2.35 History of Infection Section 1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1	
<b>General Description</b>	The History of Infection section shall contain a narrative description of any infections the patient may have contracted prior to the patient's current visit or admission.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
56838-6	R	History of infectious disease

1655

1660      <component>  
           <section>  
           <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1'>  
           <id root=' ' extension=' ' />  
           <code code='56838-6' displayName='History of infectious disease'  
               codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />  
           <text>  
               Text as described above  
           </text>  
           </section>  
       </component>

1665

**Figure 6.3.3.2.35-1: Specification for History of Infection Section**

1670

Add Section 6.3.3.2.36
------------------------

### 6.3.3.2.36 Coded Social History Section 1.3.6.1.4.1.19376.1.5.3.1.3.16.1

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.3.16.1	
<b>Parent Template</b>	Social History (1.3.6.1.4.1.19376.1.5.3.1.3.16)	
<b>General Description</b>	The social history section shall contain a narrative description of the person's beliefs, home life, community life, work life, hobbies, and risky habits. It shall include Social History Observations.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
29762-2	R	SOCIAL HISTORY
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.13.4	R	Social History Observation

```

1675 <component>
<section>
<templateId root='11.3.6.1.4.1.19376.1.5.3.1.3.16' />
<templateId root='11.3.6.1.4.1.19376.1.5.3.1.3.16.1' />
<id root=' ' extension=' ' />
<code code='29762-2' displayName='SOCIAL HISTORY'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
</section>
</component>

```

**Figure 6.3.3.2.36-1: Specification for Coded Social History Section***Add Section 6.3.3.2.37***6.3.3.2.37 Coded History of Infection Section 1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1.1	
<b>Parent Template</b>	History of Infection (1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1)	
<b>General Description</b>	The History of Infection section shall contain a narrative description of any infections the patient may have contracted prior to the patient's current condition. It shall include entries for problems as described in the Entry Content Modules.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
56838-6	R	History of infectious disease
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.5.2	R	<a href="#">Problem Concern Entry</a>

1690

```

1695 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1.1.1' />
<id root=' ' extension=' ' />
<code code='56838-6' displayName='History of infectious disease'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
</section>
</component>

```

**Figure 6.3.3.2.37-1: Specification for Coded History of Infection Section**

1705

*Add Section 6.3.3.2.38***6.3.3.2.38 Prenatal Events Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.2**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3..1.1.21.2.2
<b>General Description</b>	The Prenatal Events Section shall include narrative text describing pertinent prenatal information that has a direct impact on the process of labor and

		delivery. It shall also include subsections if known.
LOINC Code	Opt	Description
57073-9	R	Prenatal events
Subsections	Opt	Description
<b>Coded Results</b> This section SHOULD contain laboratory results and procedures as pertaining to the pregnancy , e.g., amniocentesis, cordocentesis, chorionic villus sampling.	R2	1.3.6.1.4.1.19376.1.5.3.1.3.28
<b>Procedures and Interventions</b> This section SHOULD contain procedures that took place during the prenatal period (i.e. prenatal care, prenatal complications, prenatal surgeries)	R2	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11
<b>Event Outcomes</b> This section contains event outcomes related to prenatal events e.g. miscarriage, infection.	R2	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.9

1710	<component> <section> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.2' /> <id root=' ' extension=' ' /> <code code='57073-9' displayName='Prenatal events' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' /> <text> Text as described above </text>
1715	</component> <component> <section> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.28' /> <!-- Required if known Coded Results Section --&gt;</td>
1720	</section> </component> <component> <section> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11' /> <!-- Required if known Procedures and Interventions Section --&gt;</td>
1725	</section> </component> <component> <section> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.9' /> <!-- Required if known Event Outcomes Section --&gt;</td>
1730	</section> </component>
1735	</component> </section> </component>

**Figure 6.3.3.2.38-1: Specification for Prenatal Events Section****Add Section 6.3.3.2.39**

1740

**6.3.3.2.39 Labor and Delivery Events Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.3**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.3	
<b>Parent Template</b>		
<b>General Description</b>	The Labor and Delivery Events Section SHALL include a narrative text containing relevant information collected during the labor and delivery process.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57074-7	R	Labor and delivery process
<b>Subsections</b>	<b>Opt</b>	<b>Description</b>
<b>Procedures and Interventions</b>  The subsection SHALL contain procedures and interventions specific to labor and delivery events. These may include induction, the delivery type (e.g. vaginal, vaginal birth after cesarean section or cesarean section along with incision type), electronic fetal monitoring, etc.	R2	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11
<b>Coded Event Outcomes</b>  This section SHOULD contain outcomes related to the labor and delivery process such as live birth or stillborn.  The subsection shall include coded event outcomes such as live birth or stillborn and also including maternal death with date/time. Furthermore, Coded Event Outcomes section shall contain a simple Observation using LOINC Code 11636-8 that reports the number of births live or dead that occurred during the delivery event.,	R2	1.3.6.1.4.1.19376.1.7.3.1.1.13.7

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.3'>
 <id root=' ' extension=' ' />
 <code code='57074-7' displayName='Labor and delivery process'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
 </component>
 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11' />
 <!-- Required if known Procedures and Interventions Section -->
 </section>
 </component>
 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.13.7' />
 <!-- Required if known Coded Event Outcomes Section -->
 </section>
 </component>
</component>

```

**Figure 6.3.3.2.39-1: Specification for Labor and Delivery Process Section**

<i>Add Section 6.3.3.2.40</i>
-------------------------------

**6.3.3.2.40 Newborn Delivery Information Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.4**

Template ID	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.4	
General Description	The Newborn Delivery Information Section SHALL include a narrative text containing information collected at the birth and up to the transfer of the infant from the birthing room to a post-natal unit.	
LOINC Code	Opt	Description
57075-4	R	Newborn delivery information from newborn
Subsections	Opt	Description
<b>Coded Detailed Physical Examination Section</b>  This section SHALL include information about the newborn genitalia; weight; length; head circumference, size (AGA, SGA or LGA); Apgar score assessment ; vital signs, physical exam findings	R	1.3.6.1.4.1.19376.1.5.3.1.1.9.15.1
<b>Active Problems</b>  This section SHALL describe problems that the newborn might have had during or immediately prior to delivery.	R2	1.3.6.1.4.1.19376.1.5.3.1.3.6
<b>Procedures and Interventions</b>  This section SHALL include the procedures and interventions received by the newborn such as suction or resuscitation.	R2	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11

<b>Medications Administered</b> This section SHALL include the medication that was administered to the newborn while in the birthing suite such as: Vitamin K (Aquamephyton) injection; erythromycin eye ointment; and resuscitation medications (if any) including date, time, and route of administration.	R2	1.3.6.1.4.1.19376.1.5.3.1.3.21
<b>Event Outcomes</b> This section SHALL include the outcomes of the procedures and interventions such as a resuscitation event.	R2	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.9
<b>Coded Event Outcomes</b>	C	1.3.6.1.4.1.19376.1.7.3.1.1.13.7
<b>Coded Results</b>	C	1.3.6.1.4.1.19376.1.5.3.1.3.28
<b>Intake and Output</b> This section SHALL include any intake and output while the newborn is in the delivery suite (excluding estimated blood loss) such as: first urine/void; stool; gastric output	C	1.3.6.1.4.1.19376.1.5.3.1.1.20.2.3

1770

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.4' />
 <id root=' ' extension=' ' />
 <code code='57075-4' displayName='Newborn delivery information from newborn'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.15.1' />
 <!-- Required Coded Detailed Physical Examination Section -->
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.6' />
 <!-- Required if known Active Problems Section -->
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11' />
 <!-- Required if known Procedures and Interventions Section -->
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.21' />
 <!-- Required if known Medications Administered Section -->
 </section>
</component>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.9' />
 <!-- Required if known Event Outcomes Section -->
 </section>
</component>
</component>

```

**Figure 6.3.3.2.40-1: Specification for Newborn Delivery Information Section**

<i>Add Section 6.3.3.2.41</i>
-------------------------------

1815

**6.3.3.2.41 Postpartum Hospitalization Treatment Section****1.3.6.1.4.1.19376.1.5.3.1.1.21.2.7**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.7	
<b>Parent Template</b>		
<b>General Description</b>	The Postpartum Treatment Section shall include a narrative description of the treatment delivered to the mother subsequent to the delivery.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57076-2	R	Postpartum hospitalization TREATMENT
<b>Subsections</b>	<b>Opt</b>	<b>Description</b>
<b>Immunizations</b> This section SHOULD contain the immunization given to the mother prior to the discharge from the birthing facility.	O	1.3.6.1.4.1.19376.1.5.3.1.4.12
<b>Medications Administered</b> This SHOULD include commonly prescribed maternal medications including contraceptive medication.	R2	1.3.6.1.4.1.19376.1.5.3.1.3.21
<b>Procedures and Interventions</b> This section SHALL include the procedures and interventions received by the mother during the immediate post-partum period e.g., transfusion or curettage.	R	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11
<b>Coded Results</b> This section SHOULD contain laboratory results and procedures as pertaining to the mother while discharged such as the hemoglobin or the hematocrit level.	R2	1.3.6.1.4.1.19376.1.5.3.1.3.28
<b>Care plan</b> This section SHOULD include the plan of care for the mother upon her discharge such as the feeding method or the contraceptive plan	O	1.3.6.1.4.1.19376.1.5.3.1.3.31
<b>Discharge Diet</b> This section SHALL include the diet that the mother was recommended upon her discharge.	R	1.3.6.1.4.1.19376.1.5.3.1.3.33

```

1820 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.7' />
<id root=' ' extension=' ' />
<code code='57076-2' displayName='POST PARTUM HOSPITALIZATION TREATMENT'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
1825 <text>
 Text as described above
</text>
<component>
<section>
1830 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.6' />
 <!-- Required Active Problems Section -->
</section>
</component>
<component>
<section>
1835 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.12' />
 <!-- Optional Immunizations Section -->
</section>
</component>
<component>
<section>
1840 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.22' />
 <!-- Required if known Hospital Discharge Medication Section -->
</section>
</component>
<component>
<section>
1845 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.13.2.11' />
 <!-- Required Procedures and Interventions Section -->
</section>
</component>
<component>
<section>
1850 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.28' />
 <!-- Required if known Coded Results Section -->
</section>
</component>
<component>
<section>
1855 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.31' />
 <!-- Optional Care Plan Section -->
</section>
</component>
<component>
<section>
1860 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.33' />
 <!-- Required Discharge Diet Section -->
</section>
</component>
</component>
1865 </component>
</section>
</component>
1870 </component>

```

**Figure 6.3.3.2.41-1: Specification for Postpartum Treatment Section**

<i>Add Section 6.3.3.2.42</i>
-------------------------------

1875

**6.3.3.2.42 Event Outcomes Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.9**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.9	
<b>Parent Template</b>		
<b>General Description</b>	The Event Outcome Section shall include a narrative description of the outcomes following a procedure, an intervention or a problem.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
42545-4	R	EVENT OUTCOME

1880

<pre>&lt;component&gt;   &lt;section&gt;     &lt;templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.9' /&gt;     &lt;id root=' ' extension=' '/&gt;     &lt;code code='42545-4' displayName='EVENT OUTCOME'       codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' /&gt;     &lt;text&gt;       Text as described above     &lt;/text&gt;   &lt;/section&gt; &lt;/component&gt;</pre>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

1885

**Figure 6.3.3.2.42-1: Specification for Event Outcomes Section**

<i>Add Section 6.3.3.2.43</i>
-------------------------------

1890

**6.3.3.2.43 Newborn Status at Maternal Discharge 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.8**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.8	
<b>Parent Template</b>		
<b>General Description</b>	The Newborn Status and Maternal Discharge section shall contain a narrative description of the status and disposition of the newborn at the time of maternal discharge.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57077-0	R	Newborn status at maternal discharge from newborn

1895

<pre>&lt;component&gt;   &lt;section&gt;     &lt;templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.8' /&gt;     &lt;id root=' ' extension=' '/&gt;     &lt;code code='57077-0' displayName='Newborn status at maternal discharge from newborn'       codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' /&gt;     &lt;text&gt;       Text as described above     &lt;/text&gt;   &lt;/section&gt; &lt;/component&gt;</pre>
-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

1900

**Figure 6.3.3.2.43-1: Specification for Newborn Status at Maternal Discharge Section**

*Add Section 6.3.3.2.44*

1905 **6.3.3.2.44 History of Surgical Procedures Section 1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2	
<b>Parent Template</b>		
<b>General Description</b>	The History of Surgical Procedures Section shall contain a narrative description of the surgical procedures performed on the patient.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
10167-5	R	History of surgical procedures

1910      <component>  
           <section>  
           <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.16.2.2' />  
           <id root=' ' extension=' '/>  
           <code code='10167-5' displayName='History of surgical procedures'  
               codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />  
           <text>  
             Text as described above  
           </text>  
         </component>

1915

**Figure 6.3.3.2.44-1: Specification for History of Surgical Procedures Section**

1920 *Add Section 6.3.3.2.45*

1920 **6.3.3.2.45 Operative Note Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.6**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.6	
<b>Parent Template</b>		
<b>General Description</b>	The Operative Note Section shall contain a narrative description of the current operation or surgical procedure in detail.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
10223-6	R	Surgical operation note surgical procedure

1925      <component>  
           <section>  
           <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.6' />  
           <id root=' ' extension=' '/>  
           <code code='10223-6' displayName='Surgical operation note surgical procedure'  
               codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />  
           <text>  
             Text as described above  
           </text>  
         </component>

1930

**Figure 6.3.3.2.45-1: Specification for Operative Note Section**

1935

Add Section 6.3.3.2.46
------------------------

### 6.3.3.2.46 Child Functional Status Assessment 1.3.6.1.4.1.19376.1.7.3.1.1.13.3

<b>Template ID</b>	1.3.6.1.4.1.19376.1.7.3.1.1.13.3	
<b>General Description</b>	This section provides a description of the child's status of normal functioning at the time the document was created. This section includes the psychomotor and the eating and sleeping assessments. This section shall include the Psychomotor Test Observation entry.	
LOINC Code	Opt	Description
47420-5	R	Functional Status Assessment
Subsections	Opt	Description
1.3.6.1.4.1.19376.1.7.3.1.1.13.4	O	Psychomotor Development
1.3.6.1.4.1.19376.1.7.3.1.1.13.5	O	Eating and sleeping assessment

#### Example

1940

```
<component>
 <section>
 <templateId root="1.3.6.1.4.1.19376.1.7.3.1.1.13.3"/>
 <id root="16696797-f854-443d-8819-231ee09cad71"/>
 <code code="47420-5" displayName="Functional Status Assessment"
 codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/>
 <title/>
 <text/>
 <component>
 <section>
 <!-- Optional Psychomotor Development section -->
 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.13.4' />
 :
 </section>
 </component>
 <component>
 <section>
 <!-- Eating and sleeping assessment section -->
 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.13.5' />
 :
 </section>
 </component>
 </section>
 </component>
```

1950

1955

1960

1965

Add Section 6.3.3.2.47
------------------------

### 6.3.3.2.47 Psychomotor Development Section 1.3.6.1.4.1.19376.1.7.3.1.1.13.4

<b>Template ID</b>	1.3.6.1.4.1.19376.1.7.3.1.1.13.4
<b>General Description</b>	This section describes a test battery in order to evaluate the psychomotricity of the newborn.

<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
xx-MCH-PsychoMDev	R	Psychomotor development
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
Simple Observation	R	1.3.6.1.4.1.19376.1.5.3.1.4.13

1970	<component> <section> <!-- Psychomotor Development section templateId --> <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.13.4' /> <id root=' ' extension=' '/> <code code='47420-5' displayName=' Functional Status Assessment ' codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" /> <text> : </text> <entry> : <!--Required simple Observation element --> <templateId root=1.3.6.1.4.1.19376.1.5.3.1.4.13> : </entry> </section> </component>
1975	
1980	
1985	

Add Section 6.3.3.2.48
------------------------

1990    **6.3.3.2.48 Eating and Sleeping Assessment Section 1.3.6.1.4.1.19376.1.7.3.1.1.13.5**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.7.3.1.1.13.5	
<b>General Description</b>	This section describes a test battery in order to evaluate the psychomotricity of the newborn.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
47420-5	R	Functional Status Assessment
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
Simple Observation	R	1.3.6.1.4.1.19376.1.5.3.1.4.13

1995	<component> <section> <!--Eating and Sleeping assessment section templateId --> <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.13.5' /> <id root=' ' extension=' '/> <code code='47420-5' displayName=' Functional Status Assessment ' codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" /> <text> : </text> <entry> : <!-- Required Simple Observation element --> <templateId root=1.3.6.1.4.1.19376.1.5.3.1.4.13/> : </entry> </section> </component>
2000	
2005	
2010	

*Add Section 6.3.3.2.49***6.3.3.2.49 Coded Event Outcomes 1.3.6.1.4.1.19376.1.7.3.1.1.13.7**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.7.3.1.1.13.7	
<b>Parent Template</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.9	
<b>General Description</b>	The Coded Event Outcome Section shall include a narrative description of the outcomes following a procedure, an intervention or a problem, and outcomes related to the labor and delivery process such as live birth or stillborn. It shall include entries for observation as described in the Simple Observation entry.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
42545-4	R	EVENT OUTCOME
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.13	R	Simple Observation
TBD	R2	Patient Transfer

2015

2020

```

<component>
 <section>
 <!--Coded Event Outcomes assessment section templateId -->
 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.13.7' />
 <id root=' ' extension=' '/>
 <code code='42545-4' displayName='Event Outcome'
 codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" />
 <text>
 :
 </text>
 <entry>
 :
 <!-- Required Simple Observation element -->
 <templateId root=1.3.6.1.4.1.19376.1.5.3.1.4.13/>
 :
 </entry>
 </section>
</component>

```

2025

2030

2035

*Add Section 6.3.3.2.50 (Occupational History - removed 2011-09 at the request of QRPH)*

### 6.3.3.2.50 Intentionally blank

*Add Section 6.3.3.2.51 (Patient Status - removed 2011-09 at the request of QRPH)*

### 6.3.3.2.51 Intentionally blank

*Add Section 6.3.3.2.52 Cancer Control - removed 2011-09 at the request of QRPH)*

2040

### 6.3.3.2.52 Intentionally blank

*Add Section 6.3.3.2.53*

### 6.3.3.2.53 Notifications, Alerts, and Reminders Section

**1.3.6.1.4.1.19376.1.5.3.1.1.20.2.1.x**

<b>Template ID</b>	<b>1.3.6.1.4.1.19376.1.5.3.1.1.20.2.1.x</b>	
<b>General Description</b>	The Notifications, Reminders and Alerts section highlights areas of care plan non-conformance and directs the need for follow-up communications.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
XXX	R	Notifications, Alerts, and Reminders
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.7	C	<u>Medications</u> Medications entries shall appear for all pending medications when present. These entries shall be in intent mood.
1.3.6.1.4.1.19376.1.5.3.1.4.19	C	<u>Procedure</u> Procedure entries shall appear for all pending procedures when present. These entries shall be in intent mood.

1.3.6.1.4.1.19376.1.5.3.1.4.14	C	<u>Encounter</u> Encounter entries should appear for all pending follow-up encounters. These entries shall be in promise or appointment request mood.
1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1	C	<u>Observation Request</u> Observation request entries should appear for all pending follow-up observations. These entries shall appear in intent mood.

2045

**Add Section 6.3.3.2.54****6.3.3.2.54 Pain Assessment Panel Section 1.3.6.1.4.1.19376.1.5.3.1.1.20.2.4**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.20.2.4	
<b>General Description</b>	This contains a narrative description of the patient's pain, including such items as severity, quality, location, time of onset, radiation, etc.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
38212-7	R	Pain Assessment Panel

2050

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.20.2.4' />
 <id root=' ' extension=' ' />
 <code code='38212-7' displayName='Pain Assessment Panel'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>

```

2055

2060

**Figure 6.3.2.54-1: Specification for Pain Assessment Panel Section****Add Section 6.3.3.2.55****6.3.3.2.55 History of Cognitive Function Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.11**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.11	
<b>General Description</b>	This contains a narrative description of a patient's mental status.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
11332-4	R	History of Cognitive Function

2065

2070

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.11'>
 <id root=' ' extension=' '/>
 <code code='11332-4' displayName='History of Cognitive Function'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'>
 <text>
 Text as described above
 </text>
 </code>
 </section>
 </component>

```

2075

**Figure 6.3.2.55-1: Specification for History of Cognitive Function Section**

2080

*Add Section 6.3.3.2.56 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)*

**6.3.3.2.56 Isolation Status Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.8**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.8	
<b>General Description</b>	The Isolation Status section describes a patient with an active infectious disease requiring additional personal protective equipment for healthcare providers.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
55017-8	R	ISOLATION OR QUARANTINE FOR ACTIVE INFECTIOUS DISEASE

2085

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.8'>
 <id root=' ' extension=' '/>
 <code code='55017-8' displayName=' ISOLATION OR QUARANTINE FOR ACTIVE INFECTIOUS DISEASE '
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'>
 <text>
 Text as described above
 </text>
 </code>
 </section>
 </component>

```

2090

2095

**Figure 6.3.3.2.56-1: Sample Isolation Status Section**

*Add Section 6.3.3.2.57 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)*

**6.3.3.2.57 Restraints Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.10**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.10	
<b>General Description</b>	The Restraints section describes the type of restraints currently in use on the patient to be transported.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
46067-5	R	DEVICES AND RESTRAINTS SET

2100

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.10' />
 <id root=' ' extension=' ' />
 <code code='46067-5' displayName='DEVICES AND RESTRAINTS SET'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>

```

2105

2110

**Figure 6.3.3.2.57-1: Sample Restraints Section**

Add Section 6.3.3.2.58. Added 2011-09 from QRPH EHCP profile

### 6.3.3.2.58 Risk Indicators for Hearing Loss

<b>Template ID</b>	1.3.6.1.4.1.19376.1.7.3.1.1.15.3.1	
<b>General Description</b>	This section SHALL include at least one entry describing hearing risk indicators for the subject	
<b>LOINC® Code</b>	<b>Opt</b>	<b>Description</b>
58232-0	R	HEARING LOSS RISK INDICATOR
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.7.3.1.1.15.5.1	R	Risk Indicators for Hearing Loss Entry

2115

```

<component>
 <section>
 <templateId root=1.3.6.1.4.1.19376.1.7.3.1.1.15.3.1' />
 <id root=' ' extension=' ' />
 <code code='58232-0' displayName='HEARING LOSS RISK INDICATOR '
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Required Risk Indicators for Hearing Loss Entry element -->
 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.15.5.1' />
 :
 </entry>
 </section>
</component>

```

2120

2125

2130

**Figure 6.3.3.2.58-1: Sample Coded Risk Indicators for Hearing Loss Section**

2135

Add Section 6.3.3.2.58. Added 2011-09 from QRPH PRPH-Ca profile

### 6.3.3.2.59 Cancer Diagnosis Section 1.3.6.1.4.1.19376.1.7.3.1.3.14.1

<b>Template ID</b>	1.3.6.1.4.1.19376.1.7.3.1.3.14.1
<b>Parent ID</b>	PCC Active Problem Section 1.3.6.1.4.1.19376.1.5.3.1.3.6 CCD 3.5 2.16.840.1.113883.10.20.1.11

<b>General Description</b>	This section contains specific detailed information about cancer diagnosis(es) that are currently being monitored for the patient. A separate entry for each cancer diagnosis SHALL be provided.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
xxxxx-x	R	Cancer Diagnosis
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.7.3.1.4.14.1	R	Cancer Diagnosis Entry

```

2140 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.3.14.1'>
 <id root=' ' extension=' '/>
 <code code='xxxxx-x' displayName='Cancer Diagnosis'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 <!-- Required Cancer Diagnosis Entry element -->
 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.4.14.1'>
 :
 </entry>
 </section>
 </component>

```

**Figure 6.3.2.59-1: Specification for Cancer Diagnosis Section**

2155

*Add Section 6.3.3.3***6.3.3.3 Medications***Add Section 6.3.3.3.1***6.3.3.3.1 Medications Section**

2160

*Add Section 6.3.3.3.2***6.3.3.3.2 Admission Medication History Section***Add Section 6.3.3.3.3***6.3.3.3.3 Medications Administered Section***Add Section 6.3.3.3.4*

2165

**6.3.3.3.4 Hospital Discharge Medications Section***Add Section 6.3.3.3.5*

### **6.3.3.3.5 Immunizations Section**

*Add Section 6.3.3.4*

### **6.3.3.4 Physical Exams**

2170

*Add Section 6.3.3.4.1*

#### **6.3.3.4.1 Physical Examination Section**

*Add Section 6.3.3.4.2*

#### **6.3.3.4.2 Detailed Physical Examination Section**

*Add Section 6.3.3.4.3*

2175

#### **6.3.3.4.3 Hospital Discharge Physical Exam Section**

*Add Section 6.3.3.4.4*

#### **6.3.3.4.4 Vital Signs Section**

*Add Section 6.3.3.4.5*

#### **6.3.3.4.5 Coded Vital Signs Section**

2180

*Add Section 6.3.3.4.29*

#### **6.3.3.4.29 Extremities**

*Add Section 6.3.3.4.30*

#### **6.3.3.4.30 Coded Detailed Physical Examination Section**

**1.3.6.1.4.1.19376.1.5.3.1.1.9.15.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.15.1	
<b>Parent Template</b>	Detailed Physical Examination (1.3.6.1.4.1.19376.1.5.3.1.1.9.15)	
<b>General Description</b>	The Coded Detailed Physical Examination section shall contain a narrative description of the patient's physical findings. It shall include subsections, if known, for the exams that are performed.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
29545-1	R	PHYSICAL EXAMINATION
<b>Subsections</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.5.3.2	R2	Coded Vital Signs Vital signs may be a subsection of the physical examination or they may stand alone
1.3.6.1.4.1.19376.1.5.3.1.1.9.16	R2	General Appearance

1.3.6.1.4.1.19376.1.5.3.1.1.9.48	R2	Visible Implanted Medical Devices
1.3.6.1.4.1.19376.1.5.3.1.1.9.17	R2	Integumentary System
1.3.6.1.4.1.19376.1.5.3.1.1.9.18	R2	Head
1.3.6.1.4.1.19376.1.5.3.1.1.9.19	R2	Eyes
1.3.6.1.4.1.19376.1.5.3.1.1.9.20	R2	Ears, Nose, Mouth and Throat
1.3.6.1.4.1.19376.1.5.3.1.1.9.21	R2	Ears
1.3.6.1.4.1.19376.1.5.3.1.1.9.22	R2	Nose
1.3.6.1.4.1.19376.1.5.3.1.1.9.23	R2	Mouth, Throat, and Teeth
1.3.6.1.4.1.19376.1.5.3.1.1.9.24	R2	Neck
1.3.6.1.4.1.19376.1.5.3.1.1.9.25	R2	Endocrine System
1.3.6.1.4.1.19376.1.5.3.1.1.9.26	R2	Thorax and Lungs
1.3.6.1.4.1.19376.1.5.3.1.1.9.27	R2	Chest Wall
1.3.6.1.4.1.19376.1.5.3.1.1.9.28	R2	Breasts
1.3.6.1.4.1.19376.1.5.3.1.1.9.29	R2	Heart
1.3.6.1.4.1.19376.1.5.3.1.1.9.30	R2	Respiratory System
1.3.6.1.4.1.19376.1.5.3.1.1.9.31	R2	Abdomen
1.3.6.1.4.1.19376.1.5.3.1.1.9.32	R2	Lymphatic System
1.3.6.1.4.1.19376.1.5.3.1.1.9.33	R2	Vessels
1.3.6.1.4.1.19376.1.5.3.1.1.9.34	R2	Musculoskeletal System
1.3.6.1.4.1.19376.1.5.3.1.1.9.35	R2	Neurologic System
1.3.6.1.4.1.19376.1.5.3.1.1.9.36	R2	Genitalia
1.3.6.1.4.1.19376.1.5.3.1.1.9.37	R2	Rectum
1.3.6.1.4.1.19376.1.5.3.1.1.16.2.1	R2	Extremities
1.3.6.1.4.1.19376.1.5.3.1.1.21.2.1 0	R2	Pelvis

2185

```

2190 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.15' />
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.15.1' />
<id root=' ' extension=' ' />
<code code='29545-1' displayName='PHYSICAL EXAMINATION'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.25' />
 <!-- Optional Vital Signs Section content -->
 </section>
</component>
</section>
</component>

```

2205 **Figure 6.3.3.4.30-1: Coded Detailed Physical Examination Section**

<i>Add Section 6.3.3.4.31</i>
-------------------------------

### 6.3.3.4.31 Pelvis Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.10

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.10	
<b>General Description</b>	The Pelvis section shall include a narrative description of any type of exam of the reproductive organs.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
10204-6	R	PELVIS
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.5	O	Problem Entry

```

2210 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.10' />
<id root=' ' extension=' ' />
<code code='10204-6' displayName='PELVIS'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
<!-- Optional Problem Entry element -->
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5' />
 :
</entry>
</section>
</component>

```

2225 **Figure 6.3.3.4.31-1: Pelvis Section**

2230

*Add Section 6.3.3.4.32***6.3.3.4.32 Admission Physical Exam Section 1.3.6.1.4.1.19376.1.5.3.1.1.22.1.1.2.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.22.1.1.2.1	
<b>General Description</b>	The Admission physical exam section shall include a narrative description of the physical exams given during the admission to a hospital or similar type of facility.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
XX-AdmissionPhysicalExam	R	Admission physical exam

2235

```
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.22.1.1.2.1' />
 <id root=' ' extension=' ' />
 <code code='XX-AdmissionPhysicalExam' displayName='Admission physical exam'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

2240

**Figure 6.3.3.4.32-1: Admission Physical Exam Section**

2245

*Add Section 6.3.3.4.33***6.3.3.4.33 Discharge Status 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.12**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.12	
<b>Parent Template</b>		
<b>General Description</b>	Discharge status should contain a narrative description of the status/condition of the patient at the time of discharge, such as stable, critical, etc.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
52523-8	R2	Discharge Status

2250

```
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.12' />
 <id root=' ' extension=' ' />
 <code code='52523-8' displayName='Discharge status'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

2255

**Figure 6.3.3.4.33-1: Discharge Status Section**

2260

**6.3.3.5 Relevant Studies***Add Section 6.3.3.5.1***6.3.3.5.1 Results***Add Section 6.3.3.5.2*2265 **6.3.3.5.2 Coded Results***Add Section 6.3.3.5.3***6.3.3.5.3 Hospital Studies Summary***Add Section 6.3.3.5.4***6.3.3.5.4 Coded Hospital Studies Summary**2270 *Add Section 6.3.3.5.5***6.3.3.5.5 Consultations 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.8**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.8	
<b>General Description</b>	The ED Consultations section shall contain a narrative description of the consultations obtained during an encounter of care.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
18693-2	R	ED CONSULTANT PRACTITIONER

2275

```
<component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.8'>
<id root=' ' extension=' ' />
<code code='18693-2' displayName='ED CONSULTANT PRACTITIONER'
codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
Text as described above
</text>
</section>
</component>
```

2280

**Figure 6.3.3.5.5-1: Specification for ED Consultations Section**

2285

*Add Section 6.3.3.5.6***6.3.3.5.6 Antenatal Testing and Surveillance Section****1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5
<b>Parent Template</b>	

<b>General Description</b>	The Antenatal Testing and Surveillance section shall contain a narrative description of reports and data from tests and surveillance performed during the pregnancy (e.g. Ultrasound, Biophysical Profile, Non-Stress Test, Contraction Stress Test)	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57078-8	R	Antenatal testing and surveillance

2290      <component>  
           <section>

      <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5' />  
       <id root=' ' extension=' ' />  
       <code code='57078-8' displayName='ANTENATAL TESTING AND SURVEILLANCE'  
           codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />  
       <text>  
           Text as described above  
       </text>  
   </component>

2295

**Figure 6.3.3.5.6-1: Specification for and Surveillance Section**

*Add Section 6.3.3.5.7*

### **6.3.3.5.7 Coded Antenatal Testing and Surveillance Section**

**1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5.1	
<b>Parent Template</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5	
<b>General Description</b>	The Antenatal Testing and Surveillance section shall contain a narrative and coded description of reports and data from tests and surveillance performed during the pregnancy (e.g. Ultrasound, Biophysical Profile, Non-Stress Test, Contraction Stress Test). It shall contain an Antenatal Testing and Surveillance Battery.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57078-8	R	ANTENATAL TESTING AND SURVEILLANCE
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.21.3.10	R	Antenatal Testing and Surveillance Battery

2305

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.5.1' />
 <id root=' ' extension=' ' />
 <code code='57078-8' displayName='ANTENATAL TESTING AND SURVEILLANCE'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Required Antenatal Testing and Surveillance Battery -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.3.10' />
 :
 </entry>
 </component>

```

**Figure 6.3.3.5.7-1: Specification for Coded Antenatal Testing and Surveillance Section**

2325 *Add Section 6.3.3.5.8 (Diagnosis - Removed 2011-09 at the request of QRPH)*

**6.3.3.5.8 Intentionally blank**

*Add Section 6.3.3.5.9 (TNM Stage – removed 2011-09 at the request of QRPH)*

**6.3.3.5.9 Intentionally blank**

2330 *Add Section 6.3.3.5.10 (Cancer Supporting Documentation - removed 2011-09 at the request of QRPH)*

**6.3.3.5.10 Intentionally blank**

*Add Section 6.3.3.5.11. (Added 2011-09 from QRPH EHCP profile)*

**2335 6.3.3.5.11 Hearing Screening Coded Results**

The Hearing Screening Coded Results section SHALL contain the hearing screening results of pass or refer for the right ear and pass or refer for the left ear, expressed as LOINC® codes as well as the coded methodology to complete the screening. Coded methodology includes (LOINC 54106-0) Automated Auditory Brainstem Response, Auditory Brainstem Response, Otoacoustic Emissions, Transient Otoacoustic Emissions, and Distortion Product Otoacoustic Emissions. If the methodology is unknown, the coded result of unknown method SHALL be used. Where the screening results are not available, the reason the results are not available SHALL be present. This could include unsuccessful, technical fail; not performed , not performed, medical exclusion. The Hearing Screening Coded Results section is required.

2345

<b>Template ID</b>	1.3.6.1.4.1.19376.1.7.3.1.1.15.3.2
--------------------	------------------------------------

<b>Parent Template</b>	Coded Results (1.3.6.1.4.1.19376.1.5.3.1.3.28)	
<b>General Description</b>	The Hearing Screening Code Results section SHALL include at least one observation entry describing the hearing screening results as described in the Entry Content Module. Where there are no hearing screening results performed, then the reason SHALL be indicated	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
30954-2	R	Relevant diagnostic tests/laboratory data
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.19	R	Procedure Entry
1.3.6.1.4.1.19376.1.5.3.1.4.4	R2	References Entry
1.3.6.1.4.1.19376.1.5.3.1.4.13	R	Simple Observation

### 6.3.3.5.11.1 Parent Template

The parent of this template is Coded Results.

```

2350 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.28' />
 <id root=' ' extension=' ' />
 <code code='30954-2' displayName='Relevant diagnostic tests/laboratory data'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Required Procedure Entry element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.19' />
 :
 </entry>
 <entry>
 :
 <!-- Required if known References Entry element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.4' />
 :
 </entry>
 <entry>
 :
 <!-- Optional Simple Observation element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
 </entry>
 </section>

```

**Figure 6.3.3.5.11-1: Hearing Screening Coded Results Section**

2380    Add Section 6.3.3.6

### 6.3.3.6 Plans of Care

*Add Section 6.3.3.6.1*

#### **6.3.3.6.1 Care Plan**

2385

*Add Section 6.3.3.6.2*

#### **6.3.3.6.2 Assessment and Plan**

*Add Section 6.3.3.6.3*

#### **6.3.3.6.3 Discharge Disposition**

*Add Section 6.3.3.6.4*

2390

#### **6.3.3.6.4 Discharge Diet**

*Add Section 6.3.3.6.5*

#### **6.3.3.6.5 Advance Directives**

*Add Section 6.3.3.6.6*

#### **6.3.3.6.6 Coded Advance Directives**

2395

*Add Section 6.3.3.6.7*

#### **6.3.3.6.7 Transport Mode**

*Add Section 6.3.3.6.8*

#### **6.3.3.6.8 Procedure Care Plan Status Report Section**

2400

##### **1.3.6.1.4.1.19376.1.5.3.1.1.9.45**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.45	
<b>Parent Template</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.40	
<b>General Description</b>	The procedure care plan status report section shall contain a description of the progress towards completing expectations for care including actions completed in fulfillment of proposals, goals, and order requests for monitoring, tracking, or improving the condition of the patient prior to the procedure.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
18776-5	R	TREATMENT PLAN

### Sample Procedure Care Plan Status Report Section

```

2405 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.40' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.45' />
 <id root=' ' extension=' ' />
 <code code='18776-5' displayName='TREATMENT PLAN'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>

```

2420     Add Section 6.3.3.6.9

### 6.3.3.6.9 Health Maintenance Care Plan Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.50

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.50	
<b>Parent Template</b>	1.3.6.1.4.1.19376.1.5.3.1.3.31	
<b>General Description</b>	The health maintenance care plan section shall contain a description of the expectations for wellness care including proposals, goals, and order requests for monitoring, tracking, or improving the lifetime condition of the patient with goals of educating the patient on how to reduce the modifiable risks of the patient's genetic, behavioral, and environmental pre-conditions and otherwise optimizing lifetime outcomes.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
18776-5	R	TREATMENT PLAN

### Sample Health Maintenance Care Plan Section

```

2425 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.31' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.50' />
 <id root=' ' extension=' ' />
 <code code='18776-5' displayName='TREATMENT PLAN'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>

```

2440     Add Section 6.3.3.6.10

### 6.3.3.6.10 Health Maintenance Care Plan Status Report Section

#### 1.3.6.1.4.1.19376.1.5.3.1.1.9.41

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.41	
<b>Parent Template</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.50	
<b>General Description</b>	The health maintenance status report section shall contain a description of the progress towards completing expectations for care including actions completed in fulfillment of proposals, goals, and order requests for monitoring, tracking, or improving the condition of the patient.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
18776-5	R	TREATMENT PLAN

### Sample Health Maintenance Care Plan Status Report Section

2445

```
<component>
 <section> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.50' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.41' />
 <id root=' ' extension=' '/>
 <code code='18776-5' displayName='TREATMENT PLAN'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

2450

2455

2460

Add Section 6.3.3.6.11
------------------------

### 6.3.3.6.11 Provider Orders Section 1.3.6.1.4.1.19376.1.5.3.1.1.20.2.1

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.20.2.1	
<b>General Description</b>	The provider orders shall contain a list of all pertinent orders from healthcare providers.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
46209-3	R	PROVIDER ORDERS
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.7	C	<u>Medications</u> Medications entries shall appear for all ordered medications when present. These entries shall be in intent mood.
1.3.6.1.4.1.19376.1.5.3.1.4.19	C	<u>Procedure</u> Procedure entries shall appear for all ordered procedures when present. These entries shall be in intent mood.
1.3.6.1.4.1.19376.1.5.3.1.4.14	O	<u>Encounter</u> Encounter entries should appear for all ordered encounters. These entries shall be in promise or appointment request mood.

1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1	O	<a href="#">Observation Requests</a> Observation request entries should appear for all ordered observations. These entries shall appear in intent mood.
------------------------------------	---	------------------------------------------------------------------------------------------------------------------------------------------------------------

**Sample Provider Orders Section**

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.20.2.1' />
 <id root=' ' extension=' ' />
 <code code='46209-3' displayName='PROVIDER ORDERS'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Required if known Medications element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.7' />
 :
 </entry>
 <entry>
 :
 <!-- Required if known Procedure element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.19' />
 :
 </entry>
 <entry>
 :
 <!-- Optional Encounter element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.14' />
 :
 </entry>
 <entry>
 :
 <!-- Optional Observation Requests element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1' />
 :
 </entry>
 </section>
</component>

```

2500      Add Section 6.3.3.6.12

**6.3.3.6.12 Birth Plan Section 1.3.6.1.4.1.19376.1.5.3.1.1.21.2.1**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.21.2.1	
<b>Parent Template</b>		
<b>General Description</b>	The Birth Plan section shall contain a narrative description of the patient's requests and expectations with respect to care she is expecting during the labor and delivery process.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57079-6	R	Birth plan

2505	<pre>&lt;component&gt;   &lt;section&gt;     &lt;templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.2.1' /&gt;     &lt;id root=' ' extension=' ' /&gt;     &lt;code code='57079-6' displayName='Birth plan'       codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' /&gt;     &lt;text&gt;       Text as described above     &lt;/text&gt;   &lt;/component&gt;</pre>
2510	

**Figure 6.3.3.6.12-1: Specification for Birth Plan Section**

2515	Add Section 6.3.3.6.13
------	------------------------

**6.3.3.6.13 Immunization Recommendations 1.3.6.1.4.1.19376.1.5.3.1.1.18.3.1**

Template Id	1.3.6.1.4.1.19376.1.5.3.1.1.18.3.1	
General Description	The Immunization Recommendation section shall be present to document the recommended vaccinations for the patient. It shall include Immunization entries in proposal mood describing the immunization plan to be developed. It may include a reference to a specific guideline in definition mood to indicate the guideline that should be conformed to, and may also include references to patient education information.	
LOINC Code	Opt	Description
18776-5	R	TREATMENT PLAN
Entries	Opt	Description
1.3.6.1.4.1.19376.1.5.3.1.4.12.2	R	Immunization Recommendation Entry At least one Immunization Plan Entry shall be present in Proposal mood to indicate what the proposed care is for the patient. Other Immunization Plan entries may appear in intent mood to indicate the current plan.

Add Section 6.3.3.6.14
------------------------

**6.3.3.6.14 Patient Education Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.38**

Template Id	1.3.6.1.4.1.19376.1.5.3.1.1.9.38	
General Description	The patient education section shall contain a description of the patient education the patient received as well as the results of the education.	
LOINC Code	Opt	Description
34895-3	R	EDUCATION NOTE
Entries	Opt	Description
1.3.6.1.4.1.19376.1.5.3.1.4.12.2	R	Immunization Recommendation Entry At least one Immunization Plan Entry shall be present in Proposal mood to indicate what the proposed care is for the patient. Other Immunization Plan entries may appear in intent mood to indicate the current plan.

2520

2525	<pre>&lt;component&gt; &lt;section&gt; &lt;templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.38' /&gt; &lt;id root=' ' extension=' ' /&gt; &lt;code code='34895-3' displayName='EDUCATION NOTE'       codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' /&gt; &lt;text&gt;   Text as described above &lt;/text&gt; &lt;/section&gt; &lt;/component&gt;</pre>
2530	

**Figure 6.3.3.6.14-1: Specification for Patient Education and Consents Section**

2535	<i>Add Section 6.3.3.6.15</i>
------	-------------------------------

**6.3.3.6.15 Care Plan Section 1.3.6.1.4.1.19376.1.5.3.1.3.31**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.3.31	
<b>Parent Template</b>	2.16.840.1.113883.10.20.1.10	
<b>General Description</b>	The care plan section shall contain a narrative description of the expectations for care including proposals, goals, and order requests for monitoring, tracking, or improving the condition of the patient.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
18776-5	R	TREATMENT PLAN
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1	O	<p><a href="#">Observation Requests</a>            The care plan may include observation requests in intent, goal or proposal mood to identify intended observations that are part of the care plan, goals of the plan, or proposed observations (e.g., from clinical decision support).</p>
1.3.6.1.4.1.19376.1.5.3.1.4.7	O	<p><a href="#">Medication</a>            The care plan may include medication entries to identify those medications that are or are proposed to be part of the care plan.</p>
1.3.6.1.4.1.19376.1.5.3.1.4.12	O	<p><a href="#">Immunization</a>            The care plan may include immunization entries to identify those immunizations that are or are proposed to be part of the care plan.</p>
1.3.6.1.4.1.19376.1.5.3.1.4.19	O	<p><a href="#">Procedure</a>            The care plan may include procedure entries to identify those procedures that are or are proposed to be part of the care plan.</p>
1.3.6.1.4.1.19376.1.5.3.1.4.14	O	<p><a href="#">Encounter</a>            The care plan may include encounter entries in to identify those encounters that are or are proposed to be part of the care plan.</p>

```

<component>
 <section>
 <templateId root='2.16.840.1.113883.10.20.1.10' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.31' />
 <id root=' ' extension=' ' />
 <code code='18776-5' displayName='TREATMENT PLAN'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Optional Observation Requests element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1' />
 :
 </entry>
 <entry>
 :
 <!-- Optional Medication element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.7' />
 :
 </entry>
 <entry>
 :
 <!-- Optional Immunization element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.12' />
 :
 </entry>
 <entry>
 :
 <!-- Optional Procedure element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.19' />
 :
 </entry>
 <entry>
 :
 <!-- Optional Encounter element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.14' />
 :
 </entry>
 </section>
</component>

```

Figure 6.3.3.6.15-1: Specification for Care Plan Section

**6.3.3.6.16 Diet and Nutrition Section 1.3.6.1.4.1.19376.1.5.3.1.1.20.2.2**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.20.2.2	
<b>General Description</b>	This shall contain a narrative description of the diet restrictions necessary due to disease.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
XX-DietAndNutrition	R	Diet and nutrition

2585

```
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.20.2.2' />
 <id root=' ' extension=' ' />
 <code code='XXDiet-Restrictions' displayName='Diet and nutrition'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

2590

**Figure 6.3.3.6.16-1: Specification for Diet Restrictions Section****6.3.3.6.17 Intake and Output Section 1.3.6.1.4.1.19376.1.5.3.1.1.20.2.3**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.20.2.3	
<b>General Description</b>	This section shall contain a narrative description of specific fluid inputs or fluid outputs for the patient.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
XX-IntakeAndOutput	R	Intake and output

2600

```
<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.20.2.3' />
 <id root=' ' extension=' ' />
 <code code='XX-FluidManagement' displayName='Intake and output'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

2605

**Figure 6.3.3.6.17-1: Specification for Fluid Management Section**

2610

*Add Section 6.3.3.6.18 (Cancer Course of Treatment – Remove 2011-09 at the request of QRPH)***6.3.3.6.18 Intentionally blank***Add Section 6.3.3.6.19 (Cancer Treatment Plan – removed 2011-09 at the request of QRPH)*2615    **6.3.3.6.19 Intentionally blank***Add section 6.3.3.6.20***6.3.3.6.20 Procedure Care Plan Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.40**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.40
--------------------	----------------------------------

<b>Parent Template</b>	1.3.6.1.4.1.19376.1.5.3.1.3.31 (1.3.6.1.4.1.19376.1.5.3.1.3.31)	
<b>General Description</b>	The procedure care plan section shall contain a description of the expectations for care including proposals, goals, and order requests for monitoring, tracking, or improving the condition of the patient prior, during and after a procedure with goals of educating the patient, reducing the modifiable risks of the procedure and anesthesia and otherwise optimizing the outcomes. The care plan will often be updated immediately following the addition of new impressions during the course of pre-procedure evaluation.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
18776-5	R	TREATMENT PLAN
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1	O	<p><a href="#">Observation Requests</a>            The care plan may include observation requests in intent, goal or proposal mood to identify intended observations that are part of the care plan, goals of the plan, or proposed observations (e.g., from clinical decision support).</p>
1.3.6.1.4.1.19376.1.5.3.1.4.7	O	<p><a href="#">Medication</a>            The care plan may include medication entries to identify those medications that are or are proposed to be part of the care plan.</p>
1.3.6.1.4.1.19376.1.5.3.1.4.12	O	<p><a href="#">Immunization</a>            The care plan may include immunization entries to identify those immunizations that are or are proposed to be part of the care plan.</p>
1.3.6.1.4.1.19376.1.5.3.1.4.19	O	<p><a href="#">Procedure</a>            The care plan may include procedure entries to identify those procedures that are or are proposed to be part of the care plan.</p>
1.3.6.1.4.1.19376.1.5.3.1.4.14	O	<p><a href="#">Encounter</a>            The care plan may include encounter entries in to identify those encounters that are or are proposed to be part of the care plan.</p>

2620	<pre>&lt;component&gt;   &lt;section&gt;     &lt;templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.31' /&gt;     &lt;templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.40' /&gt;     &lt;id root=' ' extension=' '/&gt;     &lt;code code='18776-5' displayName='TREATMENT PLAN'       codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' /&gt;     &lt;text&gt;       Text as described above     &lt;/text&gt;   &lt;/section&gt; &lt;/component&gt;</pre>
2625	
2630	
2635	

**Figure 6.3.3.6.20-1: Sample Care Plan Section**

Add section 6.3.3.6.21 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)

### 6.3.3.6.21 Protocols Used Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.5

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.5
--------------------	------------------------------------

<b>General Description</b>	The Protocols Used section describes the protocol used by EMS personnel to direct the clinical care of the patient.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
52019-7	R	DESCRIPTION OF SERVICES PERFORMED TO SUPPORT LEVEL OF SERVICE

2640

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.5' />
 <id root=' ' extension=' ' />
 <code code='52019-7' displayName='DESCRIPTION OF SERVICES PERFORMED TO SUPPORT LEVEL OF SERVICE'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

2645

2650

**Figure 6.3.3.6.21-1: Sample Protocols Used Section**

2655

Add section 6.3.3.6.22 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)

**6.3.3.6.22 Invasive Airway Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.7**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.7	
<b>General Description</b>	The Invasive Airway section describes if, and what type, of advanced airway used.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
NEMSIS EProtocols.01 (1)	R	PROTOCOLS USED AIRWAY

2660

2665

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.7' />
 <id root=' ' extension=' ' />
 <code code='NEMSIS EProtocols.01(1)' displayName='PROTOCOLS USED AIRWAY'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='NEMSIS' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

**Figure 6.3.3.6.22-1: Sample Invasive Airway Section**

2670

Add section 6.3.3.6.23 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)

**6.3.3.6.23 Ventilator Usage Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.11**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.11
<b>General Description</b>	The Ventilator Usage section describes

LOINC Code	Opt	Description
20124-4	R	VENTILATION MODE

2675      <component>  
           <section>  
           <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.11' />  
           <id root=' ' extension=' ' />  
           <code code='2012404' displayName='VENTILATION MODE'  
               codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />  
           <text>  
               Text as described above  
           </text>  
           </section>  
       </component>

2680

2685      **Figure 6.3.3.6.23-1: Sample Ventilator Usage Section**

Add Section 6.3.3.7

### 6.3.3.7 Administrative and Other Information

Add Section 6.3.3.7.1

2690      **6.3.3.7.1 Payers**

Add Section 6.3.3.7.2

#### 6.3.3.7.2 Referral Source

Add Section 6.3.3.7.3

#### 6.3.3.7.3 Transport Mode

2695      Add Section 6.3.3.7.4

#### 6.3.3.7.4 ED Disposition

Add Section 6.3.3.7.5 (Cancer Payers – Removed 2011-09 at the request of QRPH)

#### 6.3.3.7.5 Intentionally blank

2700      Add Section 6.3.3.7.6 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)

#### 6.3.3.7.6 Sending Facility Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.1

Template ID	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.1
General Description	The Sending Facility section contains the name and address of the healthcare facility

	that is sending the patient for transport.	
LOINC Code	Opt	Description
52023-9	R	ORIGINATION SITE NAME AND ADDRESS

2705           <component>  
                 <section>  
                 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.1' />  
                 <id root=' ' extension=' ' />  
                 <code code='52023-9' displayName='ORIGINATION SITE NAME AND ADDRESS'  
                   codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />  
                 <text>  
                   Text as described above  
                 </text>  
                 </section>  
       </component>

2710

**Figure 6.3.3.7.6-1: Sample Sending Facility Section**

Add Section 6.3.3.7.7 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)

**6.3.3.7.7 Receiving Facility Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.2**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.2	
<b>General Description</b>	The Receiving Facility section contains the name and address of the healthcare facility that is receiving the transported patient.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
52026-2	R	DESTINATION SITE NAME & ADDRESS

2720

2725           <component>  
                 <section>  
                 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.2' />  
                 <id root=' ' extension=' ' />  
                 <code code='52026-2' displayName='DESTINATION SITE NAME & ADDRESS'  
                   codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />  
                 <text>  
                   Text as described above  
                 </text>  
                 </section>  
       </component>

2730

**Figure 6.3.3.7.7-1: Sample Receiving Facility Section**

2735 Add Section 6.3.3.7.8 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)

**6.3.3.7.8 Mass Casualty Incident Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.3**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.3
<b>General Description</b>	The Mass Casualty Incident Section indicates if this event would be considered a mass

	casualty incident overwhelming existing EMS and ED resources.	
LOINC Code	Opt	Description
NA – NEMSIS EScene.07	R2	MASS CASUALTY INCIDENT

2740

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.3' />
 <id root=' ' extension=' '/>
 <code code='NEMSIS Escene.07' displayName='MASS CASUALTY INCIDENT'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='NEMSIS' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

2745

**Figure 6.3.3.7.8-1: Sample Mass Casualty Incident Section**

2750

Add Section 6.3.3.7.9 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)

**6.3.3.7.9 Unit Response Level Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.4**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.4	
<b>General Description</b>	The Unit Response Level section describes the level of service provided for this transport.	
LOINC Code	Opt	Description
51995-9	R	RATIONALE FOR TYPE OF TRANSPORT

2755

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.4' />
 <id root=' ' extension=' '/>
 <code code='51995-9' displayName='RATIONALE FOR TYPE OF TRANSPORT'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 </section>
</component>
```

2760

2765

**Figure 6.3.3.7.9-1: Sample Unit Response Level Section**

Add Section 6.3.3.7.10 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)

2770

**6.3.3.7.10 Extra Attendants Information Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.6**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.6
<b>General Description</b>	The Protocols Used section describes the protocol used by EMS personnel to direct the

	clinical care of the patient.	
LOINC Code	Opt	Description
52074-2	R2	EXTRA ATTENDANTS INFORMATION

2775	<component> <section> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.6'> <id root=' ' extension=' '/> <code code='52074-2' displayName='EXTRA ATTENDANTS INFORMATION' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' /> <text> Text as described above </text> </section> </component>
2780	

**Figure 6.3.3.7.10-1: Sample Extra Attendants Information Section**

2785     *Add Section 6.3.3.7.11 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)*

### 6.3.3.7.11 Provider Level Section 1.3.6.1.4.1.19376.1.5.3.1.1.25.2.9

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.25.2.9	
<b>General Description</b>	The Provider Level section describes the certification or licensure level of the healthcare provider.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
NEMSIS DConfiguration.02	R	STATE CERTIFICATION LICENSURE LEVELS

2790	<component> <section> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.25.2.9'> <id root=' ' extension=' '/> <code code='NEMSIS DConfiguration.02' displayName='STATE CERTIFICATION LICENSURE LEVELS' codeSystem='2.16.840.1.113883.6.1' codeSystemName='NEMSIS' /> <text> Text as described above </text> </section> </component>
2795	

**Figure 6.3.3.7.11-1: Sample Provider Level Section**

*Add Section 6.3.3.8*

### 6.3.3.8 Interventions

2805     This section contains section content modules that describe interventions, procedures, therapeutic treatments, et cetera, performed on the patient.

Add Section 6.3.3.8.3
-----------------------

### 6.3.3.8.3 Procedures and Interventions Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11	
<b>General Description</b>	The Procedures and Interventions section shall contain a narrative description of the actions performed by a clinician.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
29544-3	R	PROCEDURES
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.19	R	Procedure This entry provides coded values for procedures performed during the encounter.

2810	<component> <section> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.11'> <id root=' ' extension=' ' /> <code code='X-PROC' displayName='PROCEDURES PERFORMED' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' /> <text> Text as described above </text> <entry> : <!-- Required Procedure element --> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.19'> : </entry> </section> </component>
2815	
2820	
2825	

**Figure 6.3.3.8.3-1: Specification for Procedures and Interventions Section**

2830	Add Section 6.3.3.8.4 (Added 2011-06 from the PCC Transport Record Summary Profiles supplement)
------	-------------------------------------------------------------------------------------------------

### 6.3.3.8.4 Intravenous Fluids Administered Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.6

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.6	
<b>General Description</b>	The intravenous fluids administered section shall contain a narrative description of fluids administered to a patient during the course of an encounter. It may include entries for IV fluid administration as described in the Entry Content Module.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57072-1	R	Intravenous fluids administered
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.13.3.2	R	<a href="#">Intravenous Fluids</a>

```

2835 <component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.6' />
 <id root=' ' extension=' ' />
 <code code='57072-1' displayName='Intravenous fluids administered'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Required Intravenous Fluids element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.2' />
 :
 </entry>
 </section>
</component>

```

**Figure 6.3.3.8.4-1: Specification for Intravenous Fluids Administered Section***Add Section 6.3.3.9***2855 6.3.3.9 Impressions**

This section contains section content modules that describe assessments, impressions, diagnoses, or other reporting of clinical opinions or judgment.

*Add Section 6.3.3.9.1***2860 6.3.3.9.1 Pre-procedure Impressions Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.42  
(Deprecated)***Add Section 6.3.3.9.2***6.3.3.9.2 Pre-procedure Risk Assessment Section 1.3.6.1.4.1.19376.1.5.3.1.1.9.44**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.9.44	
<b>General Description</b>	The pre-procedure risk section shall contain a description of the risks the patient faces because of the planned procedure and associated anesthesia, especially in the context of modifiable risks identified by patient findings. It shall include entries for patient risks as described in the Entry Content Module.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
11450-4	R	PROBLEM LIST

2865

The parent of this template is 1.3.6.1.4.1.19376.1.5.3.1.3.6

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.3.6' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.9.44' />
 <id root=' ' extension=' ' />
 <code code='11450-4' displayName='PROBLEM LIST'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Required Problem Entry element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5' />
 :
 </entry>
 </section>
</component>

```

**Figure 6.3.3.9.2-1: Specification for Pre-procedure Risk Assessment Section**

2890     Add Section 6.3.3.9.3

### 6.3.3.9.3 Antepartum Visit Summary Flowsheet Section

#### 1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.2

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.2	
<b>General Description</b>	This section is a running history of the most important elements noted for a pregnant woman.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
57059-8	R	Pregnancy visit summary
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.13	R	<a href="#">Simple Observation</a> The flowsheet contains one simple observation to represent the Prepregnancy Weight. This observation SHALL be valued with the LOINC code 8348-5, BODY WEIGHT^PRE PREGNANCY-MASS-PT-QN-MEASURED. The value SHALL be of type PQ. The units may be either "lb_av" or "kg".
1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.2	R	<a href="#">Antepartum Flowsheet Panel</a> Other entries on the flowsheet are "batteries" which represent a single visit.

```

2895 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.2.2' />
<id root=' ' extension=' ' />
<code code='57059-8' displayName='Pregnancy visit summary'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
 <!-- Required Simple Observation element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
</entry>
<entry>
 :
 <!-- Required Antepartum Flowsheet Panel element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.2' />
 :
</entry>
</section>
</component>

```

**Figure 6.3.3.9.3-1: Specification for Antepartum Visit Summary Flowsheet Section**

2920

Add Section 6.3.3.9.4
-----------------------

**6.3.3.9.4 Progress Note Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.7**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.7	
<b>General Description</b>	The Progress Note section shall contain a narrative description of the sequence of events from initial assessment to discharge for an encounter.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
18733-6	R	SUBSEQUENT EVALUATION NOTE (ATTENDING PHYSICIAN)

```

2925 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.7' />
<id root=' ' extension=' ' />
<code code='18733-6' displayName='SUBSEQUENT EVALUATION NOTE (ATTENDING PHYSICIAN)'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
</section>
</component>

```

**Figure 6.3.3.9.4-1: Specification for Progress Note Section**

2940

*Add Section 6.3.3.9.5***6.3.3.9.5 ED Diagnosis Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.9**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.9	
<b>General Description</b>	The ED diagnosis section shall contain a narrative description of the conditions that were diagnosed or addressed during the ED course, as well as those active conditions that modify the complexity of the patient encounter. It should include entries for patient conditions as described in the Entry Content Module.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
11301-9	R	ED DIAGNOSIS
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.4.5	R	<a href="#">Problem Entry</a>

2945

```

<component>
 <section>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.9' />
 <id root=' ' extension=' ' />
 <code code='11301-9' displayName='ED DIAGNOSIS'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text>
 Text as described above
 </text>
 <entry>
 :
 <!-- Required Problem Entry element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5' />
 :
 </entry>
 </section>
</component>

```

2950

2955

2960

**Figure 6.3.3.9.5-1: Specification for ED Diagnosis Section***Add Section 6.3.3.9.6***6.3.3.9.6 Acuity Assessment Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.2**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.2	
<b>General Description</b>	The Acuity Assessment section contains a description of the acuity of the patient upon presentation to the Emergency department.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
11283-9	R	ACUITY ASSESSMENT
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.13.3.1	R	<a href="#">Acuity</a> This entry provides coded values giving the triage acuity.

2965

```

2970 <component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.2' />
<id root=' ' extension=' ' />
<code code='11283-9' displayName='ACUITY ASSESSMENT'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
 <!-- Required Acuity element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.1' />
 :
</entry>
</section>
</component>

```

**Figure 6.3.3.9.6-1: Specification for Acuity Assessment Section**

2985

Add Section 6.3.3.9.7
-----------------------

**6.3.3.9.7 Assessments Section 1.3.6.1.4.1.19376.1.5.3.1.1.13.2.4**

<b>Template ID</b>	1.3.6.1.4.1.19376.1.5.3.1.1.13.2.4	
<b>General Description</b>	The assessments section contains narrative assessments of the patient status.	
<b>LOINC Code</b>	<b>Opt</b>	<b>Description</b>
51848-0	R	ASSESSMENT
<b>Entries</b>	<b>Opt</b>	<b>Description</b>
1.3.6.1.4.1.19376.1.5.3.1.1.13.3.4	O	<a href="#">Nursing Assessments Battery</a>

2990

```

<component>
<section>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.2.4' />
<id root=' ' extension=' ' />
<code code='51848-0' displayName='ASSESSMENT'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text>
 Text as described above
</text>
<entry>
 :
 <!-- Optional Nursing Assessments Battery element -->
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.4' />
 :
</entry>
</section>
</component>

```

**Figure 6.3.3.9.7-1: Specification for Assessments Section**

### 6.3.4 CDA Entry Content Modules

3010

Add Section 6.3.4.25

#### 6.3.4.25 Family History Observation 1.3.6.1.4.19376.1.5.3.1.4.13.3

3015 A family history observation is a [Simple Observation](#) that uses a specific vocabulary, and inherits constraints from CCD. Family history observations are found inside [Family History Organizers](#).

##### 6.3.4.25.1 Standards

CCD [ASTM/HL7 Continuity of Care Document](#)

##### 6.3.4.25.2 Parent Template

The parent of this template is [Simple Observation](#). This template is compatible with the ASTM/HL7 Continuity of Care Document template: 2.16.840.1.113883.10.20.1.22

3020

##### 6.3.4.25.3 Specification

```

<observation typeCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 <templateId root='2.16.840.1.113883.10.20.1.22' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13.3' />
 <id root=' ' extension=' ' />
 <code code=' ' displayName=' ' codeSystem=' ' codeSystemName=' ' />
 <text><reference value='#xxx' /></text>
 <statusCode code='completed' />
 <effectiveTime value=' ' />
 <repeatNumber value=' ' />
 <value xsi:type='CD' .../>
 <interpretationCode code=' ' codeSystem=' ' codeSystemName=' ' />
 <methodCode code=' ' codeSystem=' ' codeSystemName=' ' />
 <targetSiteCode code=' ' codeSystem=' ' codeSystemName=' ' />
</observation>
```

**Figure 6.3.4.25.3-1: Family History Specification**

#### 6.3.4.25.4 <templateId root='2.16.840.1.113883.10.20.1.22'> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13.3'>

3040 The <templateId> elements identify this observation as a family history observation, and shall be present as shown above.

#### 6.3.4.25.5 <code code=' ' displayName=' ' codeSystem=' ' codeSystemName=' ' />

The <code> indicates the type of observation made (e.g., Diagnosis, et cetera). See the code element in the Problem Entry entry for suggested values.

3045 **6.3.4.25.6 <value xsi:type='CD' code=' ' displayName=' ' codeSystem=' ' codeSystemName=' '/>**

The <value> element indicates the information (e.g., diagnosis) of the family member. See the value element in the Problem Entry for suggested values.

*Add Section 6.3.4.26*

3050

*Add Section 6.3.4.27*

**6.3.4.27 EDD Observation 1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.1**

3055 The EDD observation reflects the clinician's best judgment about the estimated delivery date of the patient. It can be supported by patient history (e.g., last menses or quickening), physical examination findings (uterine size), or Ultrasound. The observation is a Simple Observation with a supporting entryRelation of another Observation. The supporting observation may in turn have an entryRelation that gives the original observation as a gestational age or date from which the estimated due date is calculated.

### 6.3.4.27.1 Specification

```

3060 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.1' />
 <statusCode code='completed' />
 <effectiveTime value=' ' />
 <author typeCode='AUT'>
 <time value=' ' />
 <assignedAuthor>
 <id root=' ' extension=' ' />
 </assignedAuthor>
 </author>
 <id root=' ' extension=' ' />
 <code code='11778-8' displayName='DELIVERY DATE-TMSTP-PT-^PATIENT-QN-CLINICAL.ESTIMATED' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<3075 text><reference value='id-foo' /></text>
 <value xsi:type='TS' value=' ' />
 <entryRelationship typeCode='SPRT'>
 <observation classCode='OBS' moodCode='EVN'>
 <id root=' ' extension=' ' />
 <statusCode code='completed' />
 <effectiveTime value=' ' />
 <author typeCode='AUT'>
 <time value=' ' />
 <assignedAuthor classCode=' ' />
 <id root=' ' extension=' ' />
 </assignedAuthor>
 </author>
 <code code='[11779-6] (xx-EDD-by-PE) | 11781-2 | (xx-EDD-by-Qck) | (xx-EDD-by-Fund)' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<3090 value type='TS' value=' ' />
 <entryRelationship typeCode='DRIV'>
 <observation classCode='OBS' moodCode='EVN'>
 <id root=' ' extension=' ' />
 <statusCode code='completed' />
 <effectiveTime value=' ' />
 <author typeCode='AUT'>
 <time value=' ' />
 <assignedAuthor>
 <id root=' ' extension=' ' />
 </assignedAuthor>
 </author>
 <informant typeCode='INF'>
 <relatedEntity classCode=' ' />
 <id root=' ' extension=' ' />
 </relatedEntity>
 </informant>
 <code code='[8655-2] (xx-ga-by-pe) | 11888-5 | (xx-date-of-qck) | (xx-date-of-fund-umb)' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<3095 <value type='[PQ|TS]' value=' ' units='week' />
 </observation>
 </entryRelationship>
 </observation>
</entryRelationship>
</observation>

```

3115

### 6.3.4.27.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.1' />

The <templateId> identifies the observation as a type of Estimated Delivery Date Observation. The root attribute SHALL be valued with '1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.1'.

**6.3.4.27.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13'>**

- 3120 EDD observation SHALL comply with the restrictions of the Simple Observation entry. The observation SHALL NOT include repeatNumber, interpretationCode, methodCode, or targetSiteCode as listed below.

**6.3.4.27.4 <code code='11778-8' codeSystem='2.16.840.1.113883.6.1'>**

- 3125 The <code> element indicates that this is a "clinically estimated" estimated delivery date (for example, this code is used to represent the field on the last line of the EDD section of the ACOG form). This code SHALL be the LOINC code 11778-8. It is good style to include the displayName and codeSystemName to help debugging.

**6.3.4.27.5 <value xsi:type='TS' value=' '>**

The value of the EDD SHALL be represented as a point in time.

- 3130 **6.3.4.27.6 <author typeCode='AUT'><assignedAuthor><id root=' ' extension=' /></assignedAuthor></author>**

3135 There may be multiple clinicians following the patient and authoring the overall document, however the EDD observation has an individual author. For CDA based content, this author SHALL be listed in the CDA header and referenced from the entry by including the id element of the assignedAuthor. For HL7 Version 3 Messages based content, the author SHALL be included in full through this element.

**6.3.4.27.7 <author typeCode='AUT'><time value=' '/></author>**

The author.time is used to record the time that the author recorded the observation. It SHALL be included.

- 3140 **6.3.4.27.8 <entryRelationship typeCode='SPRT'>**

3145 The <entryRelationship> element binds the clinicians estimated EDD to supporting observations by different methods. Supporting observations SHOULD be included. If included, the typeCode SHALL be 'SPRT'. For HL7 Version 3 Messages based content, the element name is <sourceOf> rather than <entryRelationship>, however the semantics, typeCode, and nested elements remain unchanged.

**6.3.4.27.9 <observation>**

3150           **<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13'>**  
               :  
               **</observation> [1st nesting]**

3155           Observations that support the clinical observation SHALL be included if known. These observations are the supporting calculated dates from various methods such as ultrasound dates or dates calculated from LMP (i.e., the left column of fields on the ACOG form). Supporting observations SHALL also conform to the simple observation template. Supporting observations MAY include a different effectiveTime, author, or informant. Supporting observations SHALL NOT include repeatNumber, interpretationCode, methodCode, or targetSiteCode. (Method is implied by the LOINC code). The templateId SHALL be valued as  
 3160           ‘1.3.6.1.4.1.19376.1.5.3.1.4.13’

**6.3.4.27.10 <code code=' ' codeSystem='2.16.840.1.113883.6.1'> [1st nesting]**

Supporting observations SHALL include one of following LOINC values to indicate the method used to calculate the EDD.

<b>Code</b>	<b>Description</b>
11779-6	Delivery date Estimated from last menstrual period
(xx-EDD-by-PE)	DELIVERY DATE-TMSTP-PT-^PATIENT-QN-ESTIMATED FROM CLINICIANS PHYSICAL EXAM
11781-2	Delivery date composite estimate
57063-0	Delivery date Estimated from quickening date
57064-8	Delivery date Estimated from date fundal height reaches umb

**3165 6.3.4.27.11 <entryRelationship typeCode='DRIV'>**

Observations of supporting EDD should provide observations from which they were derived such as the patient’s last menses, or gestational age value at a point in time.

For HL7 Version 3 Messages based content, the element name is <sourceOf> rather than <entryRelationship>, however the semantics, typeCode, and nested elements remain unchanged.

3170           **6.3.4.27.12 <observation>**  
               **<templateId root=' '/>**  
               :  
               **</observation> [2st nesting]**

3175           Observations that support the calculation of supporting observation SHALL be included if known. These observations are the supporting dates or ages from various methods such as ultrasound gestational age or the date of last Menses (for example, the right column of fields on the ACOG form). Supporting observations SHALL also conform to the simple observation template. Supporting observations MAY include a different effectiveTime, author, or informant.

3180 Supporting observations SHALL NOT include repeatNumber, interpretationCode, methodCode, or targetSiteCode. (Method is implied by the LOINC code)

#### **6.3.4.27.13 <code code=' ' codeSystem='2.16.840.1.113883.6.1'> [2nd nesting]**

This code is used to represent either the relevant date, or the gestational age observation from which the EDD is derived. The following table lists the relevant LOINC codes for methods used. For observations that record the gestational age the value is recorded as a physical quantity (PQ) with the units of weeks and the activity time should be recorded to indicate the date at which the gestational age was observed. For observations that simply record a date (e.g., LMP) the observation value is recorded as a point in time (TS).

Code	Description	Type
8655-2	DATE LAST MENSTRUAL PERIOD-TMSTP-PT-^PATIENT-QN-REPORTED	TS
11884-4	GESTATIONAL AGE-TIME-PT-^FETUS-QN-ESTIMATED FROM CLINICIANS PHYSICAL EXAM M	PQ
11888-5	Gestational age composite estimate	PQ
57065-5	Quickening date	TS
57066-3	Date fundal height reaches umbilicus	TS

3190 **6.3.4.27.14 <repeatNumber value=' '/> <interpretationCode code=' ' codeSystem=' '/> <targetSiteCode code=' ' codeSystem=' '/>**

The <repeatNumber> <interpretationCode>, and <targetSiteCode> elements should not be present in an EDD observation.

3195 *Add Section 6.3.4.28*

#### **6.3.4.28 Antepartum Visit Summary Battery 1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.2**

This entry describes a single row in the Antepartum Visit Summary Flowsheet. The single observation date and provider is applied to all other observations.

### 6.3.4.28.1 Specification

```

3200 <entry>
 <organizer classCode='BATTERY' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.2' />
 <id root=' ' extension=' ' />
 <code code='57061-4' displayName='Antepartum flowsheet panel'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <statusCode code='completed' />
 <author>
 <time value=' ' />
 <assignedAuthor>
 <id root=' ' extension=' ' />
 </assignedAuthor>
 </author>
 <component>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
 </observation>
 </component>
 <component>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
 </observation>
 </component>
 :
 </organizer>
 </entry>

```

### 6.3.4.28.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.2' />

3230 The <templateId> element specifies that this organizer entry conforms to the APS profile Antepartum Visit Summary Flowsheet battery. The root attribute SHALL contain the value "1.3.6.1.4.1.19376.1.5.3.1.1.11.2.3.2"

### 6.3.4.28.3 <organizer classCode='BATTERY' moodCode='EVN' >

3235 Each row in the visit Summary flowsheet of the Antepartum Summary SHALL be represented by an organizer with the classCode of 'BATTERY' and the moodCode of 'EVN'

### 6.3.4.28.4 <id root=' ' extension=' ' />

Each battery SHALL have a globally unique identifier.

### 6.3.4.28.5 <code code='(57061-4)' codeSystem='2.16.840.1.113883.6.1' />

3240 The <code> element specifies the LOINC code that represents the content of the battery. The codeSystem attribute SHALL contain the value '2.16.840.1.113883.6.1'. The code attribute SHALL contain the value='(57061-4)'. It is good practice to include displayName and codeSystemName for clarity and debugging. The corresponding values are 'ACOG VISIT SUMMARY BATTERY--PT--' and 'LOINC' respectively.

**6.3.4.28.6 <author><time><assignedAuthor><id></assignedAuthor></author>**

- 3245 The <author> relation element points at the author that records the visit battery. This assignedAuthor may be different than the author of the document. The time element is used to record when the assigned author recorded the battery.

**6.3.4.28.7 <statusCode code='completed' />**

The status code for all batteries SHALL be 'completed'

- 3250 **6.3.4.28.8 <component>**

The battery is made of several component simple observations. The following table lists the allowable LOINC codes, displayNames, and observation types, and unit of measures for these observations.

<b>LOINC Code</b>	<b>displayName</b>	<b>xsi:type</b>	<b>units</b>	<b>value set</b>
11884-4	Gestational age Clinical.estimate	PQ	week	
57067-1 or 11727-5 (by US)	Fetal Body weight Estimated by palpation or Fetal weight estimated by US	PQ	g, kg, lb_av, or oz_av	
11881-0	Uterus Fundal height Tape measure	PQ	cm	
11876-0 (by PE) or 11877-8 (by US)	Fetal presentation by palpitation or Fetal presentation US	CD		SNOMED CT Vertex (70028003) Breech (6096002) Transverse (73161006) Oblique (63750008) Compound (124736009) Brow (8014007) Face (21882006)
11948-7 or 57068-9	Fetal Heart rate US or Fetal Heart rate Auscultation	PQ	/min	
57088-7	Fetal Movement - Reported	CO		SNOMED CT fetal movement activity (finding) CID 364755008 baby kicks a lot (finding) CID 276368003 baby not moving (finding) CID 276370007 reduced fetal movement (finding) CID 276369006 fetal movements present (finding) CID 289431008 fetal movements felt (finding)

<b>LOINC Code</b>	<b>displayName</b>	<b>xsi:type</b>	<b>units</b>	<b>value set</b>
				CID 268470003 fetal movements seen (finding) CID 169731002
57069-7	Preterm labor symptoms	BL		
11709-7 or 11785-3	DILATION-LEN-PT-CERVICAL CANAL.external os -QN-PALPATION or DILATION-LEN-PT-CERVICAL CANAL.external os-QN-US	PQ	cm	
11867-9	Effacement Cervix by palpitation	PQ	percent	
11961-0	Cervix [Length] US	PQ	cm	
8480-6	Systolic blood pressure	PQ	mmHg	
8462-4	Diastolic blood pressure	PQ	mmHg	
3141-9	Body weight Measured	PQ	g, kg, lb_av, or oz_av	
1753-3	Albumin [Presence] in Urine	CO		SNOMED CT Negative (finding) CID 167273002 Trace (finding) CID 167274008 1+ (finding) CID 167275009 2+ (finding) CID 167276005 3+ (finding) CID 167277001 4+ (finding) CID 167278006
2349-9 or 25428-4(test strip)	Glucose [Presence] in Urine or Glucose [Presence] in Urine by Test strip	CO		SNOMED CT Negative (finding) CID 167261002 Trace (finding) CID 167262009 1+ (finding) CID 167264005 2+ (finding) CID 167265006 3+ (finding) CID 167266007 4+ (finding) CID 167267003
44966-0	Edema	CO		SNOMED CT Trace 4496-0 1+ pitting edema 420829009 2+ pitting edema 421605005 3+ pitting edema 421346005 4+ pitting edema 421129002
38208-5	Pain severity - Reported	CO		0 (no pain) :

<b>LOINC Code</b>	<b>displayName</b>	<b>xsi:type</b>	<b>units</b>	<b>value set</b>
				10 (worst possible pain) Note: This observation should correspond to the functional status <a href="#">pain score observation</a>
57070-5	Date next clinic visit	PQ	day,week,mo	
48767-8	Annotation comment	ED		

3255

Add Section 6.3.4.29

**6.3.4.29 Advance Directive Observation 1.3.6.1.4.1.19376.1.5.3.1.4.13.7**

An advance directive observation is a simple observation that uses a specific vocabulary, and inherits constraints from CCD.

3260

**6.3.4.29.1 Standards****CCD** [ASTM/HL7 Continuity of Care Document](#)**6.3.4.29.2 Specification**

3265

```
<observation typeCode='OBS' moodCode='EVN'>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
<templateId root='2.16.840.1.113883.10.20.1.17' />
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13.7' />
<id root=' ' extension=' '/>
<code code=' ' codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT' />
<text><reference value='#xxx' /></text>
<statusCode code='completed' />
<effectiveTime value=' ' />
<value xsi:type='BL' value='true|false' />
<reference typeCode='REFR'>
<templateId root='2.16.840.1.113883.10.20.1.36' />
<externalDocument classCode='DOC' moodCode='EVN'>
<id root=' ' extension=' '/>
<text><reference value=' ' /></text>
</externalDocument>
</reference>
</observation>
```

3270

3275

3280

An advanced directive <observation> shall be represented as shown above. They shall not contain any <repeatNumber>, <interpretationCode>, <methodCode> or <targetSiteCode> elements.

**6.3.4.29.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />  
<templateId root='2.16.840.1.113883.10.20.1.17' />  
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13.7' />**

3285

The <templateId> elements shown above shall be present, and indicated that this is an Advance Directive entry.

**6.3.4.29.4 <code code=' ' codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT'>**

- 3290 The <code> element records the type of advance directive. It should use one of the following SNOMED codes in the table below.

Code	Description	Data Type
304251008	Resuscitation	BL
52765003	Intubation	
225204009	IV Fluid and Support	
89666000	CPR	
281789004	Antibiotics	
78823007	Life Support	
61420007	Tube Feedings	
116859006	Transfusion of blood product	
71388002	Other Directive	
		<value> not permitted

**6.3.4.29.5 <value xsi:type='BL' value='true|false'>**

- 3295 The advance directive observation may include a <value> element using the Boolean (xsi:type='BL') data type to indicate simply whether the procedure described is permitted. Absence of the <value> element indicates that an advance directive of the specified type has been recorded, and must be examined to determine what type of treatment should be performed. The value element is not permitted when the <code> element describes an Other directive.

- 3300 **6.3.4.29.6 <reference typeCode='REFR'>**  
**<templateId root='2.16.840.1.113883.10.20.1.36'>**  
**<externalDocument classCode='DOC' moodCode='EVN'>**  
**<id root=' ' extension=' '/>**  
**<text><reference value=' '/></text>**

- 3305 The advanced directive observation may contain a single reference to an external document. That reference shall be recorded as shown above. The <id> element shall contain the appropriate root and extension attributes to identify the document. The <text> element may be present to provide a URL link to the document in the value attribute of the <reference> element. If the <reference> element is present, the Advance Directive in the narrative shall contain a <linkHTML> element to the same URL found in the value attribute.
- 3310

3315 *Add Section 6.3.4.30***6.3.4.30 Blood Type Observation 1.3.6.1.4.1.19376.1.5.3.1.4.13.6**

The blood type observation is a Simple Observation of the patient's blood type. It conforms to the CCD Result observation template.

**6.3.4.30.1 Standards**

CCD [ASTM/HL7 Continuity of Care Document](#)

3320 **6.3.4.30.2 Specification**

```

<observation typeCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13.6' />
 <templateId root='2.16.840.1.113883.10.20.1.31' />
 <id root=' ' extension=' '/>
 <code code='882-1' displayName='ABO+RH GROUP'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <text><reference value='#xxx' /></text>
 <statusCode code='completed' />
 <effectiveTime value=' ' />
 <repeatNumber value=' ' />
 <value xsi:type='CE' code=' ' displayName=' ' codeSystem=' ' codeSystemName=' ' />
 <interpretationCode code=' ' codeSystem=' ' codeSystemName=' ' />
 <methodCode code=' ' codeSystem=' ' codeSystemName=' ' />
 <targetSiteCode code=' ' codeSystem=' ' codeSystemName=' ' />
<observation>
```

**6.3.4.30.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />**

3340           **<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13.6' />**  
                 **<templateId root='2.16.840.1.113883.10.20.1.31' />**

These <templateId> elements identify this as a blood type observation. They shall be present in the <observation> element as shown above.

3345           **6.3.4.30.4 <code code='882-1' displayName='ABO+RH GROUP'**  
                  **codeSystem='2.16.840.1.113883.6.1'**  
                  **codeSystemName='LOINC' />**

The <code> element shall be present to represent this as a finding of the patient's composite blood type. It shall use the code and codeSystem attributes shown above.

3350           **6.3.4.30.5 <repeatNumber value=' ' />**

The <repeatNumber> element should not be present in a blood type observation.

**6.3.4.30.6 <value xsi:type='CE' code=' ' displayName=' ' codeSystem=' ' codeSystemName=' '/>**  
3355

The <value> element shall be present and shall use the CE data type. The code attribute should be valued using a vocabulary that supports encoding of blood types. The table below shows some coding systems that may be used to encode blood type.

Coding System	OID
ISBT 128	2.16.840.1.113883.6.18
SNOMED CT	2.16.840.1.113883.6.96

3360

**.3.4.30.7 <interpretationCode code=' ' codeSystem=' ' codeSystemName=' '/>  
~~<methodCode code=' ' codeSystem=' ' codeSystemName=' '/>~~  
~~<targetSiteCode code=' ' codeSystem=' ' codeSystemName=' '/>~~**

3365 The <interpretationCode>, <methodCode>, and <targetSiteCode> should not be present in a blood type observation.

<i>Add Section 6.3.4.31</i>
-----------------------------

**6.3.4.31 Encounters 1.3.6.1.4.1.19376.1.5.3.1.4.14**

3370 An Encounter is an interaction between a patient and care provider(s) for the purpose of providing healthcare-related service(s). Healthcare services include health assessment. Examples: outpatient visit to multiple departments, home health support (including physical therapy), inpatient hospital stay, emergency room visit, field visit (e.g., traffic accident), office visit, occupational therapy, or telephone call.

3375 **6.3.4.31.1 Standards**

CCD [ASTM/HL7 Continuity of Care Document](#)

**6.3.4.31.2 Specification**

```

3380 <encounter classCode='ENC' moodCode='PRMS|ARQ|EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.14' />
 <templateId root='2.16.840.1.113883.10.20.1.21' />
 <templateId root='2.16.840.1.113883.10.20.1.25' />
 <id root='' extension='' />
 <code code='' codeSystem='2.16.840.1.113883.5.4' codeSystemName='ActEncounterCode' />
 <text><reference value="#xxx"/></text>
 <effectiveTime>
 <low value='' />
 <high value='' />
 </effectiveTime>
 <priorityCode code='' />
 <performer typeCode='PRF'>
 <time><low value='' /><high value='' /></time>
 <assignedEntity>...</assignedEntity>
 </performer>
 <author />
 <informant />
 <participant typeCode='LOC'>
 <participantRole classCode='SDLOC'>
 <id />
 <code />
 <addr>...</addr>
 <telecom value='' use='' />
 <playingEntity classCode='PLC' determinerCode='INST'>
 <name></name>
 </playingEntity>
 </participantRole>
 </participant>
</encounter>

```

#### **6.3.4.31.2.1 <encounter classCode='ENC' moodCode='APT|ARQ|EVN'>**

- 3410 This element is an encounter. The classCode shall be 'ENC'. The moodCode may be PRMS to indicate a scheduled appointment, ARQ to describe a request for an appointment that has been made but not yet scheduled by a provider, or EVN, to describe an encounter that has already occurred.

#### **6.3.4.31.2.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.14'>**

- 3415 The templateId indicates that this <encounter> entry conforms to the constraints of this content module. NOTE: When the encounter is in event mood (moodCode='EVN'), this entry conforms to the CCD template 2.16.840.1.113883.10.20.1.21, and when in other moods, this entry conforms to the CCD template 2.16.840.1.113883.10.20.1.25.

#### **6.3.4.31.2.3 <id root=" extension="/">**

- 3420 This required element shall contain an identifier for the encounter. More than one encounter identifier may be present.

#### **6.3.4.31.2.4 <code code="" codeSystem='2.16.840.1.113883.5.4' codeSystemName='ActEncounterCode' />**

- 3425 This required element should contain a code from the HL7 ActEncounterCode vocabulary describing the type of encounter (e.g., inpatient, ambulatory, emergency, et cetera). Developers should take care to check that rational combinations of encounter.code and encounter.moodCode are used , but this Technical Framework does not restrict any combination.

**6.3.4.31.2.5 <text><reference value="#xxx'/'></text>**

The <text> element shall contain a reference to the narrative text describing the encounter.

**3430 6.3.4.31.2.6 <effectiveTime><low value="/"><high value="/"></effectiveTime>**

This element records the time over which the encounter occurred (in EVN mood), or the desired time of the encounter in ARQ or APT mood. In EVN or APT mood, the effectiveTime element should be present. In ARQ mood, the effectiveTime element may be present, and if not, the priorityCode may be present to indicate that a callback is required to schedule the appointment.

**3435 6.3.4.31.2.7 <priorityCode code='CS'>**

This element may be present in ARQ mood to indicate a callback is requested to schedule the appointment.

**6.3.4.31.2.8 <performer>**

For encounters in EVN mood, at least one performer should be present that identifies the provider of the service given during the encounter. More than one performer may be present. The <time> element should be used to indicate the duration of the participation of the performer when it is substantially different from that of the effectiveTime of the encounter. In ARQ mood, the performer may be present to indicate a preference for a specific provider. In APT mood, the performer may be present to indicate which provider is scheduled to perform the service.

**3445 6.3.4.31.2.9 <participant typeCode='LOC'>  
<participantRole classCode='SDLOC'>**

A <participant> element with typeCode='LOC' may be present to provide information about the location where the encounter is to be or was performed. This element shall have a <participantRole> element with classCode='SDLOC' that describes the service delivery location.

**3450 6.3.4.31.2.10 <id>**

The <id> element may be present to identify the service delivery location.

**6.3.4.31.2.11 <code/>**

The <code> element may be present to classify the service delivery location.

**6.3.4.31.2.12 <addr>...</addr>**

3455 The <addr> element should be present, and gives the address of the location.

**6.3.4.31.2.13 <telecom value=" use="/">**

The <telecom> element should be present, and gives the telephone number of the location.

**6.3.4.31.2.14 <playingEntity classCode='PLC'>**

3460      <name>...</name>  
           </playingEntity>

The <playingEntity> shall be present, and gives the name of the location in the required <name> element.

*Add Section 6.3.4.32*

3465      **6.3.4.32 Update Entry 1.3.6.1.4.1.19376.1.5.3.1.4.16**

The update entry shall contain references to the entries or sections which are being replaced or updated. This reference shall not be present when the update entry is adding a new entries or sections.

3470      Entries and sections can be added, updated, or removed from a PHR. An update entry indicates the entry in the original PHR Extract that should be replaced or updated with new information contained within the entry. Only one organizer of this type is allowed in a section, and if present, it must be the first entry in the section.

**6.3.4.32.1 Specification**

3475      <entry>  
           <organizer classCode='BATTERY' moodCode='EVN'>  
             <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.16' />  
             <reference typeCode='RPLC'>  
               <externalAct classCode='ACT' moodCode='EVN'>  
                 <id root='' extension='' />  
               </externalAct>  
             </reference>  
           </organizer>  
</entry>

**6.3.4.32.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.16'>**

This templateId indicates that the organizer is used to update a PHR Extract.

3490      **6.3.4.32.3 <reference typeCode='RPLC'>**

A reference element shall be present with typeCode RPLC. The reference element lists the acts that are affected by the update. It indicates that any referenced act is being replaced with new information. This element must be present, and may be repeated to replace more than one act at a time.

3495

#### **6.3.4.32.4 <externalAct classCode='ACT' moodCode='EVN'>**

This element must appear as shown above. It indicates that the reference is to an external act (a section or entry contained in the parent document).

3500

#### **6.3.4.32.5 <id root=' ' extension=' '/>**

This element identifies the information being replaced or updated. The identifier is of the entry or section being replaced. If the identifier is to a section being replaced, only one reference element is permitted.

3505

*Add Section 6.3.4.33*

#### **6.3.4.33 Procedure Entry 1.3.6.1.4.1.19376.1.5.3.1.4.19**

The procedure entry is used to record procedures that have occurred, or which are planned for in the future.

3510

#### **6.3.4.33.1 Standards**

CCD [ASTM/HL7 Continuity of Care Document](#)

### 6.3.4.33.2 Specification

```

<procedure classCode='PROC' moodCode='EVN|INT'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.19' />
 <templateId root='2.16.840.1.113883.10.20.1.29' /><!-- see text of section 0 -->
 <templateId root='2.16.840.1.113883.10.20.1.25' /><!-- see text of section 0 -->
 <id root='' extension=''/>
 <code code='' codeSystem='2.16.840.1.113883.5.4' codeSystemName='ActCode' />
 <text><reference value="#xxx"/></text>
 <statusCode code='completed|active|aborted|cancelled' />
 <effectiveTime>
 <low value=''/>
 <high value=''/>
 </effectiveTime>
 <priorityCode code=''/>
 <approachSiteCode code='' displayName='' codeSystem='' codeSystemName='' />
 <targetSiteCode code='' displayName='' codeSystem='' codeSystemName='' />
 <author />
 <informant />
 <entryRelationship typeCode='COMP' inversionInd='true'>
 <act classCode='ACT' moodCode=''>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.4.1' />
 <id root='' extension=''/>
 </act>
 </entryRelationship>
 <entryRelationship typeCode='RSON'>
 <act classCode='ACT' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.4.1' />
 <id root='' extension=''/>
 </act>
 </entryRelationship>
</procedure>

```

#### 6.3.4.33.2.1 <procedure classCode='PROC' moodCode='EVN|INT'>

This element is a procedure. The classCode shall be 'PROC'. The moodCode may be INT to indicate a planned procedure or EVN, to describe a procedure that has already occurred.

3545 **6.3.4.33.2.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.19' />**

The templateId indicates that this <procedure> entry conforms to the constraints of this content module. NOTE: When the procedure is in event mood (moodCode='EVN'), this entry conforms to the CCD template 2.16.840.1.113883.10.20.1.29, and when in intent mood, this entry conforms to the CCD template 2.16.840.1.113883.10.20.1.25.

3550 **6.3.4.33.2.3 <id root=" extension="/">**

This required element shall contain an identifier for the procedure. More than one procedure identifier may be present.

#### 6.3.4.33.2.4 <code code="" displayName="" codeSystem="" codeSystemName="" />

This element shall be present, and should contain a code describing the type of procedure.

3555 **6.3.4.33.2.5 <text><reference value="#xxx"/></text>**

The <text> element shall contain a reference to the narrative text describing the procedure.

**6.3.4.33.2.6 <statusCode code='completed|active|aborted|cancelled'>**

3560 The <statusCode> element shall be present when used to describe a procedure event. It shall have the value 'completed' for procedures that have been completed, and 'active' for procedures that are still in progress. Procedures that were stopped prior to completion shall use the value 'aborted', and procedures that were cancelled before being started shall use the value 'cancelled'.

**6.3.4.33.2.7 <effectiveTime><low value="/"><high value="/"></effectiveTime>**

This element should be present, and records the time at which the procedure occurred (in EVN mood), or the desired time of the procedure in INT mood.

3565 **6.3.4.33.2.8 <priorityCode code="/">**

This element shall be present in INT mood when effectiveTime is not provided, it may be present in other moods. It indicates the priority of the procedure.

**6.3.4.33.2.9 <approachSiteCode code=" displayName=" codeSystem=" codeSystemName="/">**

3570 This element may be present to indicate the procedure approach.

**6.3.4.33.2.10 <targetSiteCode code=" displayName=" codeSystem=" codeSystemName="/">**

This element may be present to indicate the target site of the procedure.

**6.3.4.33.2.11 <entryRelationship typeCode='COMP' inversionInd='true'>**

3575 This element may be present to point the encounter in which the procedure was performed, and shall contain an internal reference to the encounter. See PCC TF-2: 6.3.4.10 Internal References for more details.

**6.3.4.33.2.12 <entryRelationship typeCode='RSON'>**

3580 A <procedure> act may indicate one or more reasons for the procedure. These reasons identify the concern that was the reason for the procedure via an Internal Reference (see PCC TF-2: 6.3.4.10 Internal References) to the concern. The extension and root of each observation present must match the identifier of a concern entry contained elsewhere within the CDA document.

Add Section 6.3.4.34

3585 **6.3.4.34 Transport 1.3.6.1.4.1.19376.1.5.3.1.1.10.4.1**

A transport entry indicates the intended or actual mode of transport and time of departure and/or arrival of the patient.

### 6.3.4.34.1 Specification

```

3590 <entry>
 <act classCode='ACT' moodCode='INT|EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.10.4.1' />
 <id root='' extension=''/>
 <code code='' displayName=''
 codeSystem='2.16.840.1.113883.6.102.4.2'
 codeSystemName='DEEDS4.02'>
 <originalText><reference value='#(ID of text coded)' /></originalText>
 </code>
 <text><reference value='#text' /></text>
 <high value=/>
 </effectiveTime>
 </act>
 </entry>

```

#### 6.3.4.34.1.1 <act classCode='ACT' moodCode='INT|EVN'>

This element indicates that the entry is an act (of transporting the patient, as indicated by the code below). This entry records the mode, and intended or actual ending time of transportation. In intent mood (moodCode='INT') this is how the estimated time of departure or arrival is indicated. In event mood (moodCode='EVN') this is how the actual departure or arrival of the patient is recorded.

#### 6.3.4.34.1.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.10.4.1'>

The <templateId> element identifies this <act> as about the transportation of the patient. The templateId must have root='1.3.6.1.4.1.19376.1.5.3.1.1.10.4.1'.

#### 6.3.4.34.1.3 <id root=" extension=""/>

The entry must have an identifier.

#### 6.3.4.34.1.4 <code code="" displayName=""

3615 **codeSystem='2.16.840.1.113883.6.102.4.2' codeSystemName='DEEDS4.02'>**

The code describes the intended mode of transport. For transport between facilities, IHE recommends the use of a code system based on the DEEDS Mode of Transportation data element value set. However, the vocabulary used within an affinity domain should be determined by a policy agreement within the domain.

3620 **6.3.4.34.1.5 <originalText><reference value='#xxx' /><originalText>**

This is a reference to the narrative text within the section that describes the mode of transportation.

#### 6.3.4.34.1.6 <text><reference value='#text' /></text>

This is a reference to the narrative text corresponding to the transport act.

3625 **6.3.4.34.1.7 <effectiveTime>**

The effectiveTime element shall be sent. It records the interval of time over which transport occurs.

**6.3.4.34.1.8 <low value="/">**

This element records the time of departure. This element shall be sent using the TS data type, as shown above.

**6.3.4.34.1.9 <high value="/">**

This element records the time of arrival. If unknown, it must be recorded using a flavor of null. This element shall be sent using the TS data type as shown above.

3635 *Add Section 6.3.4.35*

**6.3.4.35 Encounter Disposition 1.3.6.1.4.1.19376.1.5.3.1.1.10.4.2**

This element records the intended or actual disposition for the patient (e.g., admit, discharge home after treatment, et cetera).

### 6.3.4.35.1 Specification

```

3640 <act classCode='ACT' moodCode='INT|EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.10.4.2' />
 <id root='' extension=''/>
 <code code='' displayName='' codeSystem='' codeSystemName='' />
 <text><reference value='#xxx' /></text>
 <statusCode code='normal|completed' />
 <effectiveTime value=''/>
 <performer typeCode='PRF'>
 <assignedEntity>
 <id root='' extension=''/>
 <addr></addr>
 <telecom value='' use=''/>
 <assignedPerson>
 <name></name>
 </assignedPerson>
 </assignedEntity>
 </performer>
 <participant typeCode='RCV'>
 <time value=''/>
 <participantRole classCode='ROL'>
 <id root='' extension=''/>
 <addr></addr>
 <telecom value='' use=''/>
 <playingEntity>
 <name></name>
 </playingEntity>
 </participantRole>
 </participant>
 <entryRelationship typeCode='COMP'>
 <act classCode='ACT'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.10.4.1' />
 :
 </act>
 </entryRelationship>
 </act>

```

#### 6.3.4.35.1.1 <act classCode='ACT' moodCode='INT|EVN'>

The disposition is recorded in an act element, to describe the disposition action taken during the encounter<sup>1</sup>. In intent mood (moodCode='INT'), this records the expected disposition of the patient. In event mood (moodCode='EVN'), this records the actual disposition.

<sup>1</sup> The HL7 RIM allows this portion of the encounter to be recorded in the dischargeDispositionCode RIM Attribute of the Encounter class, but the Encounter class is constrained within CDA. To record the disposition act therefore requires the use of the Act class.

#### 6.3.4.35.1.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.10.4.2' />

3680 The templateId indicates that this <encounter> entry conforms to the constraints of this content module.

#### 6.3.4.35.1.3 <id root=" extension="/">

This required element shall contain an identifier.

**6.3.4.35.1.4 <code code="" displayName="" codeSystem="" codeSystemName="" />**

- 3685 This required element indicates the disposition of the patient. The code shall come from a coding system that is able to record common patient dispositions (e.g., Discharged, Transferred, Admitted). The "Administrative Procedure" concept (14734007) of SNOMED CT contains several code values that cover a wide variety of dispositions routinely recorded. Other vocabularies that are commonly in use to describe discharge disposition codes are DEEDS (See section 8.02), and in the US, the Uniform National Billing Code.
- 3690

**6.3.4.35.1.5 <text><reference value="#xxx"/></text>**

- The <text> element shall contain a reference to the narrative text describing the disposition of the patient. <statusCode code='normal|completed' /> When the disposition act has occurred (moodCode='EVN'), the statusCode element shall be present, and shall contain the value 'completed'. When the disposition act is intended (moodCode='EVN') the statusCode element shall contain the value 'normal'.
- 3695

**6.3.4.35.1.6 <effectiveTime><low value="/" /><high value="/" /><effectiveTime/>**

- 3700 When the disposition has occurred, this element shall be sent, and indicates the effective time for the disposition process. This element may be sent to record when the disposition act is intended to occur. The <low> element records the time at which the disposition process was started. The <high> value records the time at which the disposition process was completed.

**6.3.4.35.1.7 <performer typeCode='PRF' />**

- 3705 The <performer> element provides information about the person that performs the discharge, admission or transfer of the patient. When the disposition is in intent mood, this element describes any expectations with respect to the performer, and is optional. When the disposition is in event mood, this element is required.

**6.3.4.35.1.8 <assignedEntity>**

The <assignedEntity> element identifies the performer of the disposition.

**6.3.4.35.1.9 <id root="" extension="/" />**

- 3710 The <id> element shall be sent when the disposition has occurred, and identifies the performer of the act.

**6.3.4.35.1.10 <addr></addr>**

The <addr> element may be sent to provide a contact postal address for the performer of the disposition.

- 3715 **6.3.4.35.1.11 <telecom value="" use="/" />**

The <telecom> element may be sent to provide a contact postal address for the performer of the disposition.

**6.3.4.35.1.12 <assignedPerson><name/></assignedPerson>**

3720 The <assignedPerson> element shall be sent to identify the person who performed the disposition of the patient.

**6.3.4.35.1.13 <participant typeCode='RCV'>**

```
<time value="/" />
<participantRole classCode='ROL'>
<id root=" extension="/" />
<addr></addr>
<telecom value=" use="/" />
<playingEntity><name/></playingEntity>
```

3725

This element identifies the person or organization that is receiving the patient. =====

3730

<entryRelationship typeCode='COMP'>

<act classCode='ACT'>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.10.4.1' /> If the disposition of the patient requires transport to another location, this information shall be recorded in a subordinate act that conforms to the Transport template described above.

3735

*Add Section 6.3.4.36*

**6.3.4.36 Coverage Entry 1.3.6.1.4.1.19376.1.5.3.1.4.17**

Payers shall be recorded as described in CCD: 3.1.2.1.1.

**6.3.4.36.1 Standards**

CCD [ASTM/HL7 Continuity of Care Document](#)

**6.3.4.36.2 Specification**

3740

**Coverage Entry Example**

3745

```
<act classCode='ACT' moodCode='DEF'>
<templateId root='2.16.840.1.113883.10.20.1.20' />
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.17' />
<id root='' extension=''/>
<code code='48768-6' displayName='Payment Sources'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<statusCode code='completed' />
<entryRelationship typeCode='COMP'>
 <sequenceNumber value=''/>
 :
</entryRelationship>
</act>
```

3750

3755

#### **6.3.4.36.3 <act classCode='ACT' moodCode='DEF'>**

Coverage shall be recorded in an <act> that groups all patient coverage together, and defines (moodCode='DEF') the payers.

- 3760    **6.3.4.36.4 <templateId root='2.16.840.1.113883.10.20.1.20'>  
                <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.17'>**

The <act> conforms to CCD: 3.1.2.1.1 as well as this specification. This shall be reflected by including the <templateId> elements shown above.

- 3765    **6.3.4.36.5 <id root=' ' extension=' '/>**

The <id> element shall be present.

#### **6.3.4.36.6 <code code='48768-6' displayName='PAYMENT SOURCES' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'>**

- 3770    The <code> element shall be recorded exactly as shown above.

#### **6.3.4.36.7 <statusCode code='completed'>**

The <statusCode> element shall be present exactly as shown above.

- 3775    **6.3.4.36.8 <entryRelationship typeCode='COMP'>**

The coverage <act> shall have one or more <entryRelationship> elements. These elements define the coverage. The entry relationships must contain Payer Entries.

#### **6.3.4.36.9 <sequenceNumber value=' '/>**

- 3780    The <sequenceNumber> element may be present. If present, it shall contain a value attribute that indicates the priority of the payment source.

*Add Section 6.3.4.37*

#### **6.3.4.37 Payer Entry 1.3.6.1.4.1.19376.1.5.3.1.4.18**

- 3785    The payer entry allows information about the patient's sources of payment to be recorded.

### 6.3.4.37.1 Standards

CCD [ASTM/HL7 Continuity of Care Document](#)

### 6.3.4.37.2 Specification

#### Payer Entry Example

```

3790 <act classCode='ACT' moodCode='EVN'>
 <templateId root='2.16.840.1.113883.10.20.1.26' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.18' />
 <id root='' extension='' />
3795 <code code='' displayName='' codeSystem='' codeSystemName='' />
 <statusCode code='completed' />
 <performer typeCode='PRF'><!-- payer -->
 <assignedEntity classCode='ASSIGNED'>
 <id root='' extension='' />
 <code code='PAYOR|GUAR|PAT' displayName=''
 codeSystem='2.16.840.1.113883.5.110' codeSystemName='RoleClass' />
 <addr></addr>
 <telecom value='' use='' />
 <representedOrganization typeCode='ORG'>
 <name></name>
 </representedOrganization>
 </assignedEntity>
 </performer>
3800 <participant typeCode='COV'><!-- member -->
 <participantRole classCode='PAT'>
 <id root='' extension='' />
 <code code='' displayName=''
 codeSystem='2.16.840.1.113883.5.111' codeSystemName='RoleCode' />
 <addr></addr>
 <telecom value='' use='' />
 <playingEntity><name></name></playingEntity>
 </participantRole>
 </participant>
3805 <participant typeCode='HLD'><!-- subscriber -->
 <participantRole classCode='PAT'>
 <id root='' extension='' />
 <playingEntity><name></name></playingEntity>
 </participantRole>
 </participant>
3810 <entryRelationship typeCode='REFR'>
 <act classCode='ACT' moodCode='DEF'>
 <id root='' extension='' />
 <code code='' displayName='' codeSystem='' codeSystemName='' />
 <text><reference value='/'></text>
 </act>
 </entryRelationship>
3815 </act>
</act>
```

### 3835 6.3.4.37.3 <act classCode='ACT' moodCode='EVN'>

The policy entry <act> describes the policy or program that has agreed to pay (moodCode='EVN') for the patient's treatment.

**6.3.4.37.4 <templateId root='2.16.840.1.113883.10.20.1.26'>  
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.18'>**

3840 The <act> conforms to CCD: 3.1.2.1.2 and this guide. This shall be reflected by including the <templateId> elements shown above.

**6.3.4.37.5 <id root=' ' extension=' '/>**

3845 The <act> shall contain at least one <id> element that represents the policy or group number of the coverage. That identifier shall appear in the extension attribute.

**6.3.4.37.6 <code code=' ' displayName=' '  
 codeSystem=' ' codeSystemName=' '/>**

3850 The <code> element should be present, and represents the type of coverage provided by the payer. Potential vocabularies to use include:

**Table 6.3.4.37.2-1: Payer Type Vocabularies**

Vocabulary	Description	OID
HL7 ActCoverageType	The HL7 ActCoverageType vocabulary describes payers and programs. Note that HL7 does not have a specific code to identify an individual payer, e.g., in the role of a guarantor or patient.	2.16.840.1.113883.5.4
X12 Data Element 1336	The X12N 271 implementation guide includes various types of payers. This code set does include a code to identify individual payers.	2.16.840.1.113883.6.255.1336

**6.3.4.37.7 <statusCode code='completed'>**

3855 The <statusCode> element shall be present, and should be recorded exactly as shown above.

**6.3.4.37.8 <performer typeCode='PRF'>  
 <assignedEntity classCode='ASSIGNED'>**

The <performer> element shall be present to represent the payer of the coverage.

3860

**6.3.4.37.9 <id root=' ' extension=' '/>**

The identity of the performer should be recorded in the <id> element.

3865      **6.3.4.37.10 <code code='PAYOR|GUAR|PAT' displayName=' '  
              codeSystem='2.16.840.1.113883.5.110'  
              codeSystemName='RoleClass' />**

The <code> element describes the role of the payer. It shall contain one of the values listed in the table below.

**Payer Role Codes**

Coding System	OID
ISBT 128	2.16.840.1.113883.6.18
SNOMED CT	2.16.840.1.113883.6.96

3870      **6.3.4.37.11 <addr></addr>**

The <addr> element shall be used to record the address of the payer. This information will usually come from the back of an insurance card.

**6.3.4.37.12 <telecom value=' ' use=' '/>**

3875      The <telecom> element shall be used to record the phone number of the payer. This information will usually come from the back of an insurance card.

**6.3.4.37.13 <representedOrganization typeCode='ORG'>  
<name></name>**

3880      The name of the payer organization shall be provided in the <name> element contained within the <representedOrganization> element.

**6.3.4.37.14 <participant typeCode='COV'>  
<participantRole classCode='PAT'>**

3885      Information about the patient with respect to the policy or program shall be recorded in the <participantRole> element shown above. This element shall be present when the patient is a member of a policy or program.

**6.3.4.37.15 <id root=' ' extension=' '/>**

3890      The <id> element should contain the identifier of the patient with respect to the payer (the subscriber or member id).

**6.3.4.37.16 <code code= displayName= codeSystem='2.16.840.1.113883.5.111' codeSystemName='RoleCode'/>**

3895 The <code> element shall indicate the covered party's relationship to the subscriber, and should come from the HL7 CoverageRoleType value set.

**6.3.4.37.17 <addr></addr>**

3900 The <addr> element should be used to record the address of the patient as known to the payer when different from that recorded in the <patientRole> element.

**6.3.4.37.18 <telecom value=' ' use=' '/>**

The <telecom> element should be used to record the phone number of the patient when different from that recorded in the <patientRole> element.

3905

**6.3.4.37.19 <playingEntity><name></name></playingEntity>**

The <name> element should be used to record the member name when it is different from that recorded in the <patient> element.

3910

**6.3.4.37.20 <participant typeCode='HLD'>  
                  <participantRole classCode='IND'>**

Information about subscriber to the policy or program shall be recorded in the <participantRole> element shown above. This element shall be present when the subscriber is different from the patient.

3915

**6.3.4.37.21 <id root=' ' extension=' '/>**

The <id> element shall contain the identifier of the subscriber when the subscriber is not the patient.

3920

**6.3.4.37.22 <addr></addr>**

The <addr> element shall be used to record the address of the subscriber when the subscriber is not the patient.

**6.3.4.37.23 <telecom value=' ' use=' '/>**

3925 The <telecom> element shall be used to record the phone number of the subscriber when the subscriber is not the patient.

**6.3.4.37.24 <playingEntity><name></name></playingEntity>**

The name of the subscriber shall be recorded in the <name> element of the <playingEntity>.

3930

**6.3.4.37.25 <entryRelationship typeCode='REFR'>  
    <act classCode='ACT' moodCode='DEF'>**

The plan information may be provided in the elements described above.

3935

**6.3.4.37.26 <id root=' ' extension=' '/>**

The health plan identifier is recorded in the <id> element.

**6.3.4.37.27 <text><reference value=' '/></text>**

3940 This <reference> element shown above should be present and the value attribute should point to the name of the plan contained in the narrative of the document.

*Add Section 6.3.4.38*

**6.3.4.38 Pain Score Observation 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.1**

3945 The pain score observation is a [Simple Observation](#) that records the patient's assessment of their pain on a scale from 0 to 10.

**6.3.4.38.1 Parent Template**

The parent of this template is [Simple Observation](#).

### 6.3.4.38.2 Specification

```

3950 <observation typeCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 <id root=' ' extension=' '/>
 <code code='38208-5|38221-8|38214-3' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'>
 <translation code='406127006' displayName='Pain intensity'
 codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT' />
 </code>
 <text><reference value='#xxx' /></text>
 <statusCode code='completed' />
 <effectiveTime value=' ' />
 <repeatNumber value=' ' />
 <value xsi:type='CO|REAL' />
 <interpretationCode code= codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT' />
 <methodCode code=' ' codeSystem=' ' codeSystemName=' ' />
 <targetSiteCode code=' ' codeSystem=' ' codeSystemName=' ' />
</observation>

```

### 6.3.4.38.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13'>

The <templateId> identifies this as a Pain Score Observation, and shall be present as shown above.

3970    **6.3.4.38.4 <code code='38208 5' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC'>**  
**<translation code='406127006' displayName='Pain intensity'**  
**codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT' />**

3975    The <code> element indicates what kind of pain observation was made. It shall contain the code and codeSystem attribute values shown above. The <translation> element may be present, and provides a mapping to SNOMED CT of the observation. If present, is shall have the code and codeSystem attribute values shown above.

Code	Data Type	Description
38208-5	CO	A Pain Score made using the Numerical Rating Scale (NRS), where pain is assessed on a scale from 0 to 10. -->The code system to use for this observation<<--

3980    **6.3.4.38.5 <value xsi:type='CO' value=' ' />**

The <value> element records the assessed pain score. If using the NRS the pain is assessed using coded ordinal values that range from 0 to 10. The use of the coded ordinal type is required because while pain assessments are ordered values, and can be compared, the differences between two pain assessment values cannot be compared, and so these values are not really numbers.

**6.3.4.38.6 <interpretationCode>**

**code='301379001|40196000|76948002|67849003'**  
**codeSystem='2.16.840.1.113883.6.96'**  
**codeSystemName='SNOMED CT'>**

- 3990 The <interpretationCode> element should be present to provide an interpretation of the pain scale assessment using SNOMED CT. When the <interpretationCode> element is present, the <translation> element described above shall be present. These interpretations are provided to assist decision support systems that are making secondary use of the assessment information, and are not intended to replace the score values.

3995

Pain Score Range	Code	Description
0	301379001	No Present Pain
1-3	40196000	Mild Pain
4-6	50415004	Moderate Pain
7-9	76948002	Severe Pain
10	67849003	Excruciating Pain

**6.3.4.38.7 <methodCode code=' ' codeSystem=' ' codeSystemName=' '/>**

The <methodCode> should not be present in a Pain Score Observation, as the method is implied by the <code> element.

4000

**6.3.4.38.8 <targetSiteCode code=' ' codeSystem=' ' codeSystemName=' '/>**

The <targetSiteCode> element should be present, and shall indicate the location of the pain being assessed.

*Add Section 6.3.4.39*

4005

**6.3.4.39 Braden Score Observation 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.2***Add Section 6.3.4.40***6.3.4.40 Braden Score Component 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.3**

4010

*Add Section 6.3.4.41***6.3.4.41 Geriatric Depression Score Observation 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.4**

*Add Section 6.3.4.42*

#### **6.3.4.42 Geriatric Depression Score Component 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.5**

4015

*Add Section 6.3.4.43*

#### **6.3.4.43 Survey Panel 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.7**

A survey panel collects related survey observations.

##### **6.3.4.43.1 Parent Template**

4020 This template is compatible with the ASTM/HL7 Continuity of Care Document template:  
2.16.840.1.113883.10.20.1.32

##### **6.3.4.43.2 Uses**

See Templates using [Survey Panel](#)

##### **6.3.4.43.3 Specification**

4025

```
<organizer classCode='CLUSTER' moodCode='EVN'>
 <templateId root='2.16.840.1.113883.10.20.1.32' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.7' />
 <id root='' extension='' />
 <code code=' ' displayName=' '
 codeSystem=' ' codeSystemName=' ' />
 <statusCode code='completed' />
 <effectiveTime value=''/>
 <!-- one or more survey observations -->
 <component typeCode='COMP'>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.6' />
 :
 </observation>
 </component>
 </organizer>
```

###### **6.3.4.43.3.1 <organizer classCode='CLUSTER' moodCode='EVN'>**

The survey panel is a cluster of related survey observations.

###### **6.3.4.43.3.2 <templateId root='2.16.840.1.113883.10.20.1.32' />**

**<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.7' />**

4030

The survey panel shall have the <templateId> elements shown above to indicate that it inherits constraints from the ASTM/HL7 CCD Specification for results organizers, and the constraints of this specification.

**6.3.4.43.3.3 <id root=' ' extension=' '/>**

4035 The organizer shall have an <id> element.

**6.3.4.43.3.4 <code code=' ' displayName=' ' codeSystem=' ' codeSystemName=' '/>**

The <code> element shall be present, and identifies the survey panel.

4040 **6.3.4.43.3.5 <statusCode code='completed' />**

The observations have all been completed.

**6.3.4.43.3.6 <effectiveTime value=' '/>**

The effective time element shall be present to indicate when the survey panel was taken.

**6.3.4.43.3.7 <!-- one or more survey observations -->  
<component typeCode='COMP'>**

4045

The organizer shall have one or more <component> elements that are <observation> elements using the [Survey Observation](#) template.

4050

*Add Section 6.3.4.44*

**6.3.4.44 Survey Observation 1.3.6.1.4.1.19376.1.5.3.1.1.12.3.6**

4055

Survey observations are used to record responses to assessment instruments. They are simple observations conforming to the CCD Result template. The vocabulary and data type constraints on survey observations is specified elsewhere, either in the specializations of the survey observation template, or by the template that makes use of it.

**6.3.4.44.1 Parent Template**

4060

See [Templates using Survey Observation](#)

**6.3.4.44.3 Specification**

```
<observation classCode='OBS' moodCode='EVN'>
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
<templateId root='2.16.840.1.113883.10.20.1.31' />
<templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.6' />
<id root=' ' extension=' ' />
```

```
<code code=' ' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
<text><reference value="#xxx"/></text>
<statusCode code='completed' />
<effectiveTime value=' ' />
<repeatNumber value=' ' />
<value xsi:type='CO|CD|INT|PQ' />
<interpretationCode code=' ' codeSystem=' ' codeSystemName=' ' />
<methodCode code=' ' codeSystem=' ' codeSystemName=' ' />
<targetSiteCode code=' ' codeSystem=' ' codeSystemName=' ' />
</observation>
```

**6.3.4.44.3.1 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13'>  
    <templateId root='2.16.840.1.113883.10.20.1.31'>  
        <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.12.3.6'>**

- 4065 A survey observation shall have the <templateId> elements shown above to indicate that it inherits constraints from the ASTM/HL7 CCD Specification for results, and the constraints of this specification.

**6.3.4.44.3.2 <code code=' ' codeSystem='2.16.840.1.113883.6.1'  
    codeSystemName='LOINC' />**

- 4070 A survey observation entry shall contain a code identifying the observation made.

**6.3.4.44.3.3 <value xsi:type='CO|CD|INT|PQ' .../>**

The <value> element shall be present, and shall be of the appropriate data type specified for the observation.

**6.3.4.44.3.4 <interpretationCode code=' ' codeSystem=' ' codeSystemName=' ' />**

- 4075 An interpretation code may be present to provide an interpretation of the observation.

**6.3.4.44.3.5 <methodCode code=' ' codeSystem=' ' codeSystemName=' ' />  
    <targetSiteCode code=' ' codeSystem=' ' codeSystemName=' ' />**

The <methodCode> and <targetSiteCode> element shall not be present, as these are not relevant to survey responses.

- 4080

*Add Section 6.3.4.45*

**6.3.4.45 Acuity 1.3.6.1.4.1.19376.1.5.3.1.1.13.3.1**

An acuity entry indicates the triage acuity entry and the triage time of the patient.

### 6.3.4.45.1 Specification

```

4085 <entry>
4086 <!-- Acuity Event -->
4087 <observation classCode='OBS' moodCode='EVN'>
4088 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.1' />
4089 <id root='' extension=''/>
4090 <code code='' displayName=''
4091 <code code='273887006' displayName='Triage index'
4092 codeSystem='2.16.840.1.113883.6.96'
4093 codeSystemName='SNOMED CT' /> <!-- Triage index (assessment scale) FullySpecifiedName --
4094 ->
4095 <originalText><reference value='#{ID of text coded}' /></originalText>
4096 </code>
4097 <text><reference value='#text' /></text>
4098 <!-- effectiveTime
4099 <effectiveTime>
4100 <low value=''/> <!-- start of triage, may be sent -->
4101 <high value=''/><!-- end of triage should be sent -->
4102 </effectiveTime>
4103 </observation>
4104 </entry>

```

#### 4105 6.3.4.45.1.1 <observation classCode='OBS' moodCode='EVN'>

This element indicates that the entry is an observation regarding the event of triage assessment. This entry records the observation and the time of the observation.

#### 6.3.4.45.1.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.1' />

The <templateId> element identifies this <act> as about Acuity Assessment of the patient. The templateId must have root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.1'.

#### 6.3.4.45.1.3 <id root=" extension=" />

The entry must have an identifier.

#### 6.3.4.45.1.4 <code code=" displayName=" 4115 codeSystem='2.16.840.1.113883.6.96' codeSystemName='SNOMED CT'>

The code describes the triage acuity scale. IHE recommends the use the Emergency Severity Index (ESI). However, the vocabulary used within an affinity domain may be determined by a policy agreement within the domain.

#### 6.3.4.45.1.5 <originalText><reference value='#xxx' /><originalText>

This is a reference to the narrative text within the section that describes the acuity description.

#### 4120 6.3.4.45.1.6 <text><reference value='#text' /></text>

This is a reference to the narrative text corresponding to the Observation act.

#### 6.3.4.45.1.7 <effectiveTime>

The effectiveTime element shall be sent. It records the interval of time over which triage occurs. The use case for this information requires that the ending time of triage be recorded. However,

- 4125 the <low value=""> element may be sent by systems that capture the beginning and end of the triage process.

#### **6.3.4.45.1.8 <high value=""/>**

- This element records the time of completion of triage, and is required. If unknown, it must be recorded using a flavor of null. This element may be sent using the TS data type, as shown above. If there is uncertainty about the time of completion of triage, the sender may record the time using the IVL\_TS data type, as shown below.

4135

```
<high xsi:type='IVL_TS'>
 <low value=''/>
 <high value=''/>
</high>
```

*Add Section 6.3.4.46*

- 4140 **6.3.4.46 Intravenous Fluids 1.3.6.1.4.1.19376.1.5.3.1.1.13.3.2**

This content module describes the general structure for intravenous fluids. All intravenous fluid administration acts should be derived from this content module.

### 6.3.4.46.1 Specification

```

4145 <substanceAdministration classCode='SBADM' moodCode='INT|EVN'>
4146 <templateId root='2.16.840.1.113883.10.20.1.24' />
4147 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.7' />
4148 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.7.1' />
4149 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.2' />
4150 <id root='' extension=''/>
4151 <code code='' codeSystem='' displayName='' codeSystemName='' />
4152 <text><reference value="#med-1"/></text>
4153 <statusCode code='completed|active' />
4154 <effectiveTime xsi:type='IVL_TS'>
4155 <low value=''/>
4156 <high value=''/>
4157 </effectiveTime>
4158 <effectiveTime operator='A' xsi:type='TS|PIVL_TS|EIVL_TS|PIVL_PPD_TS|SXPR_TS'>
4159 :
4160 </effectiveTime>
4161 <routeCode code='' codeSystem='' displayName='' codeSystemName='' />
4162 <doseQuantity value='' unit='' />
4163 <approachSiteCode code='' codeSystem='' displayName='' codeSystemName='' />
4164 <rateQuantity value='' unit='' />
4165 <consumable>
4166 :
4167 </consumable>
4168 <!-- 0..* entries describing the components -->
4169 <entryRelationship typeCode='COMP' >
4170 <sequenceNumber value=''/>
4171 </entryRelationship>
4172 <!-- An optional entry relationship that indicates the reason for use -->
4173 <entryRelationship typeCode='RSON'>
4174 <act classCode='ACT' moodCode='EVN'>
4175 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.4.1' />
4176 <id root='' extension=''/>
4177 </act>
4178 </entryRelationship>
4179 <!-- An optional entry relationship that provides prescription activity -->
4180 <entryRelationship typeCode='REFR'>
4181 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.7.3' />
4182 :
4183 </entryRelationship>
4184 <precondition>
4185 <criterion>
4186 <text><reference value=''/></text>
4187 </criterion>
4188 </precondition>
4189 </substanceAdministration>

```

This content module is derived from the Medication content module to specifically and more easily describe the necessary details of intravenous fluid administration. For the purpose of EDER and other profiles employing this content module, the table below identifies and describes the fields and constraints on whether or not they are required to be sent. The fields are listed in the order that they appear in the CDA XML content.

### 6.3.4.46.1.1 Medication Fields

Field	Opt.	CDA Tag	Description
Start and Stop Date	R2	<effectiveTime>	The date and time when the fluid regimen began and is expected to finish. The first component of the <effectiveTime> encodes the lower and upper bounds over which the <substanceAdministration> occurs,

Field	Opt.	CDA Tag	Description
			and the start time is determined from the lower bound. If the fluid has been known to be stopped, the high value must be present, but expressed as a flavor of null (e.g., Unknown).
Dose	R2	<doseQuantity>	The amount of fluid given. This should be in some known and measurable fluid unit, such as milliliters, or may be measured in "administration" units (such "units" of blood or "packs" of platelets).
Site	O	<approachSiteCode>	The site where the fluid is administered (i.e. "Left Antecubital", or "Central Line").
Rate	R2	<rateQuantity>	The rate is a measurement of how fast the fluid is given to the patient over time (e.g., .5 liter / 1 hr).
Product	R	<consumable> <name> </consumable>	The name of the substance or product. This should be sufficient for a provider to identify the type of fluid. It may be a trade name (Plasmalyte®) or a generic name. This information is required in all fluid entries. The name should not include packaging, strength or dosing information.
Code	R2	<consumable> <code> </consumable>	A code describing the product from a controlled vocabulary, such as RxNorm, First DataBank, et cetera.

#### **6.3.4.46.1.2 <substanceAdministration classCode='SBADM' moodCode='INT|EVN'>**

4200 The general model is to record each fluid administered in a <substanceAdministration> intent (moodCode='INT'). Fluids that have been started but not completely administered should be recorded in a <substanceAdministration> intent (moodCode='INT'). Fluids that have been completed should be recorded as an event (moodCode='EVN').

#### **6.3.4.46.1.3 <templateId root='2.16.840.1.113883.10.20.1.24'> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.7'> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.7.1' />**

4205 All intravenous fluid entries use the <templateId> elements specified above to indicate that they are IV fluid administration acts. This element is required.

#### **6.3.4.46.1.4 <id root=" extension="/">**

4210 The <substanceAdministration> element must be uniquely identified. If there is no explicit identifier for this observation in the source EMR system, a GUID may be used for the root attribute, and the extension may be omitted. Although HL7 allows for multiple identifiers, this Technical Framework profile requires that one and only one be used.

#### **6.3.4.46.1.5 <code code="" displayName="" codeSystem="" codeSystemName="">**

4215 The <code> element is required, and is used to supply a code that describes the act of fluid administration, not the fluid being administered. This may be a procedure code, such as those found in CPT-4 (and often used for billing), or may describe the method of administration, such as by intravenous injection.

**6.3.4.46.1.6 <text><reference value="/"></text>**

4220 The URI given in the value attribute of the <reference> element points to an element in the narrative content that contains the complete text describing the fluid administration.

**6.3.4.46.1.7 <statusCode code='completed|active' />**

The status of all <substanceAdministration> elements must be "completed" or "active". If "completed", then the administration has occurred, or the request or order has been placed. If "active", then at the time recorded, the fluid was still being administered.

4225 **6.3.4.1.46.8 <effectiveTime xsi:type='IVL\_TS'>**

The first <effectiveTime> element encodes the start and stop time of the administration. This is an interval of time (xsi:type='IVL\_TS'), and must be specified as shown. This is an additional constraint placed upon CDA Release 2.0 by this Technical Framework profile, and simplifies the exchange of start/stop and frequency information between EMR systems.

4230 **6.3.4.46.1.9 <low value="/"><high value="/">**

The <low> and <high> values of the first <effectiveTime> element represent the start and stop times for the fluid administration. The <low> value represents the start time, and the <high> value represents the stop time. If either the <low> or the <high> value is unknown, this shall be recorded by setting the nullFlavor attribute to UNK. The <high> value records the end of the fluid administration according to the information provided in the initial fluid order or RN documentation. For example, if the fluid order is for one liter, and the fluid is to be delivered at 250 mL/hr, then the high value should contain a datetime that is 4 hours later than the <low> value. The rationale is that a provider, seeing a discontinued fluid could normally assume that the fluid has been stopped, even if the intent of the treatment plan is to continue the fluid continuously.

4240

**6.3.4.46.1.10 <approachSiteCode code="" codeSystem="">  
originalText><reference value="/"></originalText>  
</approachSiteCode>**

4245

The <approachSiteCode> element contains a URI in the value attribute of the <reference> that points to the text in the narrative identifying the site. It may be coded to a controlled vocabulary that lists such sites (e.g., SNOMED-CT).

**6.3.4.46.1.11 <doseQuantity><low value="" unit=""><high value="" unit="">  
</doseQuantity>**

4250

The dose is specified if the <doseQuantity> element. If a dose range is given (e.g., 125-250 mL/hr [i.e. to replace fluid losses]), then the <low> and <high> bounds are specified in their respective elements, otherwise both <low> and <high> have the same value. The unit attribute should be derived from the HL7 UnitsOfMeasureCaseSensitive vocabulary .

**6.3.4.46.1.12 <low|high value=""> <translation> <originalText><reference value="/"></originalText> </translation></low|high >**

- 4255 Any <low> and <high> elements used for <doseQuantity> or <rateQuantity> should contain a <translation> element that provides a <reference> to the <originalText> found in the narrative body of the document .

**6.3.4.46.1.13 <rateQuantity><low value="" unit="/"><high value="" unit="/"></rateQuantity>**

- 4260 The rate is specified in the <rateQuantity> element. The rate is given in units that have measure over time. In this case, the units should be specified as a string made up of a unit of measure (see doseQuantity above), followed by a slash (/), followed by a time unit (s, min, h or d) (i.e. mL/hr). Again, if a range is given, then the <low> and <high> elements contain the lower and upper bound of the range, otherwise, they contain the same value.

- 4265 **6.3.4.46.1.14 <consumable>**

The <consumable> element shall be present, and shall contain a <manufacturedProduct> entry conforming to the Product Entry template (see PCC TF-2: 6.3.4.19).

*Add Section 6.3.4.47*

- 4270 **6.3.4.47 Nursing Assessments Battery 1.3.6.1.4.1.19376.1.5.3.1.1.13.3.4**

This entry describes a single row in the Nursing Assessment flowsheet. The single observation date/time and provider is applied to all other observations.

- 4275

- 4280

### 6.3.4.47.1 Specification

```

4285 <entry>
 <organizer classCode='BATTERY' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.4' />
 <id root=' ' extension=' '/>
 <code code='XX-ASSESS' displayName='Nursing Assessments Battery'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <statusCode code='completed' />
 <author>
 <time value=' ' />
 <assignedAuthor>
 <id root=' ' extension=' '/>
 </assignedAuthor>
 </author>
 <component>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
 </observation>
 </component>
 <component>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
 </observation>
 </component>
 :
 </organizer>
 </entry>

```

#### 6.3.4.47.1.1 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.13.3.4' />

The <templateId> element specifies that this organizer entry conforms to the Nursing Interventions battery. The root attribute SHALL contain the value "1.3.6.1.4.1.19376.1.5.3.1.1.13.3.4"

#### 6.3.4.47.1.2 <organizer classCode='BATTERY' moodCode='EVN'>

Each row in the Nursing Interventions battery SHALL be represented by an organizer with the classCode of 'BATTERY' and the moodCode of 'EVN'

#### 4320 6.3.4.47.1.3 <id root=' ' extension=' '/>

Each battery SHALL have a globally unique identifier.

#### 6.3.4.47.1.4 <code code='X-ASSESS' codeSystem='2.16.840.1.113883.6.1' />

The <code> element specifies the LOINC code that represents the content of the battery. The codeSystem attribute SHALL contain the value '2.16.840.1.113883.6.1'. The code attribute SHALL contain the value='X-ASSESS'. It is good practice to include displayName and codeSystemName for clarity and debugging. The corresponding values are 'Nursing Assessments battery' and 'LOINC' respectively.

**6.3.4.47.1.5 <author><time><assignedAuthor><id></assignedAuthor></author>**

- 4330 The <author> relation element points at the author that records the visit battery. This assignedAuthor may be different than the author of the document. The time element is used to record when the assigned author recorded the battery.

**6.3.4.47.1.6 <statusCode code='completed'>**

The status code for all batteries SHALL be 'completed'

**6.3.4.47.1.7 <component>**

- 4335 The battery is made of several component Simple Observations (see PCC TF-2: 6.3.4.20). The following table lists the allowable LOINC codes, displayNames, and observation types, and unit of measures for these observations.

<b>LOINC Code</b>	<b>displayName</b>	<b>xsi:type</b>	<b>value set</b>
9269-2	GLASGOW COMA CORE.TOTAL	CO	3..15
9268-4	GLASGOW COMA SCORE.MOTOR	CO	1..6
11454-6	LEVEL OF RESPONSIVENESS	CO	ALERT VERBAL RESPONSE PAINFUL RESPONSE UNRESPONSIVE
38208-5	PAIN SEVERITY	CO	0-10
48767-8	(COMMENT FIELD)	ED	

- 4340 Add Section 6.3.4.48

**6.3.4.48 Antenatal Testing and Surveillance Battery****1.3.6.1.4.1.19376.1.5.3.1.1.21.3.10**

This entry describes a single row in the Antenatal Testing and Surveillance Section. The single observation date/time and provider is applied to all other observations.

4345 **6.3.4.48.1 Specification**

```

<entry>
 <organizer classCode='BATTERY' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.3.10' />
 <id root=' ' extension=' ' />
 <code code='XX-ANTENATALTESTINGBATTERY' displayName='ANTENATAL TESTING AND SURVEILLANCE
BATTERY'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <statusCode code='completed' />
 <author>
 <time value=' ' />
 <assignedAuthor>
 <id root=' ' extension=' ' />
 </assignedAuthor>
 </author>
 <component>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
 </observation>
 </component>
 <component>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
 </observation>
 </component>
 :
 </organizer>
</entry>
```

**6.3.4.48.1.1 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.21.3.10' />**

The <templateId> element specifies that this organizer entry conforms to the Antenatal Testing and Surveillance Battery. The root attribute SHALL contain the value "1.3.6.1.4.1.19376.1.5.3.1.1.21.3.10"

4380 **6.3.4.48.1.2 <organizer classCode='BATTERY' moodCode='EVN'>**

Each row in the Antenatal Testing and Surveillance Battery SHALL be represented by an organizer with the classCode of 'BATTERY' and the moodCode of 'EVN'

**6.3.4.48.1.3 <id root=' ' extension=' ' />**

Each battery SHALL have a globally unique identifier.

4385 **6.3.4.48.1.4 <code code='XX- XX-ANTENATALTESTINGBATTERY'
codeSystem='2.16.840.1.113883.6.1' />**

The <code> element specifies the LOINC code that represents the content of the battery. The codeSystem attribute SHALL contain the value '2.16.840.1.113883.6.1'. The code attribute SHALL contain the value='XX-ANTENATALTESTINGBATTERY'. It is good practice to include displayName and codeSystemName for clarity and debugging. The corresponding values are 'ANTENATAL TESTING AND SURVEILLANCE BATTERY' and 'LOINC' respectively.

**6.3.4.48.1.5 <author><time><assignedAuthor><id></assignedAuthor></author>**

4395 The <author> relation element points at the author that records the visit battery. This assignedAuthor may be different than the author of the document. The time element is used to record when the assigned author recorded the battery.

**6.3.4.48.1.6 <statusCode code='completed'>**

The status code for all batteries SHALL be 'completed'

**6.3.4.48.1.7 <component>**

4400 The battery is made of several component Simple Observations (see PCC TF-2: 6.3.4.20). The following table lists the allowable LOINC codes, displayNames, and observation types, and unit of measures for these observations.

<b>LOINC Code</b>	<b>displayName</b>	<b>xsi:type</b>
11630-1	Biophysical profile.amniotic fluid volume	ED
11631-9	Biophysical profile.body movement	ED
11632-7	Biophysical profile.breathing movement	ED
11633-5	Biophysical profile.heart rate reactivity	ED
11635-0	Biophysical profile.tone	ED
11634-3	Biophysical profile.sum	ED
<b>35096-7</b>	<b>Ultrasound morphologic</b>	<b>ED</b>
<b>49086-2</b>	<b>Nuchal translucency screening</b>	<b>ED</b>
<b>51659-1</b>	<b>Hbs1 Antigen</b>	<b>ED</b>

Add Section 6.3.4.49

**4405 6.3.4.49 Immunization Recommendation 1.3.6.1.4.1.19376.1.5.3.1.4.12.2**

An immunization recommendation entry is used to record intended or proposed immunization activities. Proposed activities are suggestions for care or treatment that are used in decision making (these might appear as an input to, or output from clinical decision support). Intended activities describe the currently accepted plan, and are part of the care activities expected to occur for the patient.

```

<substanceAdministrationIntent typeCode='SBADM' moodCode='INT|PRP' negationInd='true|false'>
 <templateId root='2.16.840.1.113883.10.20.1.25'/>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.12.2'/>
 <id root='' extension=''/>
 <code code='IMMUNIZ' codeSystem='2.16.840.1.113883.5.4' codeSystemName='ActCode' />
 <text><reference value="#xxx"/></text>
 <statusCode code='completed'/>
 <effectiveTime><low value=''/><high value=''/></effectiveTime>
 <routeCode code='' codeSystem='' codeSystemName='RouteOfAdministration' />
 <approachSiteCode code='' codeSystem='' codeSystemName='HumanSubstanceAdministrationSite' />
 <doseQuantity value=' ' units=''/>
 <consumable typeCode='CSM'>
 :
 .
 </consumable>
 <!-- An optional entry relationship that identifies the immunization series number -->
 <entryRelationship typeCode='SUBJ'>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='2.16.840.1.113883.10.20.1.46'/>
 <code code='30973-2' displayName='Dose Number'
 codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <statusCode code='completed'/>
 <value xsi:type='INT' value=''/>
 </observation>
 </entryRelationship>
 <!-- Optional <entryRelationship> element containing comments -->
</substanceAdministrationIntent>

```

- 4440 Note: The CCD represents the observation of a series number in EVN mood, as we have shown above. However, when the immunization is “intended” to be the second of a series, we do not believe this is the correct mood code. How should this be addressed?

#### **6.3.4.49.1 <substanceAdministrationIntent typeCode='SBADM' moodCode='INT|PRP' negationInd='true|false'>**

- 4445 This entry represents the intent or proposal to administer (when negationInd = false), or not administer (when negationInd = true) an immunization to a patient.

#### **6.3.4.49.2 <templateId root='2.16.840.1.113883.10.20.1.25'>**

This element represents a plan of care activity for the patient, and so shall conform to the CCD Plan of Care activity template.

#### **6.3.4.49.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.12.2'>**

- 4450 This element is an instance of an IHE PCC Immunization Recommendation entry, and shall indicate that conformance by inclusion of the template identifier given above.

#### **6.3.4.49.4 <id root=" extension="/">**

Each plan of care activity shall contain an identifier.

**6.3.4.49.5 <code code='IMMUNIZ' codeSystem='2.16.840.1.113883.5.4'  
codeSystemName='ActCode'>**

4455 An immunization care plan entry shall include a code identifying that this is in reference to an immunization. If no coding system is required by the source, simply record as shown above. See PCC TF-2:6.3.4.17 for additional coding systems that may appear in this element.

**6.3.4.49.6 <text><reference value="#xxx"/></text>**

4460 In a CDA document, the URI given in the value attribute of the <reference> element points to an element in the narrative content that contains the complete text describing the immunization activity. In an HL7 message, the content of the text element shall contain the complete text describing the immunization activity.

**6.3.4.49.7 <statusCode code='active'>**

4465 The status code shall be active for all active proposals or intentions.

**6.3.4.49.8 <effectiveTime><low value="" /><high value="" /></effectiveTime>**

The <effectiveTime> element should be present to indicate time interval over which the suggested activity should take place. Intervals shall be represented using the IVL\_TS data type.

**6.3.4.49.9 <routeCode code=" codeSystem="  
codeSystemName='RouteOfAdministration' />  
<approachSiteCode code=" codeSystem="  
codeSystemName='HumanSubstanceAdministrationSite' />  
<doseQuantity value=" units="" />**

4470 The <routeCode>, <approachSiteCode> and <doseQuantity> elements are used to represent additional attributes of the proposed care. When present these elements must be consistent with the rules for these elements specified in PCC TF-2:6.3.4.16 Medication Entry and PCC TF-2:6.3.4.17 Immunization Entry.

**6.3.4.49.10 <consumable typeCode='CSM'>**

4480 The <consumable> element shall be present, and shall contain a <manufacturedProduct> entry conforming to the Product Entry template found in PCC TF-2:6.14.19.

**6.3.4.49.11 <entryRelationship typeCode='SUBJ'>  
<observation classCode='OBS' moodCode='EVN'>  
<templateId root='2.16.840.1.113883.10.20.1.46' />  
<code code='30973-2'  
displayName='Dose Number'  
codeSystem='2.16.840.1.113883.6.1'  
codeSystemName='LOINC' />  
<statusCode code='completed' />**

4490                    **<value xsi:type='INT' value=""/>**  
                          **</observation>**  
                          **</entryRelationship>**

4495     The immunization plan of care entry may contain a single entry relationship identifying the immunization series number. This entry shall use the CCD template (2.16.840.1.113883.10.20.1.46) defined for that purpose.

#### **6.3.4.49.12 <!-- Optional <entryRelationship> element referencing guidelines -->**

*Add Section 6.3.4.50*

#### **6.3.4.50 Alert Entry 1.3.6.1.4.1.19376.1.5.3.1.4.12.3**

4500     The alert entry is an observation whose subject is any clinical statement. This entry provides additional information about a clinical statement that may be of relevance to the care being described. For example, some treatments may be contraindicated by other conditions or co-occurring treatment. For example, the use of aspirin and a blood thinning agent at the same time may not be recommended. The alert entry is provided to record these annotations. An example use of this entry is in a clinical decision support service that uses the alert entry to identify vaccinations that are considered to be of reduced effectiveness when making immunization recommendations for the patient. Another example of this use in a similar system would be to identify treatments that are contraindicated subsequent to an immunization.

4510                   **<entryRelationship typeCode='SUBJ' inversionInd='true'>**  
                          **<templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.12.3' />**  
                          **<observation classCode='ALRT|OBS' moodCode='EVN'>**  
                           **<id root='' id=''/>**  
                           **<code code='' displayName='' codeSystem='' codeSystemName='' />**  
                           **<text><reference value='#ref-1' /></text>**  
                          **</observation>**  
                          **</entryRelationship>**

#### **6.3.4.50.1 <entryRelationship typeCode='SUBJ' inversionInd='true'>**

The alert has a preexisting entry as its subject (typeCode=SUBJ).

4520     **6.3.4.50.2 <templateId root='TBD' />**

This alert complies with the rules specified in the PCC technical framework for alerts, and so must include the templateId specified above.

#### **6.3.4.50.3 <observation classCode='OBS|ALRT' moodCode='EVN'>**

4525     An alert is an observation that has occurred (moodCode=EVN). The HL7 classCode value of ALRT shall be used where permitted (e.g. in an HL7 Care Record message). Where not permitted, the classCode shall be OBS (e.g., in CDA Document).

**6.3.4.50.4 <id root="extension"/>**

Each alert observation may have an identifier.

**6.3.4.50.5 <code code=" displayName=" codeSystem='2.16.840.1.113883.5.4' codeSystemName='ActDetectedIssueCode'/>**

4530

Each alert observation shall have a code identifying the type of issue detected. The HL7 ActDetectedIssueCode value set (2.16.840.1.113883.1.11.16124) is one possible source of codes for these issues.

**6.3.4.50.6 <text><reference value="#ref-1"/></text>**

4535

The text of the observation should provide some human readable explanation for the alert. In a CDA document, this would appear within the narrative of the clinical document, and so would be referenced by the alert. In an HL7 Version 3 message, this text would appear in the **<text>** element of the alert entry.

4540

*Add Section 6.3.4.51*

**6.3.4.51 Antigen Dose 1.3.6.1.4.1.19376.1.5.3.1.4.12.1**

An Antigen Dose entry is used to record additional details about the patient's immunization history. These entries may be used to provide dose details about a specific antigen received during an Immunization.

4545

**6.3.4.51.1 Specification**

4550

```

<substanceAdministration typeCode='SBADM' moodCode='EVN' negationInd='false'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.12.1'>
 <id root='' extension=''/>
 <code code='IMMUNIZ' codeSystem='2.16.840.1.113883.5.4' codeSystemName='ActCode' />
 <text><reference value='#xxxx' /></text>
 <statusCode code='completed' />
 <effectiveTime value=''/>
 <routeCode code=' ' codeSystem=' ' codeSystemName='RouteOfAdministration' />
 <approachSiteCode code=' ' codeSystem=' '
 codeSystemName='HumanSubstanceAdministrationSite' />
 <doseQuantity value=' ' units=''/>
 <consumable typeCode='CSM'>
 <manufacturedProduct classCode='MANU'>
 <manufacturedLabeledDrug classCode='MMAT' determinerCode='KIND'>
 <code code=' ' codeSystem=' ' codeSystemName=' '>
 <originalText><reference value='#yyy' /></originalText>
 </code>
 </manufacturedLabeledDrug>
 </manufacturedProduct>
 </consumable>
 </substanceAdministration>

```

4555

4560

4565

**Figure 6.3.4.51.1-1: Immunizations Example**

**6.3.4.51.2 <substanceAdministration typeCode='SBADM' moodCode='EVN' negationInd='false'>**

4570 An antigen dose entry is a substance administration event.

**6.3.4.51.3 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.12.1'>**

The <templateId> element identifies this <substanceAdministration> as an antigen dose, allowing for validation of the content. The templateId must have root='1.3.6.1.4.1.19376.1.5.3.1.4.12.1'.

4575 **6.3.4.51.4 <id root=' ' extension=' '/>**

This shall be the identifier for the antigen dose event.

**6.3.4.51.5 <code/><text/><statusCode/><effectiveTime value=' '/>  
<routeCode code=' ' codeSystem=' ' codeSystemName='RouteOfAdministration'/>  
<approachSiteCode code=' ' codeSystem=' '  
codeSystemName='HumanSubstanceAdministrationSite'/>**

4580

Since the antigen dose entry is subordinate to an immunization entry, the code, text, status code, effective time, route and approach site would all be repetitive and therefore should not be present.

**6.3.4.51.6 <doseQuantity value=' ' units=' '/>**

4585 This element gives the dose quantity of the specific antigen.

**6.3.4.51.7 <consumable typeCode='CSM'>**

The <consumable> element shall be present, and shall contain a <manufacturedProduct> entry conforming to the Product Entry template found in PCC TF-2:6.3.4.19. This product entry describes the antigen to which the dose is applied.

4590

*Add Section 6.3.4.52 (Occupation Observation – removed 2011-09 at the request of QRPH)*

**6.3.4.52 Intentionally blank**

*Add Section 6.3.4.53 (Industry Observation removed 2011-09 at the request of QRPH)*

4595 **6.3.4.53 Intentionally blank**

*Add Section 6.3.4.54*

### 6.3.4.54 Observation Request 1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1

4600 The observation request entry is used to record goals, plans or intention for an observation to be performed (e.g., assessment, laboratory test, imaging study, et cetera).

#### 6.3.4.54.1 Uses

See Templates using Observation Request

#### 6.3.4.54.2 Specification

4605

```

<observation classCode='OBS' moodCode='INT|PRP|GOL'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1' />
 <templateId root='2.16.840.1.113883.10.20.1.25' />
 <id root='' extension='' />
 <code code='' displayName='' codeSystem='' codeSystemName='' />
 <!-- for CDA -->
 <text><reference value='#xxx' /></text>
 <!-- For HL7 Version 3 Messages
 <text>text</text>
 -->
 <statusCode code='active' />
 <effectiveTime value='' />
 <repeatNumber value='' />
 <value xsi:type='' ... />
 <interpretationCode code='' codeSystem='' codeSystemName='' />
 <methodCode code='' codeSystem='' codeSystemName='' />
 <targetSiteCode code='' codeSystem='' codeSystemName='' />
 <author typeCode='AUT'>
 <assignedAuthor typeCode='ASSIGNED'><id ... /></assignedAuthor> <!-- for CDA --
 >
 <!-- For HL7 Version 3 Messages
 <assignedEntity typeCode='ASSIGNED'>
 <Person classCode='PSN'>
 <determinerCode root='' />
 <name>...</name>
 </Person>
 <assignedEntity>
 -->
 </assignedEntity>
 </author>
 </observation>
```

**Figure 6.3.4.54.1-1: Observation Request Example**

#### 6.3.4.54.2.1 <observation classCode='OBS' moodCode='INT|PRP|GOL'>

These acts are observations that form the care plan or which can be used in decision support. In intent mood (moodCode='INT') these are what is intended to be performed as part of the care plan. In proposal mood (moodCode='PRP'), these observations are being proposed, for example, as the output of a clinical decision support system. In goal mood (moodCode='GOL'), these observations described the intended goal of a treatment plan.

**6.3.4.54.2.2 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.1.20.3.1'>**

4645 The <templateId> element identifies this <observation> as an observation request, allowing for validation of the content. The templateId must appear as shown above.

**6.3.4.54.2.3 <templateId root=2.16.840.1.113883.10.20.1.25'>**

The IHE Observation Request template conforms to the Plan of care activity defined by the HL7 Continuity of Care Document. This template id must be present to indicate conformance.

**6.3.4.54.2.4 <id root=' ' extension=' '/>**

4650 Each observation shall have an identifier.

**6.3.4.54.2.5 <code code=' ' displayName=' ' codeSystem=' ' codeSystemName=' '/>**

4655 Observations shall have a code describing what is to be measured. The code system used is determined by the vocabulary constraints on the types of measurements that might be recorded in a section. Modules that are derived from this one may restrict the code system and code values used for the observation.

**6.3.4.54.2.6 <text><reference value="#xxx"/></text> -OR- <text>text</text>**

4660 Each observation request entry may contain a <text> element providing the free text that provides the same information as the observation within the narrative portion of the document with a <text> element. For CDA based uses of Observation Requests, this element SHALL be present, and SHALL contain a <reference> element that points to the related string in the narrative portion of the document. For HL7 Version 3 based uses, the <text> element MAY be included.

**6.3.4.54.2.7 <statusCode code='active'>**

4665 The <statusCode> element shall be present and shall describe the current state of the observation. Goals, intents and proposals that are available for action shall have an 'active' status, but other status values are permitted.

**6.3.4.54.2.8 <effectiveTime value=' '/>**

The <effectiveTime> element shall be present in observation requests to indicate the date and time when the measurement should be taken.

4670 **6.3.4.54.2.9 <value xsi:type=' ' .../>**

The value of the observation may be recorded using a data type appropriate to the observation to indicate the desired value (e.g., in GOL or PRP mood).

**6.3.4.54.2.10 <methodCode code=' ' codeSystem=' ' codeSystemName=' '/>**

4675 The methodCode element may be used to record the specific method used to make an observation when this information is not already pre-coordinated with the observation code.

#### **6.3.4.54.2.11 <targetSiteCode code=' ' codeSystem=' ' codeSystemName=' '/>**

The targetSiteCode may be used to record the target site where the observation should be made when this information is not already pre-coordinated with the observation code.

#### **6.3.4.54.2.12 <author><assignedAuthor**

4680 **classCode='ASSIGNED'>...<assignedAuthor></author>**

In CDA uses, the observation request is assumed to be authored by the same author as the document through context conduction. However, observation requests would often be used to record orders, and in these cases, the author of the order shall be recorded in the author element.

4685 For HL7 Version 3 purposes, the <author> element SHOULD be present unless it can be

determined by conduction from organizers or higher level structures. When used for HL7 Version 3 the role element name is <assignedEntity> and the author is represented as <assignedPerson> element.

*Add Section 6.3.4.55 (Added 2011-09 from QRPH EHCP profile)*

4690 **6.3.4.55 Risk Indicators for Hearing Loss Entry 1.3.6.1.4.1.19376.1.7.3.1.1.15.5.1**

This entry describes the Risk Indicators for Hearing Loss.

#### **6.3.4.55.1 Specification**

4695

```

<entry>
 <organizer classCode='BATTERY' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.15.5.1' />
 <id root=' ' extension=' ' />
 <code code='58232-0' displayName='Hearing Loss Risk Indicators' codeSystem='2.16.840.1.113883.6.1' codeSystemName='LOINC' />
 <statusCode code='completed' />
 <author>
 <time value=' ' />
 <assignedAuthor>
 <id root=' ' extension=' ' />
 </assignedAuthor>
 </author>
 <component>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
 </observation>
 </component>
 <component>
 <observation classCode='OBS' moodCode='EVN'>
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.13' />
 :
 </observation>
 </component>
 :
 </organizer>
</entry>

```

**Figure 6.3.4.93.1-1: Sample Risk Indicators for Hearing Loss Entry****6.3.4.55.2 <templateId root='1.3.6.1.4.1.19376.1.7.3.1.1.15.5.1' />**

- 4700 The <templateId> element specifies that this organizer entry conforms to the Nursing Interventions battery. The root attribute SHALL contain the value "1.3.6.1.4.1.19376.1.7.3.1.1.15.5.1"

**6.3.4.55.3 <organizer classCode='BATTERY' moodCode='EVN'>**

- 4705 Each row in the Nursing Interventions battery SHALL be represented by an organizer with the classCode of 'BATTERY' and the moodCode of 'EVN'

**6.3.4.55.4 <id root=' ' extension=' ' />**

Each battery SHALL have a globally unique identifier.

**6.3.4.55.5 <code code='58232-0' codeSystem='2.16.840.1.113883.6.1'>**

- 4710 The <code> element specifies the LOINC® code that represents the content of the battery. The codeSystem attribute SHALL contain the value '2.16.840.1.113883.6.1'. The code attribute SHALL contain the value='58232-0'. It is good practice to include displayName and codeSystemName for clarity and debugging. The corresponding values are 'Hearing Loss Risk Indicators' and 'LOINC®' respectively.

**6.3.4.55.6 <author/><time/><assignedAuthor><id/></assignedAuthor></author>**

- 4715 The <author> relation element points at the author that records the visit battery. This assignedAuthor MAY be different than the author of the document. The time element is used to record when the assigned author recorded the battery.

**6.3.4.55.7 <statusCode code='completed'>**

The status code for all batteries SHALL be 'completed'

4720 **6.4.4.55.8 <component>**

The battery is made of several component Simple Observations. The observation values SHALL be constrained to those coded values and descriptions described by the JCIH-EHDI Risk Indicators for Hearing Loss (LOINC®) Value Set (1.3.6.1.4.1.19376.1.7.3.1.1.15.2.24).

4725 *Add Section 6.3.4.56. (Added 2011-09 from QRPH PRPH-Ca profile)*

**6.3.4.56 Cancer Diagnosis Entry 1.3.6.1.4.1.19376.1.7.3.1.4.14.1**

A Cancer Diagnosis entry collects details of the patient's cancer diagnosis, including histology, behavior, primary site, laterality, diagnosis date, TNM Stage, and Best Method of Confirmation.

**6.3.4.56.1 Parent Template**

- 4730 The parent of this template is Problem Concern Entry (1.3.6.1.4.1.19376.1.5.3.1.4.5.2).

### 6.3.4.56.2 Specification

```

<section>
 <templateId root="2.16.840.1.113883.10.20.1.11"/>
 <templateId root="1.3.6.1.4.1.19376.1.5.3.1.3.6"/>
 <templateId root="1.3.6.1.4.1.19376.1.7.3.1.3.14.1"/>
 <title>"Cancer Diagnosis"</title>
 <text>"Malignant melanoma of the left leg, Stage 1"</text>
 <entry>
 <act classCode='ACT' moodCode='EVN'>
 <templateId root='2.16.840.1.113883.10.20.1.27' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5.2' />
 <code nullFlavor='NA' />
 <statusCode code='active' />
 <effectiveTime>
 <low value='20110101' />
 <high nullFlavor="NA" />
 </effectiveTime>
 <entryRelationship typeCode="SUBJ" inversionInd="false" >
 <observation classCode='OBS' moodCode='EVN' negationInd="false">
 <templateId root='2.16.840.1.113883.10.20.1.28' />
 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5.' />
 <templateId root="1.3.6.1.4.1.19376.1.7.3.1.4.14.1" />
 <code code="282291009" codeSystem="2.16.840.1.113883.6.96"
codeSystemName="SNOMED CT" displayName="Diagnosis"/>
 <text><reference value="" ></reference></text>
 <statusCode code="completed" />
 <effectiveTime>
 <low value="20110101" />
 <high nullFlavor="NI" />
 </effectiveTime>
 <!--The <value> is the condition that was found.-->
 <value xsi:type="CD" code="8742" codeSystem="2.16.840.1.113883.3.520.3.2"
codeSystemName="NAACCR Histologic Type" displayName="Lentigo Maligna" >
 <!--Behavior Qualifier-->
 <qualifier>
 <name code="31206-6" codeSystem="2.16.840.1.113883.6.1"
codeSystemName="LOINC" displayName="Behavior ICD-O-3" />
 <value code="2" codeSystem="2.16.840.1.113883.3.520.3.14"
codeSystemName="NAACCR Behavior Code" displayName="In Situ" />
 </qualifier>
 <qualifier>
 <!--Best Method of Diagnosis Qualifier-->
 <name code="21861-0" codeSystem="2.16.840.1.113883.6.1"
codeSystemName="LOINC" displayName="Diagnostic Confirmation" />
 <value xsi:type="CD" code="2"
codeSystem="2.16.840.1.113883.3.520.3.3" codeSystemName="NAACCR Diagnostic Confirmation"
displayName="Positive cytology, no positive histology" />
 </qualifier>
 </value>
 <!--Primary Site -->
 <targetSiteCode code="C447" codeSystem="2.16.840.1.113883.6.43.1"
codeSystemName="ICD-O-3 (Topography Section)" displayName="Leg" >
 <!--Laterality-->
 <qualifier>
 <name code="20228-3" codeSystem="2.16.840.1.113883.6.1"
codeSystemName="LOINC" displayName="Anatomic part Laterality" />
 <value code="1" codeSystem="2.16.840.1.113883.3.520.3.1"
codeSystemName="NAACCR Laterality at Diagnosis" displayName="origin of primary: right" />
 </qualifier>
 </targetSiteCode>
 <entryRelationship typeCode="SUBJ" inversionInd="true" >
 <!--TNM Stage Information-->
 </entry>
 </section>

```

```

 <observation classCode="OBS" moodCode="EVN">
 <templateId
root="1.3.6.1.4.1.19376.1.7.3.1.4.14.2"/>
 <code code="21908-9" displayName="TNM Clinical Stage
Group" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/>
<!-- Narrative TNM Clinical Stage -->
 <text> Stage 0 TisN0M0 </text>
 <statusCode code="completed"/>
 <value xsi:type="CD" code="0"
codeSystem="2.16.840.1.113883.3.520.3.9" codeSystemName="NAACCR TNM Clinical Stage Group"
displayName="In Situ">
 <qualifier>
<!--TNM Clinical Stage Descriptor Observation -->
 <name code="21909-7" displayName="TNM
Clinical Stage Descriptor" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/>
 <value xsi:type="CD" code="0"
codeSystem="2.16.840.1.113883.3.520.3.10" codeSystemName="NAACCR TNM Clinical Stage Descriptor"
displayName="None"/>
 </qualifier>
<!--AJCC TNM Edition Number.-->
 <qualifier>
 <name code="21917-0" displayName="TNM
Edition Number" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/>
 <value xsi:type="CD" code="7"
codeSystem="2.16.840.1.113883.3.520.3.5" codeSystemName="NAACCR TNM Edition Number"
displayName="7th Edition"/>
 </qualifier>
 </value>
 <participant typeCode="PPRF">
 <participantRole>
 <code code="21910-5"
codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" displayName="Stager.clinical Cancer"/>
 <playingEntity nullFlavor="NA">
 <code xsi:type="CE" code="1"
codeSystem="2.16.840.1.113883.3.520.3.4" codeSystemName="TNM Clinical Staged By"
displayName="Managing Physician"/>
 </playingEntity>
 </participantRole>
 </participant>
 <entryRelationship typeCode="COMP">
<!-- 6.3.4.62 TNM Clinical Tumor Observation-->
 <observation classCode="OBS" moodCode="EVN">
 <templateId root="1.3.6.1.4.1.19376.1.5.3.1.4.13"/>
 <code code="21905-5" displayName="TNM Clinical T"
codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/>
 <statusCode code="completed"/>
 <value xsi:type="CD" code="Tis"
codeSystem="2.16.840.1.113883.3.520.3.6" codeSystemName="NAACCR TNM Clinical Tumor"
displayName="In Situ"/>
 </observation>
 </entryRelationship>
<!--6.3.4.63 TNM Clinical Nodes Observation -->
 <entryRelationship typeCode="COMP">
 <observation classCode="OBS" moodCode="EVN">
 <templateId root="1.3.6.1.4.1.19376.1.5.3.1.4.13"/>
 <code code="21906-3" displayName="TNM Clinical N"
codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/>
 <statusCode code="completed"/>
 <value xsi:type="CD" code="N0"
codeSystem="2.16.840.1.113883.3.520.3.7" codeSystemName="NAACCR TNM Clinical Nodes"
displayName="None"/>
 </observation>
 </entryRelationship>

```

```

4865 <!--6.3.4.64 TNM Clinical Metastases Observation-->
4866 <entryRelationship typeCode="COMP">
4867 <observation classCode="OBS" moodCode="EVN">
4868 <templateId root="1.3.6.1.4.1.19376.1.5.3.1.4.13"/>
4869 <code code="21907-1" displayName="TNM Clinical M"
4870 codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/>
4871 <statusCode code="completed"/>
4872 <value xsi:type="CD"
4873 codeSystem="2.16.840.1.113883.3.520.3.8" codeSystemName="NAACCR TNM Clinical Metastases"
4874 code="M0" displayName="None"/>
4875 </observation>
4876 </entryRelationship>
4877 </observation>
4878 </entryRelationship>
4879 </act>
4880 </entry>
4881 </section>

```

**Figure 6.3.4.94.1-1: Sample Cancer Diagnosis Entry****6.3.4.56.3 <act classCode='ACT' moodCode='EVN'>**

4885 All concerns reflect the act of recording (<act classCode='ACT'>) the event (moodCode='EVN') of being concerned about a problem, allergy or other issue about the patient condition.

**6.3.4.56.4 <templateId root='2.16.840.1.113883.10.20.1.27'/> <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5.2'/>**

4890 These template identifiers indicates this entry conforms to the concern content module. This content module inherits constraints from the HL7 CCD Template for problem acts, and so also includes that template identifier.

**6.3.4.56.5 <!-- 1..\* entry relationships identifying problems of concern --><entryRelationship type='SUBJ'><observation classCode='OBS' moodCode='EVN'><templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5'>...</observation>**

4895 >

This entry shall contain one or more problem entries that conform to the Problem Entry template 1.3.6.1.4.1.19376.1.5.3.1.4.5. The typeCode SHALL be “SUBJ” and inversionInd SHALL be “false”.

**6.3.4.56.6 <observation classCode="OBS" moodCode="EVN">**

4900 The <observation> classCode and moodCode SHALL be recorded as shown above.

**6.3.4.56.7 <templateId root='1.3.6.1.4.1.19376.1.5.3.1.4.5'> <templateId root='1.3.6.1.4.1.19376.1.7.3.1.4.14.1'>**

4905 These <templateId> elements identify this <entry> as a cancer diagnosis entry and its parent, Problem Entry, allowing for validation of the content. The <templateId> elements shall be recorded as shown above.

**6.3.4.56.8 <code code="282291009" codeSystem=" 2.16.840.1.113883.6.96" codeSystemName="SNOMED CT" displayName="Diagnosis"/>**

4910 The <code> element indicates that this is the Diagnosis information. This code SHALL be the SNOMED CT code “282291009” for “Diagnosis”. It is good style to include the displayName and codeSystemName to help debugging.

**6.3.4.56.9 <statusCode code='completed' />**

The status code for all Cancer Diagnosis Entries SHALL be ‘completed’.

**6.3.4.56.10 <effectiveTime value="xxx"/>**

4915 This element records the date of initial diagnosis by a recognized medical practitioner for the cancer being reported.

**6.3.4.56.11 <value xsi:type='CD' code=' ' codeSystem=' ' codeSystemName=' ' displayName=' '>**

4920 The <value> records the Histologic Type, which is the cell type of the tumor/cancer (e.g., carcinoma, melanoma, sarcoma, lymphoma, leukemia). This element is required. It is always represented using the CD datatype (xsi:type='CD'), even though the value may be a coded or uncoded string. If coded, it SHALL follow the appropriate realm constraints for vocabulary. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).

**6.3.4.56.12 <qualifier><name code="31206-6" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" displayName="Behavior ICD-O-3 Cancer"/><value code="" codeSystem="" codeSystemName="" displayName="" /> </qualifier>**

4925 This <qualifier> provides Behavior information, indicating whether the tumor is benign, in situ, malignant or metastatic. The code and codeSystem attributes SHALL be recorded exactly as shown above. If coded, it SHALL follow the appropriate realm constraints for vocabulary. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).

**6.3.4.56.13 <qualifier><name code="21861-0" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" displayName="Dx confirmed by Cancer"/><value xsi:type="CD" code="" codeSystem="" codeSystemName="" displayName="" /> </qualifier>**

4935 This <qualifier> provides Best Method of Diagnosis information, indicating the best method used to confirm the presence of the cancer being reported. The code and codeSystem attributes SHALL be recorded exactly as shown above. The <value> records the best method of diagnosis, and if coded, it SHALL follow the appropriate realm constraints for vocabulary. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).

4940 **6.3.4.56.14 <targetSiteCode code="" codeSystem="" codeSystemName="" displayName="">**

The <targetSiteCode> element SHALL be present and shall indicate the anatomic location where the primary tumor originated. Vocabulary used SHALL follow the appropriate realm constraints. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).

4945 **6.3.4.56.15 <qualifier><name code="20228-3" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" displayName="Anatomic part Laterality"/> <value code="" codeSystem="" codeSystemName="" displayName="/" /></qualifier>**

4950 This <qualifier> provides the laterality, which indicates the side of a paired organ or side of the body on which the reportable tumor originated. The code and codeSystem attributes SHALL be recorded exactly as shown above. The <value> records the laterality, if coded, it SHALL follow the appropriate realm constraints for vocabulary. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).

**6.3.4.56.16 <entryRelationship typeCode="SUBJ" inversionInd="false">**

4955 One <entryRelationship> element should be present providing information on the TNM Clinical Stage.

When present, this <entryRelationship> element SHALL contain an observation conforming to the TNM Stage Information (1.3.6.1.4.1.19376.1.7.3.1.4.14.2) template. The typeCode SHALL be “SUBJ” and inversionInd SHALL be “false”.

4960 **6.3.4.56.17 <observation classCode="OBS" moodCode="EVN"> <templateId root="1.3.6.1.4.1.19376.1.7.3.1.4.14.2"/> [1<sup>st</sup> nesting]**

Observations that describe the TNM Stage Information SHALL be included if known.

**6.3.4.56.18 <code code="xxxxx-x" displayName="TNM Clinical Stage Information" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/> [1<sup>st</sup> nesting]**

4965 The <code> element indicates that this observation is the TNM Clinical Stage Information. This code SHALL be the LOINC code xxxxx-x. It is good style to include the displayName and codeSystemName to help debugging.

**6.3.4.56.19 <statusCode code="completed"/> [1<sup>st</sup> nesting]**

The status code for all TNM Clinical Stage Information observations SHALL be ‘completed’.

4970 **6.3.4.56.20 <value xsi:type="CD" code="" codeSystem="" codeSystemName="" displayName=""> [1<sup>st</sup> nesting]**

The <value> records the TNM Clinical Stage Group, which is a detailed site-specific code for the clinical stage group as defined by AJCC and recorded by the physician. This element is required. It is always represented using the CD datatype (xsi:type='CD'), even though the value

- 4975 may be a coded or uncoded string. If coded, it SHALL follow the appropriate realm constraints for vocabulary. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).
- 6.3.4.56.21 <qualifier><name code="21909-7" displayName="Descriptor.clinical Cancer" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/><value xsi:type="CD" code="" codeSystem="" codeSystemName=" " displayName=" "/></qualifier> [1<sup>st</sup> nesting]**
- 4980 This <qualifier> provides TNM Clinical Stage Descriptor information, indicating The AJCC clinical stage prefix/suffix recorded by the physician. AJCC stage descriptors identify special cases that require separate analysis. The code and codeSystem attributes SHALL be recorded exactly as shown above. The <value> records the TNM Clinical Stage Descriptor, and if coded, it SHALL follow the appropriate realm constraints for vocabulary. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).
- 6.3.4.56.22 <qualifier><name code="21917-0" displayName="Version TNM Classification" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/><value xsi:type="CD" code="" codeSystem="" codeSystemName=" " displayName=" "/></qualifier> [1<sup>st</sup> nesting]**
- 4990 This <qualifier> provides TNM Edition Number information, indicating the edition number of the AJCC Staging Manual. The code and codeSystem attributes of <name> SHALL be recorded exactly as shown above. If coded, it SHALL follow the appropriate realm constraints for vocabulary. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).
- 5000 **6.3.4.56.23 <participant typeCode="PPRF"> <participantRole> <code code="21910-5" codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC" displayName="Stager.clinical Cancer"/><playingEntity nullFlavor="NA"> <code xsi:type="CE" code="" codeSystem="" codeSystemName=" " displayName=" "/> [1<sup>st</sup> nesting]**
- 5005 This <participant> element should specify the person who recorded the AJCC staging elements and stage group in the patient's medical record. The code and codeSystem attributes for <participantRole> SHALL be recorded exactly as shown above. The <code> attribute of <playingEntity> identifies the person who recorded the staging elements, and SHALL follow the appropriate realm constraints for vocabulary. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).
- 6.3.4.56.24 <!-- 0..3 entryRelationships identifying simple observations for TNM Clinic Tumor, TNM Clinical Nodes, and TNM Clinical Metastases--><entryRelationship typeCode="COMP" inversionInd="false"><observation classCode='OBS'moodCode='EVN'><templateIDroot='1.3.6.1.4.1.19376.1.5.3.1.4.13'>...</observation>[2nd nesting]**
- 5010 Each <entryRelationship> element should contain a simple observation that specifies the TNM Clinic Tumor, TNM Clinical Nodes, and TNM Clinical Metastases, each of which is a

5015 component of the TNM Stage Group. Simple observations that describe the TNM Clinic Tumor, TNM Clinical Nodes, and TNM Clinical Metastases SHALL be included if known and inversionInd SHALL be “false”.

#### **6.3.4.56.25 <code code="" displayName=" " codeSystem="2.16.840.1.113883.6.1" codeSystemName="LOINC"/> [2<sup>nd</sup> nesting]**

Observations SHALL include one of following LOINC values to indicate the component of TNM Stage Group represented in the Observation.

5020

LOINC Code	Display Name	Description
21905-5	TNM Clinical T	A detailed site-specific code for the clinical tumor (T) as defined by AJCC and recorded by the physician.
21906-3	TNM Clinical N	A detailed site-specific code for the clinical nodes (N) as defined by AJCC and recorded by the physician.
21907-1	TNM Clinical M	A detailed site-specific staging code for the clinical metastases (M) as defined by AJCC and recorded by the physician.

#### **6.3.4.56.26 <value xsi:type="CD" code="" codeSystem="" codeSystemName=" " displayName=" "/>**

5025

The <value> of the observation SHALL be recorded using the vocabulary appropriate to the coded observation according to the table above and SHALL follow the appropriate realm constraints for vocabulary. (See Volume 4 in the QRPH PRPH-Ca profile found at <http://www.ihe.net>).

#### **6.3.4.57 Patient Transfer TBD OID**

5030

The Patient Transfer entry shall record the transfer of the patient to an internal department or external entity such as a different hospital or skilled nursing facility.

#### **6.3.4.56.1 Parent Template**

The parent of this template is Problem Concern Entry (1.3.6.1.4.1.19376.1.5.3.1.4.5.2).

### 6.3.4.56.2 Specification

```

5035 <act classCode='ACT' moodCode='EVN'>
<templateId root='PatientTransferAct' />
<id/>
<!-- code is fixed -->
<code code='107724000' displayName='patient transfer' codeSystem='2.16.840.1.113883.6.96' />
5040 <effectiveTime value='' />
<participant typeCode='DST'>
 <templateId root='destinationLocation' />
 <participantRole classCode='SDLOC'>
 <id/>
 <code/>
 <addr/>
 <telecom/>
 <playingEntity classCode='ENT'>
 <name/>
 </playingEntity>
 </participantRole>
</participant>
</act>

```

**Figure 6.3.4.57.2-1: Sample Cancer Diagnosis Entry**

### 6.3.4.57.3 <act classCode='ACT' moodCode='INT|EVN'>

5055 The transfer is recorded in an act element, to describe a patient transfer. In intent mood (moodCode='INT'), this records the expected transfer of the patient. In event mood (moodCode='EVN'), this records the actual transfer.

### 6.3.4.57.4 <templateId root='TBD' />

5060 The templateId indicates that this transfer entry conforms to the constraints of this content module.

### 6.3.4.57.5 <id root=" extension="/" />

This required element shall contain an identifier.

### 6.3.4.57.6 <code code="" displayName="" codeSystem="" codeSystemName="" />

5065 The code shall be code='107724000' displayName='patient transfer' codeSystem='2.16.840.1.113883.6.96' />

### 6.3.4.57.7 <text><reference value="#xxx"/></text>

The <text> element shall contain a reference to the narrative text describing the transfer of the patient.

### 6.3.4.57.8 statusCode

5070 <statusCode code='normal|completed' /> When the transfer act has occurred (moodCode='EVN'), the statusCode element shall be present, and shall contain the value 'completed'. When the transfer act is intended (moodCode='EVN') the statusCode element shall contain the value 'normal'.

#### **6.3.4.57.9 <effectiveTime><low value="/"><high value="/"><effectiveTime/>**

- 5075 When the transfer has occurred, this element shall be sent, and indicates the effective time for the transfer. This element may be sent to record when the transfer act is intended to occur. The <low> element records the time at which the transfer process was started. The <high> value records the time at which the transfer was completed.

#### **6.3.4.57.10 participant**

- 5080 The <participant> element encodes the destination with a typeCode of DST  
<participant typeCode='DST'>

#### **6.3.4.57.11 templateId**

The template id identifies the facility or department which is the transfer destination.

- 5085 <templateId root='destinationLocation' />

#### **6.3.4.57.12 participantRole**

The participant role is fixed to <participantRole classCode='SDLOC'>

- 5090 **6.3.4.57.13 <id root=" extension="/>**

The <id> element shall be sent when the transfer has occurred, and identifies the performer of the act.

#### **6.3.4.57.14 <code>**

The code shall indicate the type of healthcare service location for the transfer destination.

- 5095 **6.3.4.57.15 <addr></addr>**

The <addr> element may be sent to provide a contact postal address for the performer of the disposition.

#### **6.3.4.57.16 <telecom>**

- 5100 The <telecom> element may be sent to provide a contact postal address for the performer of the disposition.

#### **6.3.4.57.17 playingEntity**

The playing entity classCode shall be ENT <playingEntity classCode='ENT'>

### 6.3.4.57.18 name

- 5105 The name element of the playing entity shall record the name of the facility or departmental destination.

*Add Section 6.4*

## 6.4 HL7 Version 2.0 Content Modules

- 5110 This section contains content modules based upon the HL7 Version 2 Standard, and related standards and/or implementation guides.

*Add Section 6.5*

## 6.5 PCC Value Sets

- 5115 This section contains value sets used by Content Modules.

Note: Although some tables in this section include a column for “Units”, units may not be applicable to all table entries and the cell will remain blank.

*Add Section 6.5.A*

### 6.5.A Antepartum History of Past Illness Value Set

#### 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.1

Name	Opt	Type	Units	SNOMED CT
Diabetes	R2	CD		73211009
Hypertension	R2	CD		38341003
Heart Disease	R2	CD		56265001
Autoimmune Disorder	R2	CD		85828009
Kidney Disease	R2	CD		90708001
UTI	R2	CD		68566005
Neurologic	R2	CD		118940003
Epilepsy	R2	CD		84757009
Psychiatric	R2	CD		74732009
Depression	R2	CD		41006004
Postpartum Depression	R2	CD		58703003
Hepatitis	R2	CD		128241005
Liver Disease	R2	CD		235856003
Varicosities	R2	CD		276504003
Phlebitis	R2	CD		61599003
Thyroid Dysfunction	R2	CD		14304000
Trauma	R2	CD		417746004

Name	Opt	Type	Units	SNOMED CT
Violence	R2	CD		225818009
History of Blood Transfusion	R2	CD		116859006
D(Rh) Sensitized	R2	CD		3885002
Pulmonary	R2	CD		19829001
Seasonal Allergies	R2	CD		367498001
Drug Allergy	R2	CD		416098002
Latex Allergy	R2	CD		300916003
Food Allergy	R2	CD		414285001
Breast	R2	CD		79604008
Hospitalizations	R2	CD		32485007
Anesthetic Complications	R2	CD		33211000
History of Abnormal Pap	R2	CD		274688009
Uterine Anomaly/DES	R2	CD		37849005
DES Exposure	R2	CD		413340008 of fetus
Infertility	R2	CD		8619003
Artificial Reproductive Therapy (ART) Treatment	R2	CD		63487001
History of Gestational Diabetes	R2	CD		
History of Incompetent Cervix	R2	CD		17382005 Code is for incompetent cervix rather than history of. Given this condition this should be okay.
History of Infant with Intrauterine Growth Restriction	R2	CD		Need Code for history of.
History of Infant with Macrosomia	R2	CD		Need Code for history of.
History of Pregnancy Induced Hypertension	R2	CD		Need code for history of.
History of Placenta Previa/Abruptio	R2	CD		Need Code for history of.
History of Preterm labor	R2	CD		441493008
History of Premature Rupture of Membranes	R2	CD		Need Code for history of.
Previous Cesarean Section	R2	CD		161805006
History of Stillbirth	R2	CD		161743003
History of Neonatal Death	R2	CD		Need code for

Name	Opt	Type	Units	SNOMED CT
				history of.
History of Postpartum Hemorrhage	R2	CD		161809000

Add Section 6.5.C

### 6.5.C Antepartum Family History and Genetic Screening Value Set

1.3.6.1.4.1.19376.1.5.3.1.1.16.5.4

Name	Opt	Type	Units	SNOMED CT	LOINC
Autism	R2	CD		408856003	
Blood Disorders	R2	CD		414022008	
Canavan Disease	R2	CD		80544005	
Chromosomal Disorder Includes any inherited genetic or chromosomal disorders	R2	CD		409709004	
Congenital Heart Defect	R2	CD		13213009	
Cystic Fibrosis	R2	CD		190905008	
Dysmorphism (Birth Defect) Patient or baby's father has a child with birth defects	R2	CD		276720006	
Down Syndrome	R2	CD		41040004	
Familial Dysautonomia	R2	CD		29159009	
Hemophilia	R2	CD		90935002	
Huntington's Chorea	R2	CD		58756001	
Maternal Metabolic Disorder	R2	CD		75934005	
Mental Retardation	R2	CD		91138005	
Muscular Dystrophy	R2	CD		73297009	
Neural Tube Defect	R2	CD		253098009	
Recurrent pregnancy loss/stillbirth	R2	CD		102878001	
Sickle Cell Disease	R2	CD		417357006	
Sickle Cell Trait	R2	CD		16402000	
Tay-Sachs	R2	CD		111385000	
Thalassemia	R2	CD		40108008	

5125

Add Section 6.5.D

## 6.5.D Antepartum Review of Systems Menstrual History Value Set

### 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.5

Name	Opt	Type	Units	SNOMED CT	LOINC
Date of Last Menstrual Period	R	TS		21840007	
Menses Monthly	R	BL		364307006	
Prior Menses Date	R	TS		21840007	
Duration of Menstrual Flow Frequency	R	PQ	days	364306002	
Frequency of Menstrual Cycles	R	PQ	days	289887006	
On Birth Control Pills at conception	R	BL		10036567	
Menarche	R	PQ		398700009	
hCG+	R	TS		250423000	

5130

*Add Section 6.5.E*

## 6.5.E Antepartum History of Infection Value Set 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.6

Name	Opt	Type	Units	SNOMED CT	LOINC
Live with someone with TB or exposed to TB	R2	CD		170464005	
History of Genital Herpes	R2	CD		402888002	
Exposed to Genital Herpes	R2	CD		240480009	
Rash since LMP	R2	CD		49882001	
Viral illness since LMP	R2	CD		34014006	
Rash or viral illness since LMP	R2	CD		49882001	
Hepatitis B	R2	CD		235871004	
Hepatitis C	R2	CD		235872006	
History of STD	R2	CD		8098009	
History of Gonorrhea	R2	CD		15628003	
History of Chlamydia	R2	CD		312099009	
History of HPV	R2	CD		302812006	
History of HIV	R2	CD		165816005	
History of Syphilis	R2	CD		76272004	

*Add Section 6.5.F*

**6.5.F Antepartum Laboratory Value Set 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.7**

Lab	LOINC Code	Comments
Antibody Screen (AB)	890-4 Ab Screen	
Blood Type (ABO Group)	883-9 ABO Group	
Rh	10331-7 Rh	
Hepatitis B virus (HBV) surface Antigen (Ag)	5196-1 HBV surface Ag (EIA) 5195-3 HBV surface Ag 5197-9 HBV surface Ag (RIA) 7905-3 HBV surface Ag (Neut)	
Hemoglobin (Hgb)/Hematocrit (Hct)	718-7 Hgb 4544-3 Hct (Automated count) 30350-3 Hgb	
Hemoglobin (Hgb) Electrophoresis	13514-5 Hemoglobin pattern [interpretation] in Blood by Electrophoresis Narrative	Appropriate code appears to be 13514-5
Aneuploidy Screening (Ultrasound)	XX-ASU Aneuploidy Screening (Ultrasound)	XX-ASU: A LOINC profile code will be requested
Pap Test/Human papilloma virus (HPV)	21440-3 HPV I/H Risk DNA Cervix (Probe) 21441-1 HPV Low Risk DNA Cervix (Probe) 10524-7 Cytology Cervix 18500-9 Thin Prep Cervix 19765-7 Cytology Cervix/Vaginal (Nominal) 19766-5 Cytology Cervix/Vaginal (Narrative)	
Rubella Virus (RUBV) Antibody (Ab)	5334-8 RUBV Ab IgG (EIA) 20458-6 RUBV Ab IgG 40667-8 RUBV Ab IgG (EIA) 8014-3 RUBV Ab IgG	
Urine Culture Screen	630-4 Bacteria Urine Culture	
Purified protein derivative (PPD)	1647-7 Purified protein derivative skin test	
Chlamydia	6347-9 Chlamydia Ag 14510-2 Chlamydia trachomatis Ag (Vaginal) 14474-1 Chlamydia trachomatis Ag (Urine) 6349-5 Chlamydia trachomatis (Unspecified specimen)	

<b>Lab</b>	<b>LOINC Code</b>	<b>Comments</b>
Gonorrhea	691-6 Neisseria Gonorrhoeae (genital specimen) 9568-7 Neisseria Gonorrhoeaea Ab	
Chlamydia Trachomatis/ Neisseria Gonorrhoeae	45067-6 Chlamydia Trachomatis Neisseria Gonorrhoeae (Cervix) 45074-2 Chlamydia Trachomatis Neisseria Gonorrhoeae (Urine)	
Ultrasound	35096-7 OB Ultrasound Panel	
Alpha-Feto Protein (Maternal) (Profile)	30525-0 Age 29463-7 Body Weight 18185-9 Gestational Age 20450-3 Alpha-1-Fetoprotein 48803-1 Neural Tube Defect Risk	
Chorionic Villus Sampling (CVS)	33774-1 Karotype	
Amniotic Fluid (Karyotype)	33773-3 Karyotype (Amino Fluid)	
Amniotic Fluid (AFP)	41273-4 Alpha-1-Fetoprotein, Amniotic Fluid Semi-Quantitative 1832-5 Alpha-1-Fetoprotein [Multiple of the median] in Amniotic Fluid 29595-6 Alpha-1-Fetoprotein [Mass/volume] in Amniotic Fluid	
Diabetes Screen	1557-8 Fasting Blood Glucose-Venous 14770-2 Fasting Blood Glucose-Capillary	
Glucose Tolerance Test (GTT)	1507-3 Glucose 1HR post 75 g glucose 14995-5 Glucose 2HR post 75 g glucose 20437-0 Glucose 3HR post 75 g glucose	
Rapid Plasma Reagin (RPR)	31147-2 Reagin Ab 20508-8 Reagin Ab by RPR	
Venereal Disease Research Laboratory (VDRL)	5292-8 Reagin Ab by VDRL	
Group B Strep	48683-7 Beta Strep Group B (PCR) 11267-2 Strep Group B	
Beta Human Chorionic Gonadotropin (HCG)	21198-7 Beta HCG	
Varicella zoster virus Ab.IgG	15410-4 Varicella zoster virus Ab.IgG (EIA)	

<b>Lab</b>	<b>LOINC Code</b>	<b>Comments</b>
	17763-4 Varicella zoster virus Ab.IgG (IF)	
Maternal Serum Triple Screen	30525-0 Age, Patient Quantitative 20450-3 Alpha-1-Fetoprotein Multiple of the Median, Serum Quantitative Calculated 20465-1 Choriogonadotropin/Choriogonatropin, Control Serum Quantitative 20466-9 Estriol/Estriol, Control Serum Quantitative	
Urinalysis (Urine Screen)	20406-5 Glucose 20505-4 Bilirubin 5797-6 Ketones 5811-5 Specific Gravity 5803-2 pH 5804-0 Protein 20405-7 Urobilinogen 20407-3 Nitrite 5794-3 Hemoglobin 5799-2 Leukocyte esterase 5767-9 Appearance 5778-6 Color 9842-6 Casts 5787-7 Epithelial cells 13945-1 Erythrocytes 5769-5 Bacteria	
First Trimester Maternal Serum Screening with Nuchal Translucency	49588-7 First trimester maternal screen with nuchal translucency [interpretation] Narrative	
Thyroid Stimulating Hormone (TSH)	11580-8 Thyrotropin (3rd generation) 3016-3 TSH 5385-0 Thyrotropin Receptor Ab 27975-2 TSH (serum)	
Triiodothyronine (T3)	3051-0 T3 Free 3052-8 T3 Reverse 3054-4 T3 True 3050-2 T3 Resin Uptake	

5135

Add Section 6.5.G
-------------------

### **6.5.G Antepartum Education Value Set 1.3.6.1.4.1.19376.1.5.3.1.1.16.5.8**

Name	Opt	Type	units	SNOMED CT	LOINC
<b>First Trimester</b>					
Risk factors identified by prenatal history	R2	CD		440047008	
Anticipated course of prenatal care	R2	CD		17629007	
Special Diet	R2	CD		171054004	
Nutrition and weight gain counseling	R2	CD		171054004	
Toxoplasmosis precautions (cats/raw meat)	R2	CD		439733009	
Sexual activity	R2	CD		162169002	
Exercise	R2	CD		171056002	
Influenza vaccine	R2	CD		xx-edu-influenza need code closest is vaccine education 171044003	
Smoking/tobacco counseling	R2	CD		171055003	
Environmental/work hazards	R2	CD		385872009	
Travel	R2	CD		439816006	
Alcohol	R2	CD		171057006	
Illicit/recreational drugs	R2	CD		425014005	
Use of any medications	R2	CD		171058001	
Indications for ultrasound	R2	CD		440227005	
Domestic violence	R2	CD		413457006	
Seatbelt use	R2	CD		440638004	
Childbirth classes/hospital facilities	R2	CD		66961001	

Name	Opt	Type	units	SNOMED CT	LOINC
<b>Second Trimester</b>					
Childbirth classes/hospital facilities	R2	CD		66961001	
Signs and symptoms of preterm labor	R2	CD		440669000	
Abnormal Lab Values	R2	CD		410299006	
Influenza vaccine	R2	CD		xx-edu-fluvaccine need code. Closest is vaccine education  171044003	
Selecting a newborn care provider	R2	CD		439908001	
Postpartum family planning	R2	CD		54070000	
Tubal sterilization	R2	CD		243064009	
<b>Third Trimester</b>					
Anesthesia/analgesia plans	R2	CD		243062008	
Intended Facility for Delivery plan				310585007	
Fetal movement monitoring	R2	CD		440309009	
Labor signs	R2	CD		440671000	
VBAC counseling	R2	CD		440073003	
Signs & Symptoms of Pregnancy-induced hypertension	R2	CD		xx-edu-sspreclampsia need to request code	
Postterm counseling	R2	CD		xx-edu-postterm need to request code	
Circumcision	R2	CD		184002001	
Bottle feeding	R2	CD		169644004	
Breast feeding	R2	CD		169643005	
Postpartum depression	R2	CD		439366005	
Newborn education (Newborn screening, jaundice, SIDS, car seat)	R2	CD		75461000	
Family medical leave or disability forms	R2	CD		40791000	
Tubal sterilization consent signed	R2	CD		408835000	

*Add section 6.5.H. (Added 2011-09 from QRPH EHCP profile)*

- 5140 The value subsets provided in this section are used both to constrain the CDA content, and to assert measure logic. These MAY be supported by the Value Set Repository actor for value set management as defined by the IHE ITI TF Sharing of Value Sets (SVS) profile.

## 6.5.H JCIH-EHDI Risk Indicators for Hearing Loss (LOINC®) Value Set

### 6.5.H.1 Metadata

5145

Metadata Element	Description	Mandatory
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.24
Name	name of the value set	JCIH-EHDI Risk Indicators for Hearing Loss (LOINC®) Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect the Risk Indicators for Hearing Loss associated with hearing loss using LOINC® concepts
Definition	A text definition describing how concepts in the value set were selected	<b>Extensional definition:</b> The value set was constructed by enumerating the codes from LOINC®
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://loinc.org">http://loinc.org</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

## 6.5.H.2 JCIH-EHDI Risk Indicators for Hearing Loss (LOINC®) Value Set Value Set

5150 *LOINC® 58232-0 Hearing Loss Risk Indicator*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.24</b>
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.1</b>
<b>Sequence</b>	<b>LOINC® Code</b>	<b>Description</b>
1	LA137-2	None
2	LA12667-4	Caregiver concern about hearing
3	LA12668-2	Family Hx of hearing loss
4	LA12669-0	NICU stay > 5 days
5	LA12670-8	ECMO
6	LA12671-6	Assisted ventilation
7	LA12672-4	Ototoxic medication use
8	LA12673-2	Exchange transfusion for Hyperbilirubinemia
9	LA12674-0	In utero infection(s)
10	LA12675-7	Craniofacial anomalies
11	LA12681-5	Physical findings of syndromes that include hearing loss
12	LA12676-5	Syndromes associated with hearing loss
13	LA12677-3	Neurodegenerative disorders
14	LA12678-1	Postnatal infections
15	LA12679-9	Head trauma
16	LA6172-6	Chemotherapy

*Add section 6.5.I. (Added 2011-09 from QRPH EHCP profile)*

## 6.5.I JCIH-EHDI Risk Indicators for Hearing Loss Codes

5155    **6.5.I.1 Metadata**

Metadata Element	Description	Mandatory
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.11
Name	name of the value set	JCIH-EHDI Risk Indicators for Hearing Loss Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect the risk indicators for hearing loss associated with hearing loss using SNOMED-CT Finding/Situation concepts
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from SNOMED-CT
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html">http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

### 6.5.I.2 JCIH-EHDI Risk Indicators for Hearing Loss Value Set

5160

#### *SNOMED-CT Risk Indicators for Hearing Loss Value Set*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.11</b>
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>
<b>Sequence</b>	<b>SNOMED-CT Code</b>	<b>Description</b>
1	439750006	Family history of hearing loss (situation)
2	441899004	History of therapy with ototoxic medication (situation)
3	276687002	Conjugated hyperbilirubinemia in infancy (disorder)
4	281610001	Neonatal hyperbilirubinemia (disorder)
5	281612009	Neonatal conjugated hyperbilirubinemia (disorder)
6	281611002	Neonatal unconjugated hyperbilirubinemia (disorder)
7	206363004	Intra-amniotic fetal infection (disorder) (Deprecated, replaced by 11618000)
8	206331005	Infections specific to perinatal period (disorder)
9	206005002	Fetus or neonate affected by maternal infection (disorder)
10	80690008	Degenerative disease of the central nervous system (disorder)
11	178280004	Postnatal infection (disorder)
12	312972009	Neonatal extracranial head trauma (disorder)
13	161653008	History of - chemotherapy (situation)
14	11618000	Intra-amniotic infection of fetus (disorder) (Replaces 206363004)

### 6.5.I.3 Pending Codes for SNOMED-CT Findings/Situation to support Risk Indicators for Hearing Loss

5165

Note that additional specificity for this value set is under way and will result in an update to this value set. Further coded values are sought to represent the following:

None
Caregiver concern about hearing
Craniofacial anomalies
Physical findings of syndromes that include hearing loss
Syndromes associated with hearing loss

*Add section 6.5.J. (Added 2011-09 from QRPH EHCP profile)*

### 6.5.J JCIH-EHDI Risk Indicators for Hearing Loss - Procedures Codes

5170

#### 6.5.J.1 Metadata

Metadata Element	Description	Mandatory

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.12
Name	name of the value set	JCIH-EHDI Risk Indicators for Hearing Loss - Procedures Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect the risk indicators for hearing loss Procedures associated with hearing loss using SNOMED-CT
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from SNOMED-CT
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html">http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

### 6.5.J.2 JCIH-EHDI Risk Indicators for Hearing Loss - Procedures Value

#### *Risk Indicators for Hearing Loss - Procedures Set*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.12</b>
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>
<b>Sequence</b>	<b>SNOMED-CT Code</b>	<b>Description</b>
1	266700009	Assisted breathing (procedure)
2	233573008	Extracorporeal membrane oxygenation (procedure)

5175

Add section 6.5.K. (Added 2011-09 from QRPH EHCP profile)

## 6.5.K Newborn Hearing Procedure Codes

### 6.5.K.1 Metadata

Metadata Element	Description	Mandatory
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.17
Name	name of the value set	JCIH-EHDI Newborn Hearing Procedure Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect the type of newborn hearing procedure identified using SNOMED-CT Procedure codes (includes both screening and other tests and examinations)
Definition	A text definition describing how concepts in the value set were selected	<b>Extensional definition:</b> The value set was constructed by enumerating the codes from SNOMED-CT
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html">http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

## 6.5.K.2 JCIH-EHDI Newborn Hearing Procedure Value Set

*Newborn Hearing Procedure Value Set:*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.17</b>
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>
<b>Sequence</b>	<b>SNOMED-CT Code</b>	<b>Description</b>
3	417491009	Neonatal hearing test (procedure)

5185

*Add section 6.5.L. (Added 2011-09 from QRPH EHCP profile)*

## 6.5.L JCIH-EHDI Newborn Hearing Screening Method Codes

### 6.5.L.1 Metadata

5190

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.4
Name	name of the value set	JCIH-EHDI Newborn Hearing Screening Method Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect the type of newborn hearing screening procedure identified using LOINC® answer codes
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from LOINC®
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://loinc.org">http://loinc.org</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY	IHE EHDI

Metadata Element	Description	Mandatory
	also have an OID assigned.	

### 6.5.L.2 JCIH-EHDI Newborn Hearing Screening Method Value Set

*Newborn Hearing Screening Method Value Set:*

*LOINC® 54106-0*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.4</b>	
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.1</b>	
Sequence	LOINC® Code	Answer Code	Description
1	LA10387-1	AABR	Automated auditory brainstem response
2	LA10388-9	ABR	Auditory brain stem response
3	LA10389-7	OAE	Otoacoustic emissions
4	LA10390-5	DPOAE	Distortion product otoacoustic emissions
5	LA10391-3	TOAE	Transient otoacoustic emissions
6	LA12406-7		Methodology unknown

5195

<i>Add section 6.5.M. (Added 2011-09 from QRPH EHCP profile)</i>
------------------------------------------------------------------

### 6.5.M JCIH-EHDI Hearing Screen Right Codes– Right

#### 6.5.M.1 Metadata

5200

Metadata Element	Description	Mandatory
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.9
Name	name of the value set	JCIH-EHDI Hearing Screen Right Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect the right ear EHDI screening using LOINC® in result type
Definition	A text definition describing how concepts in the value set were selected	<b>Extensional definition:</b> The value set was constructed by enumerating the codes from LOINC®
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://loinc.org">http://loinc.org</a>

Metadata Element	Description	Mandatory
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

### 6.5.M.2 JCIH-EHDI Hearing Screen Right Value Set

*NB hearing scn –R:Result Type*

*Hearing Screen Right Value Set*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.9</b>
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.1</b>
<b>Sequence</b>	<b>LOINC® Code</b>	<b>Description</b>
1	53109-4	Newborn Hearing Screen Right

5205

*Add section 6.5.N. (Added 2011-09 from QRPH EHCP profile)*

### 6.5.N JCIH-EHDI Hearing Screen Left Codes

#### 6.5.N.1 Metadata

Metadata Element	Description	Mandatory
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.8
Name	name of the value set	JCIH-EHDI Hearing Screen Left Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect the left ear EHDI hearing screening result type using LOINC®
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from LOINC®
Source URI	Most sources also have a URL or	<a href="http://loinc.org">http://loinc.org</a>

Metadata Element	Description	Mandatory
	document URI that provides further details regarding the value set.	
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

### 6.5.N.2 JCIH-EHDI Hearing Screen Left Value Set

5210 *Hearing Screen Left Value Set*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.8</b>
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.1</b>
<b>Sequence</b>	<b>LOINC® Code</b>	<b>Description</b>
1	53108-6	Newborn Hearing Screen Left

*Add section 6.5.O. (Added 2011-09 from QRPH EHCP profile)*

### 6.5.O JCIH-EHDI Reason for no Hearing Loss Diagnosis or Screening Codes(SNOMED)

5215 **6.5.O.1 Metadata**

Metadata Element	Description	Mandatory
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.15
Name	name of the value set	JCIH-EHDI Reason for no Hearing Loss Diagnosis or Screening Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general	To Reflect Reason for no hearing loss

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
	purpose of the value set	diagnosis coded with SNOMED-CT.
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from SNOMED-CT
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html">http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

### 6.5.O.2 JCIH-EHDI Reason for no Hearing Loss Diagnosis or Screening Value Set

#### *Reason for no Hearing Loss Diagnosis*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.15</b>	
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>	
<b>Sequence</b>	<b>SNOMED Code</b>	<b>Description</b>	<b>EHDI Concept</b>
1	397709008	Patient died (finding)	No screening or diagnosis: Infant died
2	360885002	Change of residence status (finding)	No diagnosis: Moved or gone elsewhere
3	184112005	Patient address unknown (finding)	No diagnosis: Unable to Contact Family
4	184118009	Patient telephone number unknown (finding)	No diagnosis: Unable to Contact Family
5	183638004	Follow-up refused	No screening diagnosis: Parents Declined Services - Follow-up refused
6	183946001	Procedure refused-uncooperative	No diagnosis: Parents Declined Services -Procedure refused - uncooperative
7	413319007	Persistent non-attender	No diagnosis: Unresponsive - Persistent non-attender

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.15</b>	
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>	
<b>Sequence</b>	<b>SNOMED Code</b>	<b>Description</b>	<b>EHDI Concept</b>
8	399307001	Loss to follow-up	No diagnosis: Unknown - Loss to follow-up
9	419984006	Inconclusive (qualifier value)	No diagnosis: Audiologic Diagnosis in Process
10	185332005	Appointment cancelled by patient (finding)	No diagnosis: Audiologic Diagnosis in Process - Rescheduled appointment
11	185333000	Appointment cancelled by doctor (finding)	No diagnosis: Audiologic Diagnosis in Process - Rescheduled appointment
12	281399006	Did not attend	No diagnosis: Audiologic Diagnosis in Process - Did not attend

5220

Add section 6.5.P. (Added 2011-09 from QRPH EHCP profile)

## 6.5.P JCIH-EHDI Newborn Hearing Loss Referrals Codes

### 6.5.P.1 Metadata

5225

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.16
Name	name of the value set	JCIH-EHDI Newborn Hearing Loss Referrals Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect EHDI referrals coded with SNOMED-CT and as a response to care plan recommendations (entered on a list of referrals in a medical summary)
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from SNOMED-CT
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html">http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html</a>
Version	A string identifying the specific version of the value set.	Version 1.0

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

### 6.5.P.2 JCIH-EHDI Newborn Hearing Loss Referrals Value Set

#### *EHDI Newborn Hearing Loss Referrals Value Set*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.16</b>	
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>	
<b>Sequence</b>	<b>SNOMED Code</b>	<b>Description</b>	<b>EHDI Concept</b>
1	306210008	Referral to pediatric diagnostic audiology service (procedure)	Referral to audiologist
2	415271004	Referral to education service (procedure)	Referral to Early Intervention (Part C/non Part C)

*Add section 6.5.Q. (Added 2011-09 from QRPH EHCP profile)*

- 5230    **6.5.Q JCIH-EHDI Newborn Hearing Loss Reason for no Follow-up – Patient Reason Codes**

### 6.5.Q.1 Metadata

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.7
Name	name of the value set	JCIH-EHDI Newborn Hearing Loss Reason for no Follow-up – Patient Reason Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect the reason that no follow-up is conducted in the case of hearing loss using SNOMED-CT reflected in

Metadata Element	Description	Mandatory
		negation of an intent to order the referral.
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from SNOMED-CT
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html">http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

5235    **6.5.Q2 JCIH-EHDI Newborn Hearing Loss Reason for no Follow-up – Patient Reason Value Set**

*EHDI Newborn Hearing Loss Reason for no Follow-up Value Set*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.7</b>	
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>	
Sequence	<b>SNOMED Code</b>	<b>Description</b>	<b>EHDI Concept</b>
1	397709008	Patient died (finding)	Incomplete outpatient screen: Infant died
2	360885002	Change of residence status (finding)	Incomplete outpatient screen: Moved or gone elsewhere
3	184112005	Patient address unknown (finding)	Incomplete outpatient screen: Unable to contact family
4	184118009	Patient telephone number unknown (finding)	Incomplete outpatient screen: Unable to contact family
5	183638004	Follow-up refused	Incomplete outpatient screen: Follow-up refused
6	183946001	Procedure refused-uncooperative	Incomplete outpatient screen: Procedure refused-uncooperative
7	413319007	Persistent non-attender	Incomplete outpatient screen: Persistent non-attender

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.7</b>	
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>	
<b>Sequence</b>	<b>SNOMED Code</b>	<b>Description</b>	<b>EHDI Concept</b>
8	399307001	Loss to follow-up	Incomplete outpatient screen: Loss to follow-up
9	185332005	Appointment cancelled by patient (finding)	Incomplete outpatient screen: Rescheduled appointment
10	185333000	Appointment cancelled by doctor (finding)	Incomplete outpatient screen: Rescheduled appointment
11	281399006	Did not attend	Incomplete outpatient screen: Did not attend

5240 *Add section 6.5.R. (Added 2011-09 from QRPH EHCP profile)*

## 6.5.R Joint Commission Medical Reason Codes

### 6.5.R.1 Metadata

Metadata Element	Description	Mandatory
Identifier	unique identifier of the value set	1.3.6.1.4.1.33895.1.3.0.75
Name	name of the value set	Joint Commission Medical Reason Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	The Joint Commission value set is used to reflect medical reason why a test was not performed
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from SNOMED-CT
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html">http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

### 6.5.R.2 Joint Commission Medical Reason Value Set

5245 EHDI specifies the re-use of the existing Medical Reason Value Set used by the Joint Commission measures.

	<b>Value Set :</b>	<b>1.3.6.1.4.1.33895.1.3.0.75</b>
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>
<b>Sequence</b>	<b>SNOMED-CT Code</b>	<b>Description</b>
1	397745006	Medical contraindication (finding)
2	397773008	Surgical contraindication (finding)

<i>Add section 6.5.S. (Added 2011-09 from QRPH EHCP profile)</i>
------------------------------------------------------------------

5250 **6.5.S JCIH-EHDI Inpatient Screening Results not Performed Codes****6.5.S.1 Metadata**

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.10
Name	name of the value set	JCIH-EHDI Inpatient Screening Results not Performed Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect the right ear EHDI results reported using LOINC® answer lists
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from LOINC®
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://loinc.org">http://loinc.org</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

**6.5.S.2 JCIH-EHDI Inpatient Screening Results not Performed Value Set**

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1. 1.15.2.10</b>		
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.1</b>		
<b>Sequence</b>	<b>LOINC®</b>	<b>Description</b>	<b>Global ID</b>	<b>Global ID</b>

	<b>Code</b>			<b>Code System</b>
1	LA12408-3	Attempted, but unsuccessful - technical fail	103709008	SN
2	LA7304-4	Not performed	262008008	SN
3	LA12409-1	Not performed, medical exclusion - not indicated	410534003	SN

5255

Add section 6.5.T. (Added 2011-09 from QRPH EHCP profile)

## 6.5.T JCIH-EHDI Evidence of Hearing Screening Performed Codes

### 6.5.T.1 Metadata

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.18
Name	name of the value set	JCIH-EHDI Evidence of Hearing Screening Performed Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect Evidence of Hearing Screening Performed through the result values of pass-Left, pass-Right, or Refer. This excludes unsuccessful results.
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from LOINC®
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://loinc.org">http://loinc.org</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

5260

**6.5.T.2 JCIH-EHDI Evidence of Hearing Screening Performed Value Set***Evidence of Hearing Screening Performed Value Set*

		<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.18</b>		
		<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.1</b>		
<b>Sequence</b>	<b>LOINC® Code</b>	<b>Answer Code</b>	<b>Description</b>	<b>Global ID</b>	<b>Global ID Code System</b>
1	LA10392-1	164059009	Pass		
2	LA10393-9	183924009	Refer		

*Add section 6.5.U. (Added 2011-09 from QRPH EHCP profile)*

5265

**6.5.U JCIH-EHDI Procedure Declined Value Set Codes****6.5.U.1 Metadata**

<b>Metadata Element</b>	<b>Description</b>	<b>Mandatory</b>
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.20
Name	name of the value set	JCIH-EHDI Procedure Declined Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general purpose of the value set	To Reflect that the hearing screening procedure was not performed due to the patient/parent declining the procedure
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from SNOMED-CT
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html">http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html</a>

Metadata Element	Description	Mandatory
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

### 6.5.U.2 JCIH-EHDI JCIH-EHDI Procedure Declined Value Set Value Set

*JCIH-EHDI Procedure Declined Value Set:*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.20</b>
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>
Sequence	<b>SNOMED-CT Code</b>	<b>Description</b>
1	183949008	Assessment examination refused (situation)
2	183945002	Procedure refused - religion (situation)
3	183948000	Refused procedure - parent's wish (situation)
4	397709008	Patient died (finding)

5270

<i>Add section 6.5.V. (Added 2011-09 from QRPH EHCP profile)</i>
------------------------------------------------------------------

### 6.5.V JCIH-EHDI Newborn Hearing Screening Abnormal Results Value Set Codes

#### 6.5.V.1 Metadata

Metadata Element	Description	Mandatory
Identifier	unique identifier of the value set	1.3.6.1.4.1.19376.1.7.3.1.1.15.2.23
Name	name of the value set	JCIH-EHDI Newborn Hearing Screening Abnormal Results Value Set
Source	source of the value set, identifying the originator or publisher of the information	IHE Quality Research and Public Health Domain
Purpose	Brief description about the general	To Reflect abnormal results from the

Metadata Element	Description	Mandatory
	purpose of the value set	hearing screening procedure
Definition	A text definition describing how concepts in the value set were selected	Extensional definition: The value set was constructed by enumerating the codes from SNOMED-CT
Source URI	Most sources also have a URL or document URI that provides further details regarding the value set.	<a href="http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html">http://www.nlm.nih.gov/research/umls/Snomed/snomed_main.html</a>
Version	A string identifying the specific version of the value set.	Version 1.0
Status	Active (Current) or Inactive	Active
Effective Date	The date when the value set is expected to be effective	8/1/2010
Expiration Date	The date when the value set is no longer expected to be used	N/A
Creation Date	The date of creation of the value set	8/1/2010
Revision Date	The date of revision of the value set	N/A
Groups	The identifiers of the groups that include this value set. A group MAY also have an OID assigned.	IHE EHDI

5275 **6.5.V.2 JCIH-EHDI Newborn Hearing Screening Abnormal Results Value Set**

*JCIH-EHDI Newborn Hearing Screening Abnormal Results Value Set:*

	<b>Value Set :</b>	<b>1.3.6.1.4.1.19376.1.7.3.1.1.15.2.23</b>
	<b>Vocabulary:</b>	<b>2.16.840.1.113883.6.96</b>
<b>Sequence</b>	<b>SNOMED-CT Code</b>	<b>Description</b>
1	313203003	Hearing test abnormal (finding)
2	308409008	Child hearing screening failure (finding)
3	185577006	Child hearing screening first failure (finding)
4	185579009	Child hearing screening second failure (finding)
5	185580007	Child hearing screening failure referred to specialist (finding)

<i>Add section 6.5.W. (Added 2011-09 from QRPH PRPH-Ca profile)</i>
---------------------------------------------------------------------

**6.5.W Primary Site Value Set**

LOINC = 22035-0	
Code System: ICD-O-3 2.16.840.1.113883.6.43.1	
Code	Meaning
	A code from ICD-O-3 (Topography Section)

5280

<i>Add section 6.5.X (Added 2011-09 from QRPH PRPH-Ca profile)</i>
--------------------------------------------------------------------

### 6.5.X Histologic Type Value Set

LOINC = 31205-8	
Code System: ICD-O-3 2.16.840.1.113883.6.43.1	
Code	Meaning
	An ICD-O-3 code (Morphology Section)

5285

<i>Add section 6.5.Y (Added 2011-09 from QRPH PRPH-Ca profile)</i>
--------------------------------------------------------------------

### 6.5.Y Derived AJCC Descriptor (T,N,M) Value Set

LOINC = 21908-9	
Code System: 2.16.840.1.113883.15.6	
Code	Meaning
c	clinical
p	pathologic
a	Autopsy classification
yc or yp	Posttherapy classification “y” prefix to utilize with “c” or “p” for denoting extent of cancer after neoadjuvant or primary systemic and/or radiation therapy
r	Retreatment Classification

5290

<i>Add section 6.5.Y (Added 2011-09 from QRPH PRPH-Ca profile)</i>
--------------------------------------------------------------------

### 6.5.Z TNM Edition Value Set

LOINC = 21917-0	
Code System: 2.16.840.1.113883.15.6	
Code	Meaning
5	AJCC Staging Manual, 5 <sup>th</sup> Edition
6	AJCC Staging Manual, 6 <sup>th</sup> Edition
7	AJCC Staging Manual, 7 <sup>th</sup> Edition

5295

*Add section 6.5.Z (Added 2011-09 from QRPH PRPH-Ca profile)***6.5.AA TNM Stage Group Value Set**

5300

Note: The AJCC Staging Manual TNM system is propriety and its definitions cannot be included in other documents without permission.

LOINC = 21908-9	
Code System: TNM 5. Edition: 2.16.840.1.113883.15.8 - tnm5 TNM 6. Edition: 2.16.840.1.113883.15.7 - tnm6 TNM 7. Edition: 2.16.840.1.113883.15.6 - tnm7	
<b>Code</b> <b>Description: Site specific descriptions prevent listing of text equivalents.</b>	
0	Site specific descriptions prevent listing of text equivalents.
0a	“
0is	“
I	“
IA	“
IA1	“
IA2	“
IB	“
IB1	“
IB2	“
IC	“
II	“
IIA	“
IIA1	“
IIA2	“
IIB	“
IIC	“
III	“
IIIA	“
IIIB	“
IIIC	“
IS	“
IV	“
IVA	“
IVB	“
IVC	“

<i>Add section 6.5.AA (Added 2011-09 from QRPH PRPH-Ca profile)</i>
---------------------------------------------------------------------

5305

### 6.5.BB TNM Stage Descriptor Value Set

Note: The AJCC Staging Manual TNM system is propriety and its definitions cannot be included in other documents without permission.

LOINC = 21909-7	
Code System: TNM 5. Edition: 2.16.840.1.113883.15.8 - tnm5 TNM 6. Edition: 2.16.840.1.113883.15.7 - tnm6 TNM 7. Edition: 2.16.840.1.113883.15.6 - tnm7	
Code	Meaning
0	None
1	E (Extranodal, lymphomas only)
2	S (Spleen, lymphomas only)
3	M (Multiple primary tumors in a single site)
4	Y (Classification during or after initial multimodality therapy)—pathologic staging only
5	E & S (Extranodal and spleen, lymphomas only)
6	M & Y (Multiple primary tumors and initial multimodality therapy)

5310

<i>Add section 6.5.BB (Added 2011-09 from QRPH PRPH-Ca profile)</i>
---------------------------------------------------------------------

### 6.5.CC TNM Tumor Value Set

5315

Note: The AJCC Staging Manual TNM system is propriety and its definitions cannot be included in other documents without permission.

LOINC = 21905-5	
Code System: TNM 5. Edition: 2.16.840.1.113883.15.8 - tnm5 TNM 6. Edition: 2.16.840.1.113883.15.7 - tnm6 TNM 7. Edition: 2.16.840.1.113883.15.6 - tnm7	
Code	Description: Site specific descriptions prevent listing of text equivalents.
Ta	Site specific descriptions prevent listing of text equivalents.
Tis	“
T0	“
T1	“
T1mic	“
T1a	“

T1a1	“
T1a2	“
T1b	“
T1b1	“
T1b2	“
T1c	“
T1d	“
T2	“
T2a	“
T2a1	“
T2a2	“
T2b	“
T2c	“
T2d	“
T3	“
T3a	“
T3b	“
T3c	“
T3d	“
T4	“
T4a	“
T4b	“
T4c	“
T4d	“
T4e	“
Tx	“

Add section 6.5.CC (Added 2011-09 from QRPH PRPH-Ca profile)

## 6.5.DD TNM Node Value Set

5320

Note: The AJCC Staging Manual TNM system is propriety and its definitions cannot be included in other documents without permission.

LOINC = 21906-3	
Code System: TNM 5. Edition: 2.16.840.1.113883.15.8 - tnm5 TNM 6. Edition: 2.16.840.1.113883.15.7 - tnm6 TNM 7. Edition: 2.16.840.1.113883.15.6 - tnm7	
<b>Code</b>	<b>Description: Site specific descriptions prevent listing of text equivalents.</b>
N0	Site specific descriptions prevent listing of text equivalents.
N1	“
N1mi	“
N1a	“

N1b	“
N1b1	“
N1b2	“
N1b3	“
N1b4	“
N1c	“
N2	“
N2a	“
N2b	“
N2c	“
N3	“
N3a	“
N3b	“
N3c	“
N	“

Add section 6.5.DD (Added 2011-09 from QRPH PRPH-Ca profile)

## 6.5.DD TNM Metastasis Value Set

5325 Note: The AJCC Staging Manual TNM system is propriety and its definitions cannot be included in other documents without permission.

LOINC = 21907-1	
Code System: TNM 5. Edition: 2.16.840.1.113883.15.8 - tnm5	
TNM 6. Edition: 2.16.840.1.113883.15.7 - tnm6	
TNM 7. Edition: 2.16.840.1.113883.15.6 - tnm7	
Code	Description: Site specific descriptions prevent listing of text equivalents.
M0	Site specific descriptions prevent listing of text equivalents.
M1	“
M1a	“
M1b	“
M1c	“
M1d	“
M1e	“
Mx	“

Add Appendix Q

5330

## **APPENDIX Q: Document Construction**

Describe document construction.